

00040.030240

PATENT APPLICATION

IN THE UNITED STATES PATENT AND TRADEMARK OFFICE

In re Application of:	)	
GERHARD HOEFLE, ET AL.	:	Examiner: Taofiq A. Solola
	)	
Application No.: 09/313,524	:	Group Art Unit: 1625
	)	
Filed: May 17, 1999	:	Confirmation No.: 4030
	)	
For: EPOTHILONES C, D, E, AND F,	:	
PREPARATION AND	)	
COMPOSITIONS	:	

Commissioner for Patents  
P.O. Box 1450  
Alexandria, VA 22313-1450

DECLARATION UNDER 37 C.F.R. §1.131

Sir:

Pursuant to the provisions of MPEP 715, Applicants hereby aver as follows

1. I, Gerhard Hoefle, am a named inventor of the above-identified application.
2. I actually reduced to practice various species within the scope of the claims, or supervised the actual reduction to practice various species within the scope of the claims, before September 26, 1996. Copies of various laboratory notebook pages confirming this are attached. Dates on the attached photocopies are redacted, but as to paragraphs 1- 59, they are all prior to September 26, 1996, with the balance at least occurring prior to filing, as developed in the course of Interference No. 105,298. For the Examiner's convenience, an initial table is provided as a summary sheet (labeled

“Evidence for the discovery of epothilones C and D from...” as well. Note too that the laboratory notebook pages have been labeled “Exhibit 3-1, 3-2, etc.” by me simply for easy reference herein.

3. Dr. Klaus Gerth, a microbiologist at Gesellschaft für Biotechnologische Forschung mbH (“GBF”), instructed Carmen Fischer, a laboratory technician in the Chemistry Group at GBF, to commence cultivation and screening using a culture sample of *Sorangium cellulosum*, Socel198 45/30.

4. Accordingly, Ms. Fischer took a sample of this Socel198 45/30 strain, and inoculated it in heso; she recorded her work on a Sample Table, Exhibit 3-1, and in particular, she recorded the growth of this strain to be “good.” Id.

5. Likewise, Ms. Fischer took a second sample of this Socel198 45/30 strain, inoculated it in Probion (this is a protein supplied by, at that time, Hoechst AG as nutrient component), and took the resultant product and used it to inoculate a 250 ml shaking flask with medium S (S is the lab term for a culture medium used for the screening of *Sorangium* strains). Id.

6. Ms. Fischer used that material in turn to inoculate three 250 ml flasks with medium S, which she recorded was harvested. Exhibits 3-1 and 3-2.

7. Ms. Fischer recorded that two (“second” is a misinterpretation of “2 K.” in the English translation) of the three 250 ml flasks were “good.” Exhibit 3-1.

8. The product harvested by Ms. Fischer was then given to Dr. Gerth, who proceeded to conduct an HPLC-UV absorbance analysis on the product – this analysis

was a “production check”; i.e., to check whether the strain was producing the desired materials.

9. In particular, Dr. Gerth injected the product onto an HPLC column and ran a UV absorbance analysis on the eluent, and the resultant spectrum is found in Exhibit 3-3.

10. On the spectrum found in Exhibit 3-3, Dr. Gerth wrote “Epo A” above the UV absorbance peak for the material eluted at 13.282 minutes, and wrote “Epo B” above the absorbance peak for the material eluted at the 14.584 minutes. Id. at page 1.

11. Dr. Gerth was able to identify these materials as epothilone A and epothilone B based upon his previous work with these materials.

12. His identification was also confirmed by the two UV peaks observed at the 13.282 and 14.584 time slices, characteristic of the presence of the thiazole side-chain (only thiazoles with the adjacent double bond show this peak pattern) found in epothilones A and B. See id. at page 4.

13. Additionally, Dr. Gerth noted at the same day the presence of an additional material or materials, which he recorded as “epo unbekannt” (meaning “unknown epothilones”) on the spectrum. Id. at page 1.

14. Dr. Gerth identified this material in this manner because it also exhibited the characteristic UV bands of epothilone A and epothilone B, but plainly was not those materials, since it eluted at a different time. [See id. at page 6.]

15. C. Fischer summarized these results on the Sample Table (Exhibit 3-1), where she recorded the names “epothilon A,” “epothilon B,” and “epo. Unbekannt” adjacent their respective peaks.

16. Ms. Fischer provided Mr. Steinmetz with a 20 microliter (erroneously “ml” in the english translation) into sample of So cel 198 45/30 in a high pressure liquid chromatography (“HPLC”) tube for analysis. Exhibit 3-4.

17. Mr. Steinmetz in turn gave the sample to Ms. Antje Ritter, a laboratory technician who also worked in the Chemistry group, and asked her to analyze the sample of Socel 198 using HPLC/UV and on-line mass spectrometry (“MS”).

18. Ms. Ritter conducted the requested HPLC-UV analysis. Exhibit 3-5. In the same HPLC run she also conducted the MS analysis on the sample. Id.

19. The same day that the HPLC/UV/MS analysis was conducted, Ms. Pohlan reviewed the print-outs from the HPLC/UV/MS analysis and wrote down “Epo A” and “Epo B” next to the peaks in the chromatogram (page 1) and next to their corresponding UV spectra (pages 2 and 3); she additionally wrote down “Epo neu” (“new epothilones”) next to the two peaks in the chromatogram as well as UV spectra corresponding to the newly identified, but as yet unisolated and uncharacterized, materials. Id.

20. Mr. Steinmetz reviewed the results of the ESI-MS analysis, and found that the two new materials showed protonated molecular ion ( $M+H^+$ ) of 478.2 and



492.5 respectively, as compared to those of epothilone A, 494.4, and epothilone B, 508.6. Exhibit 3-5.

21. Mr. Steinmetz therefore realized that the new materials each differed from epothilones A and B by the atomic mass 16, which is the mass of an oxygen atom.

22. In addition, Mr. Steinmetz reviewed with Ms. Pohlan the HPLC/UV/MS results, and noted that the new materials, which exhibited the characteristic UV bands of epothilones, each had an atomic mass 16 less than epothilones A and B; accordingly, Ms. Pohlan then recorded "Epo neu" (meaning "new epothilone") adjacent the UV absorbance curves taken for the 29.25 and 30.57 time slices. Exhibit 3-5.

23. Mr. Steinmetz discussed with Ms. Pohlan the 16 amu mass difference, and she recorded below the "Epo new" entries the equation " $mz = -16$ ," meaning that these new materials differed in mass from epothilones A and B by 16, respectively. Exhibit 3-5, at page 3.

24. Upon analyzing the data, Dr. Hoefle and Mr. Steinmetz concluded that the new materials had the same structures as epothilones A and B, except that they were missing an oxygen atom, presumably the epoxide group (one oxygen atom, atomic mass 16).

25. A departmental meeting was attended by Drs. Hoefle, Gerth, Reichenbach, Mr. Steinmetz and others, and was recorded in Meeting Minutes prepared by Dr. Reichenbach. Exhibit 3-6.

26. At the departmental meeting, Dr. Gerth reported his finding that the *Sorangium cellulosum* Soce 1198 strain produced epothilones A and B, and additionally produced small quantities of two unknown compounds exhibiting the characteristic ultraviolet (UV) spectrum of epothilone A and B, and that the new compounds were more lipophilic than epothilones A and B.

27. The “two peaks” referred to by Dr. Gerth at the departmental meeting were the thiazole double peaks exhibited by the UV absorbance spectrum.

28. In the Minutes of the meeting Dr. Höfle wrote:

According to HPLC/MS studies carried out by Herr Steinmetz, the substance is composed of homologues ( $\Delta 14$ ) possessing one oxygen atom less than epothilone A and B. There were ca. 1-2 mg of the new epothilones present in the shaken culture.

Exhibit 3-6 (English translation).

29. Dr. Höfle’s reference to “homologues” was intended to indicate that the new compounds had a similar structure to the known compounds epothilone A and epothilone B, except for the absence of an oxygen atom. Id.

30. Dr Höfle recorded that the next step was to isolate the new compounds individually, or as a mixture. Id.

31. Ms. Pohlen received from Dr. Gerth a methanol extract of adsorber resin collected by Dr. Gerth from the shaking flasks in the screening, which was labeled “So cel 1198-45/30, Screening.” Exhibit 3-7.

32. Ms. Pohlan evaporated the MeOH off, and recorded that she had obtained 198 mg ( "g", gram is an error in the English translation) of material. On the same day she run an analytical tlc (on the left hand side of Exhibit 3-7) comparing the extract with authentic epothilone A. Exhibit 3-7.

33. Ms. Pohlan fractionated the material using a Sephadex LH-20 column 1.5 cm in diameter and 70 cm long, and obtained six fractions, which she recorded as "LH-1" through "LH-6," respectively. Exhibit 3-7.

34. Ms. Pohlan spotted these fractions onto different positions of a thin layer plate, and she recorded that the developed chromatogram showed spots in fraction "LH-2," which indicated that the epothilones A and B as well as the new compounds would be present in that fraction. Exhibit 3-7.

35. Ms. Pohlan conducted an analytical check on fraction LH-2 by subjecting a small sample to reverse phase HPLC separation and UV detection. Exhibit 3-9.

36. The UV absorption trace of the analytical check on fraction LH-2 displayed in the 2.43 time slice the characteristic UV spectrum of epothilones which confirmed the presence of epothilone A and epothilone B. Exhibit 3-9.

37. The UV absorption trace of the analytical check on fraction LH-2 further displayed in the 4.26 and 5.04 time slices the characteristic UV spectrum of epothilones, thus confirming that this sample also contained the new compounds. Exhibit 3-9.

38. Ms. Pohlen recorded in her notebook that she had conducted a reverse phase ("RP") chromatographic fractionation on sample LH-2, using a Nucleosil 100 column (20 x 250 mm) and a solvent consisting of 73 parts methanol and 27 parts water. Exhibit 3-10.

39. Ms. Pohlen also obtained a UV absorbance trace of the eluent, and on that trace she noted that fractions 35 to 37 exhibited the peaks corresponding to epothilone A and epothilone B, which she labeled "epo A" and "epo B" respectively on the trace. Exhibit 3-11.

40. However, Ms. Pohlen also noted a UV absorption peak spanning fractions [46 and 47], which she labeled "RP-1" on the trace, and another UV absorption peak spanning fractions 50 and 51, which she labeled "RP-2" on the trace.

41. These UV absorbance peaks were believed to have been produced by the new compounds, and thus eluted fractions 46 and 47 were believed to contain one of the new compounds, and eluted fractions 50 and 51 were believed to contain the other of the new compounds.

42. Ms. Pohlen used a thin layer chromatograph technique to analyze a small amount of each of the eluents corresponding to peaks RP-1 and RP-2, and she recorded that the resultant spots had a violet color, which is the color that epothilones A and B were known to develop after spraying with vanillin/sulfuric acid. Exhibit 3-10.

43. Dr. Hoefle instructed Ms. Pohlen to submit the fraction corresponding to peak RP-2 for NMR analysis. Id.

44. Ms. Pohlen submitted the sample with an NMR Request Form, requesting a proton analysis (standard spectrum and COSY, 1D and 2D, respectively); the NMR analyses were conducted on the same day, and the resultant standard spectrum was given Spectrum no. 2550. Exhibit 3-12.

45. Spectrum no. 2550 for sample RP-2 was given to Dr. Hoefle; he reviewed it and found it to have characteristics of epothilone B, such as the five methyl singlets in the range of 1 to 2.7 ppm, and an olefinic singlet around 6.6 ppm. Exhibit 3-12.

46. However, Dr. Hoefle noted the presence in the NMR spectrum for sample RP-2 of a singlet at about 1.7 ppm; if the substance were epothilone B, this singlet would have been present at about the 1.2 ppm position.

47. Thus the singlet location was shifted, relative to epothilone B.

48. Given the previously recognized mass difference of 16, which is the weight of an atom of oxygen, Dr. Hoefle attributed the singlet shift in the RP-2 sample, relative to epothilone B, to the presence of a double bond, which had replaced the epoxide group.

49. After reviewing the NMR print-out, Dr. Hoefle sketched on it the CH<sub>3</sub> structure that he attributed to the singlet on the NMR print-out. See Exhibit 3-12.

50. Dr. Hoefle also drew a more complete picture of the molecular structure of the RP-2 material on the COSY NMR spectrum, which shows in addition to the structure of the -CH=CHCH<sub>3</sub>- group that replaced the epoxide group, the surrounding partial structures of epothilone. Exhibit 3-12.

51. The RP-2 eluent material, which Dr. Hoefle structurally characterized in the manner explained above, is the material that Dr. Hoefle and Mr. Steinmetz named "epothilone D."

52. As to the eluent material identified as RP-1, Dr. Hoefle noticed from the chromatograph (Exhibit 3-11) that it sat on a broader peak, which would suggest that it was mixed with other material.

53. Accordingly, Dr. Hoefle directed that the material RP-1 be further purified, and Ms. Pohlan subjected sample RP-1 to separation on silica gel plate. Exhibit 3-13.

54. Ms. Pohlan then used a thin layer chromatographic technique to analyze the resultant RP-1/DC, (DC means purified by thin-layer chromatography), and recorded that it exhibited a single band only, indicating a purified product. Id.

55. The next day, Ms. Pohlan submitted the purified sample of RP-1 for NMR proton analysis, and the resultant NMR was given Spectrum no. 2630. Exhibit 3-15.

56. Ms. Pohlan showed the NMR print-out for sample RP-1 to Dr. Hoefle, and he was immediately able to characterize the structure of the material, which he drew on the right-hand side of the NMR print-out. Exhibit 3-15.

57. The RP-1 eluent material, which Dr. Hoefle structurally characterized in the manner explained above, is the material that Dr. Hoefle and Mr. Steinmetz named "epothilone C."

58. A departmental meeting was held which was attended by a number of people, including Drs. Hoefle, Reichenbach, Sasse and Mr. Steinmetz and was recorded in Meeting Minutes prepared by Dr. Reichenbach. Exhibit 3-16.

59. At the meeting, Dr. Gerth, Mr. Steinmetz and Dr. Sasse reported to the attendees the following:

The strains So ce1198, So ce1275 and So ce1294 form two new epothilones as well as epothilone, but with the epoxide missing (Gerth, Steinmetz). They had considerably reduced action, but were not abolished: the IC50 for L929 cells was 150 ng/ml for RP1 (from So ce1198), and 100 ng/ml for RP2. Noticeable effect on Tubulin could be detected in cell cultures. (Sasse)  
Perhaps patenting is possible?

Exhibit 3-16 (English translation).

60. Notably, the above, initial work was completed using strains So ce1198, 1275 and 1294. Further isolation work of epothilone C and D was then conducted using a variant or mutant strain of *Sorangium cellulosum*, So ce90. The wild version of So ce90 had previously been deposited with the German Collection for Microorganisms ("Deutsche Sammlung von Mikroorganismen") as DSM 6773.

61. A number of cultures were prepared from DSM 6773. These cultures, which were called "clones" generally do not have the same population mixture or production profile as DSM 6773.

62. In particular, So ce90 A3, which based on earlier work was known to be a good producer of epothilone A and epothilone B, was used for this further isolation work since its production profile was similar to that of So ce1198. Ms. Fischer recorded

the preparation of cultures medium for 15 L and 150 L fermentors (having working volumes of 10 L and 100 L respectively). Exhibits 4-10 to 4-12. The details and monitoring of these fermentations were also recorded, including a description of the fermentation medium used. Exhibits 4-13 to 4-26.

63. The product of these fermentations was then used to charge a 750 L fermenter. The details and monitoring of this fermentation was also recorded, including a description of the fermentation medium used. Exhibit 4-27 to 4-34.

64. The harvest of this fermentation was then undertaken by recovering the XAD absorber resin from the 750 L fermenter, filtering the absorber resin, (Exhibit 4-35 to 4-36); eluting the absorber resin with methanol, (Exhibit 4-37 to 4-38); concentrating the eluent by evaporation to a 20 L concentrate, (Exhibit 4-39 to 4-40); performing an ethyl acetate extraction, (Exhibit 4-41 to 4-42); and then subjecting the extract to rotary evaporation to yield crude extract, (Exhibit 4-43 to 4-44).

65. The crude extract was next tested for the presence of epothilone A, B, C and D. See Exhibits 4-45 to 4-49, particularly the mass spectrometer results shown in Exhibit 4-49, which depict the peaks indicative of the presence of these species.

66. The crude extract was then dried, distributed between methanol and heptane (the heptane was discarded), (Exhibit 4-50 to 4-51); and passed through a Sephadex LH20 chromatographic column, (Exhibit 4-52). Fractions were collected and utilized for further analysis. Exhibits 4-53 to 4-57. Again, the mass spectrometer results shown in Exhibit 4-57 exhibited peaks indicating epothilone A, B, C and D were present.



67. A reverse phase chromatography was next performed using fractions 6-12. Exhibits 4-58 to 4-60. A UV absorbance analysis indicated that fractions 8-12 contained epothilone A and B. Exhibits 4-61 to 4-63. Fraction 24, on the other hand, was subject to HPLC MS analysis, and was found to exhibit a peak indicating the presence of epothilone C. Exhibit 4-64. Similarly, fraction 28 was found to exhibit a peak indicating the presence of epothilone D. Exhibit 4-66. Mass spectrometer analyses confirmed these results. Exhibits 4-65, 4-67.

68. The epothilone C and D in the fractions referenced above were next purified using reverse phase RP-18 chromatography, and then analyzed. Exhibits 4-68 to 4-75. In particular, fraction RP-1 was subjected to UV absorbance analysis and exhibited a clean peak, indicating that it contained pure epothilone D. Exhibit 4-74. Fraction RP-2, which was epothilone C, was subject to TLC analysis. That analysis showed it to be free of trace contaminants. Exhibit 4-71.

69. There followed an NMR analysis of fraction RP-2, which confirmed the peaks as those of epothilone C. Exhibit 4-78 to 4-91. An NMR analysis was also conducted of fraction RP-1, which confirmed the peaks as those of epothilone D. Exhibit 4-92 to 4-105.

70. The data from the tests run with the So ce90 A3 clone were then used to prepare Example 1 in the subject application. Thus, the data reported in Example 1 were not generated with DSM 6773, and DSM 6773 was erroneously listed in the application as the starting material.

71. To demonstrate that wild strain DSM 6773 produces epothilones C and D, the strain DSM 6773 was ordered in 2005 and the production process as reported in the subject application was followed to generate and isolate epothilones C and D. These experiments are described at page 5 of the accompanying document, titled "Reply to the Opposition Statement against EP-B-1186606." (For completeness of the record, "epothilones A and B" at page 6, line 6 therein should read --epothilones C and D--.) The experiment produced 1.4 mg epothilone C and 0.5 mg epothilone D as reported in the attached Reply.

I hereby declare that all statements made herein of my own knowledge are true and that all statements made on information and belief are believed to be true; and further that these statements were made with the knowledge that willful false statements and the like so made are punishable by fine or imprisonment, or both, under Section 1001 of Title 18 of the United States Code and that such willful false statements may jeopardize the validity of the application or any patent issued thereon.

  
Gerhard Hoefle

Date: June 10, 2008

Evidence for the discovery of epothilones C and D from *Sorangium cellulosum* at  
the GBF in early [REDACTED]

No.	Date of document	Type of document	Operation, Results
1.	[REDACTED]	Screening record	<i>Sorangium cellulosum</i> , strain So cell198 produces two new epothilones
2.	[REDACTED]	Screening record	HPLC sample preparation
3.	[REDACTED]	HPLC/DAD chromatogram	Two minor components identified as new lipophilic epothilones
4.	[REDACTED]	Sampling record	Sample given to H. Steinmetz for HPLC/MS analysis
5.	[REDACTED]	HPLC/MS chromatogram	The new epothilones contain one oxygen less than epothilones A and B
6.	[REDACTED]	Record of project meeting No. 230	Two new epothilone homologues with one oxygen less than epothilones A and B
7.	[REDACTED]	Isolation record	TLC and Sephadex LH 20 chromatography, enriched fraction
8.	[REDACTED]	LH 20 chromatogram	Separation of crude extract
9.	[REDACTED]	Analytical HPLC	The two new epothilones localised (X)
10.	[REDACTED]	Separation record	Fractions RP1 (0.7 mg) and RP2 (1.0 mg) isolated
11.	[REDACTED]	RP18 chromatogram	Separate peaks for RP1 and RP2
12.	[REDACTED]	NMR Spectra	<sup>1</sup> H and COSY spectra prove that RP2 is an epothilone with a methyl substituted 12,13 double bond later named epothilone D
13.	[REDACTED]	TLC	Purification of fraction RP1 to give RP1/DC
14.	[REDACTED]	Flow diagram	Origin of the two new epothilones (C and D)
15.	[REDACTED]	NMR spectra	<sup>1</sup> H spectrum proves that RP1/DC is an epothilone with a 12,13 double bond later named epothilone C
16.	[REDACTED]	Record of project meeting No. 231	The two new epothilones show reduced cytotoxicity and tubulin activity

# Exhibit 3-1

7


EXHIBIT A

DECLARATION OF TRANSLATOR

I, Ronald Richards, of the Technical Translation Agency, 2136  
Laa/Thaya, Austria,

do hereby avow and declare that I am conversant with the English and  
German languages and am a competent translator of German into English. I  
declare further that to the best of my knowledge and belief the following is a true  
and correct translation prepared and reviewed by me of the document in the  
German language attached hereto.

I hereby declare that all statements made herein of my own knowledge are  
true and that all statements made on information and belief are believed to be true;  
and further that these statements were made with the knowledge that willful false  
statements and the like so made are punishable by fine or imprisonment, or both,  
under Section 1001 of Title 18 of the United States Code and that such willful  
false statements may jeopardize the validity of U.S. Patent Application Serial No.  
09/313,524 or any patent issued thereon.

  
\_\_\_\_\_

Date: 5. August 2003

03/4/169

1

Table 1

Strain: 1198 – 45/30	Aim: Screening	Growth	Liquid culture	Growth
Received on: [REDACTED]	Inoculated on heso: [REDACTED]	good	Inoculated in 50 ml heso: [REDACTED]	good
Preserved:	Inoculated on propion: [REDACTED]		Inoculated in 50 ml propion (German assumed to be misspelled): [REDACTED]	
			Inoculated in 250 ml S: [REDACTED]	moderate
			Inoculated in 3 x 250 ml am (formate?): [REDACTED]	
			Medium 360: 1st flask very dark, Medium propion (German assumed to be misspelled): [REDACTED]	Harvest: 30.4 2nd flask good
			Medium heso:	
Screening on:				
	Inhibition dil. step:			
E.coli			HPLC	Substances
Micrococcus			Method screen 1	
Staph				
Nocardia				
Mucor				
Hansenula			Method screen 2	
Candida			13.282 HP	Epothilone A
Schizo			14.584 NP	Epothilone B
Rhodotorula			20.204 + 21.458 NP	Unknown epo
			Special comments:	

No b. scr.  
sample was  
prepared by  
Herr  
Steinmetz!

Tabella 1

Stamm: 1198 - 45/30	Ziel: <i>Surelevierung</i>	Wachstum	Flüssigkultur	Wachstum
Erhalten am: [redacted]	Überimpft auf Heso: [redacted]	gut	Überimpft in 50 ml Heso: [redacted]	gut
Konserviert:	Überimpft auf Propion: [redacted]		" in 250 ml S: [redacted]	wäßig
			Überimpft in 3 x 250 ml am: [redacted]	
			Medium 360: <i>A. kelben</i> sehr dunkel	Ernte: 30.4.
			Medium Probion:	2k. gut
			Medium Heso:	
Screening am: [redacted]				
E. coli	Hemmung Verd. Stufe:		HPLC	Substanzen
Micrococcus			Methode Screen. 1	
Staph			13.282 HP	Epothilon A
Nocardia			14.584 NP	Epothilon B
			20.204 + 21.458 NP	Epo. unbekannt
Mucor				
Hansenula			Methode Screen. 2	
Candida				
Schizo				
Rhodotorula				
			Besonderheiten:	

1

# Exhibit 3-2



2.

EXHIBIT

DECLARATION OF TRANSLATOR

I, Ronald Richards, of the Technical Translation Agency, 2136  
Laa/Thaya, Austria,

do hereby avow and declare that I am conversant with the English and  
German languages and am a competent translator of German into English. I  
declare further that to the best of my knowledge and belief the following is a true  
and correct translation prepared and reviewed by me of the document in the  
German language attached hereto.

I hereby declare that all statements made herein of my own knowledge are  
true and that all statements made on information and belief are believed to be true;  
and further that these statements were made with the knowledge that willful false  
statements and the like so made are punishable by fine or imprisonment, or both,  
under Section 1001 of Title 18 of the United States Code and that such willful  
false statements may jeopardize the validity of U.S. Patent Application Serial No.  
09/313,524 or any patent issued thereon.

  
\_\_\_\_\_

Date: 5. August 2003

[page of handwritten notes]

Monday

- Sample from F100 fetched, microscope inspection → OK  
Can be inoculated
- Protocol for fermentation comp. 55 taken to Herr Schüler,  
Flask order to Frau Heiber.
- Medium + XAD from Ratjadon – fermentation taken to Herr Ebert
- Strain cultures inoculated + for fermentation
- Antifoam flask + alkali flask autoclaved, antifoam sterilised in drying cupboard →  
transferred/decanted (sterile)
- Sample from F900 checked → OK
- 1 litre P medium boiled ; E medium → thorax treatment → autoclaved.

2

100% gluc - 0% H<sub>2</sub>O  
    \        /  
      35% gluc  
      /        \  
35 ml gluc      65 ml H<sub>2</sub>O      =      100 ml  
  
3 times amount added      —      300 µl / 50 ml P medium

=====

Tue

- Screening - strains : Harvest!  
1198 – 45/30 , 1230 , 1233 , 1235 :
  - 7.30 - sample removed , preparation as for HPLC →  
take methanol flask to Herr Steinmetz →  
carries out analysis
  - check other fermentation protocols
  - HPLC of Soce90 clone (medium with skimmed milk from KS)
  - 2 new screening - strains prepared in 10 ml H medium : Soce1266 + 1257
  - For fermentor : Soce360A1 further inoculated / 6 flasks available  
                              : Soce 1149      "      "      / 3 flasks available
- =====

Thu

- Fermentor sample : HPLC – preparation (Herr Steinmetz)
- Screening – strains -- analysis (1198 – 45/30 ; 1227; 1230; 1233; 1235; 1251)
- Protocol 44 – first evaluation
- Fermentor protocols : sterility check , further inoculation
- Mon, Soce 1149 3 flasks, further inoculation → 6 flasks
- Fri, Soce 360A1 1,5 litres in inoculation flask → inoculation deadline 10.15



# Exhibit 3-3

Epo unknown  
Epo unbekannt

Injection Date : 05-13-2000 13:38:30  
Sample Name : 1198-45/30  
Acq. Operator : Gerth


Seq. Line : 1  
Vial : 0  
Inj : 1  
Inj Volume : 10 µl

Sequence File : C:\HPCHEM\1\SEQUENCE\DEF\_LC.S  
Method : C:\HPCHEM\1\METHODS\SCREEN1.M  
Last changed : 05-13-2000 12:36:13 by Gerth  
Screening 1 Methode

Instrument Conditions: At Start At Stop  
Temperature: 39.8 39.8 °C  
Pressure: 190.0 213.0 bar  
Flow: 0.500 0.500 ml/min

Solvent Description :  
LC, Solvent A : Wasser  
LC, Solvent B : Methanol

DAD1 A, Sig=300,200 Ref=598,4 of 000-0101.D



### Area Percent Report

Sorted by Signal

Multiplier : 1.000000  
Dilution : 1.000000


Signal 1: DAD1 A, Sig=300,200 Ref=598,4

Peak #	RT [min]	Type	Width [min]	Area [mAU*s]	Height [mAU]	Area %
1	2.544	BB	0.208	580.48407	37.79762	2.0091
2	3.447	BB	0.380	6948.34668	277.56702	24.0491
3	4.922	BB	0.443	3222.51660	111.79380	11.1535

RT [min]	Type	Width [min]	Area [mAU*s]	Height [mAU]	Area %
8.715	BB	0.689	1983.70227	39.05926	6.8658
10.388	BB	0.513	3327.00781	91.31287	11.5152
11.921	BB	0.400	344.29031	11.95683	1.1916
13.282	BB	0.466	6934.54688	215.95573	24.0013
14.584	BB	0.355	3090.42578	125.53922	10.6963
17.392	BB	0.289	759.20789	37.43296	2.6277
18.577	BB	0.231	320.88089	20.30108	1.1106
20.207	BB	0.265	701.43530	37.26771	2.4278
21.458	BB	0.228	453.19839	28.68863	1.5686
42.724	BB	0.213	226.30806	17.73971	0.7833

Totals : 28892.35156 1052.41248


Peak :1 at 2.544 min Name : ?



-> The purity factor exceeds the thres

Purity factor : 796.557 (100%  
of spectra)  
Threshold : 990 (Set by user)  
Reference : Peak Apex  
(integrated) (2.545167)  
Spectra : 2 (Selection  
automatic, 3)


Peak :2 at 3.447 min Name : ?



-> The purity factor exceeds the thres

Purity factor : 842.491 (100%  
of spectra)  
Threshold : 990 (Set by user)  
Reference : Peak Apex  
(integrated) (3.444667)  
Spectra : 2 (Selection  
automatic, 3)


Peak :3 at 4.922 min Name : ?



-> Not enough data for purity calculat

Purity factor : Not available  
Threshold :  
Reference : Peak Apex  
(integrated) (4.925)  
Spectra : 1 (Selection  
automatic, 3)


Peak :4 at 8.715 min Name : ?



-> Not enough data for purity calculat

Purity factor : Not available  
Threshold :  
Reference : Peak Apex  
(integrated) (8.719333)  
Spectra : 1 (Selection  
automatic, 3)

Peak :5 at 10.388 min Name : ?




-> Not enough data for purity calculat

Purity factor : Not available  
Threshold :  
Reference : Peak Apex  
(integrated) (10.38833)  
Spectra : 1 (Selection  
automatic, 3)




Peak :6 at 11.921 min Name : ?



-> Not enough data for purity calculation

Purity factor : Not available  
 Threshold :  
 Reference : Peak Apex  
 (integrated) (11.91916)  
 Spectra : 1 (Selection  
 automatic, 3)


Peak :7 at 13.282 min Name : ?



-> The purity factor is within the threshold

Purity factor : 999.671 (100%  
 of spectra)  
 Threshold : 990 (Set by user)  
 Reference : Peak Apex  
 (integrated) (13.28033)  
 Spectra : 2 (Selection  
 automatic, 3)  
 Warning : Spectral  
 absorbances > 1000 mAU


Peak :8 at 14.584 min Name : ?



-> The purity factor is within the threshold

Purity factor : 999.896 (100%  
 of spectra)  
 Threshold : 990 (Set by user)  
 Reference : Peak Apex  
 (integrated) (14.5725)  
 Spectra : 2 (Selection  
 automatic, 3)


Peak :9 at 17.392 min Name : ?



-> The purity factor is within the threshold

Purity factor : 997.584 (100%  
 of spectra)  
 Threshold : 990 (Set by user)  
 Reference : Peak Apex  
 (integrated) (17.38966)  
 Spectra : 2 (Selection  
 automatic, 3)


Peak :10 at 18.577 min Name : ?



-> The purity factor exceeds the threshold

Purity factor : 922.261 (100%  
 of spectra)  
 Threshold : 990 (Set by user)  
 Reference : Peak Apex  
 (integrated) (18.57533)  
 Spectra : 2 (Selection  
 automatic, 3)


Peak :11 at 20.207 min Name : ?



-&gt; The purity factor is within the thr

Purity factor : 994.925 (100%  
of spectra)  
Threshold : 990 (Set by user)  
Reference : Peak Apex  
(integrated) (20.2055)  
Spectra : 2 (Selection  
automatic, 3)


Peak :12 at 21.458 min Name : ?



-&gt; Not enough data for purity calculat

Purity factor : Not available  
Threshold :  
Reference : Peak Apex  
(integrated) (21.45866)  
Spectra : 1 (Selection  
automatic, 3)

Peak :13 at 42.724 min Name : ?



-&gt; Not enough data for purity calculat

Purity factor : Not available  
Threshold :  
Reference : Peak Apex  
(integrated) (42.72333)  
Spectra : 1 (Selection  
automatic, 3)

\*\*\* End of Report \*\*\*

# Exhibit 3-4


EXHIBIT

DECLARATION OF TRANSLATOR

I, Ronald Richards, of the Technical Translation Agency, 2136  
Laa/Thaya, Austria,

do hereby avow and declare that I am conversant with the English and  
German languages and am a competent translator of German into English. I  
declare further that to the best of my knowledge and belief the following is a true  
and correct translation prepared and reviewed by me of the document in the  
German language attached hereto.

I hereby declare that all statements made herein of my own knowledge are  
true and that all statements made on information and belief are believed to be true;  
and further that these statements were made with the knowledge that willful false  
statements and the like so made are punishable by fine or imprisonment, or both,  
under Section 1001 of Title 18 of the United States Code and that such willful  
false statements may jeopardize the validity of U.S. Patent Application Serial No.  
09/313,524 or any patent issued thereon.



Date: 9. August 2003

[page of handwritten notes]

4

- Soce 660 Epo

H (only survivors inoculated) very good

(illegible comment in margin) H very good

H very good

HPLC, Tues. 21.5

Tues. [REDACTED]

- HPLC of different samples or cultures
- New : prepare in liquid culture + plate :  
Soce 1241 sci.  
1198 - 45/30 several plates  
1199  
471 Epo  
523 Epo  
613 Epo

- Check Soce's on plate - mould?
- Plates from protocol 44 cleared away into cool room
- 20 ml Soce 1198 - 45/30 given to Herr Steinmetz in HPLC tube for analysis

PROTOCOL 45 STILL HAS TO BE CARRIED OUT USING SOCE 1198 - 45/30 !

ALSO PROTOCOL 37 ! ALSO PROTOCOL 44 !

- Protocol 37 :

Boil 1.5 litres E medium

E medium for 1500 ml :

0.4%	skimmed milk KS	6g
0.2%	yeast extract	3g
1%	starch	15g
0.1%	CaCl <sub>2</sub>	1.5g
0.1%	MgSO <sub>4</sub>	1.5g
50mM	hepes	17.85g
8 mg/l	Fe EDTA	12mg

After autoclaving in the thorax the runny skimmed milk homogenises.

PH 7.4

Divide up into 30 x 50 ml portions in 250ml flasks + 1 ml XAD per flask each time.

Inoculate Soce 90 clone , Soce 950 Epo + Soce 660 Epo  
25 ml of each culture are required !

Addition of malonic acid diamide (malonamide) + succinate

• Socé 660 Epo  
 4 (nur Überstand überimpft) sehr gut  
 optisch 4 sehr gut  
 11 sehr gut

4

• HPLC, Di 21.5.

• HPLC von verschiedenen Proben bzw. Kulturen

• neu aussetzen in Flüssigkultur + Platte

Socé 11241 Scr.  
 1198-45/30 mehrere Platten  
 1199  
 471 Epo  
 523 Epo  
 613 Epo

• Kontrolle der Socé's auf Platte — Pilze?

• Platten vom Protokoll 44 in den Kühlraum geräumt

• Dpl in HPLC-Pörlchen gegeben von Socé 1198-45/30 für H-Steuerietz zur Analyse

PROTOKOLL 45 MUSS NOCH MIT SOCE 1198-45/30 DURCHGE-  
 FÜHRT WERDEN!  
 PROTOKOLL 37 AUCH! PROTOKOLL 44 AUCH!

• Protokoll 37

1.5 l E-Med. kochen

E-Medium für 1500 ml

0.4 l	Magermilch KS	6 g	} PH 7.4
0.2 l	Yeast extract	3 g	
1 l	Starke	15 g	
0.1 l	CaCl <sub>2</sub>	1.5 g	
0.1 l	MgSO <sub>4</sub>	1.5 g	
50 ml	Hepes	17.85 g	
8 mg/l	Fe-EDTA	12 mg	

nach dem Autoklavieren  
 im Thonrx die gesamte  
 Magermilch homogenisieren  
 auf 3 250 ml Kolben  
 a 50 ml aufteilen +  
 jeweils 1 ml XAD/Kolbe

Socé 90klou, Socé 950 Epo + Socé 660 Epo animpfen.  
 Von jeder Kultur werden 25 ml benötigt.

(Malonat)  
 Zugabe von Malonsäurediamid + Succinat (Zerstickensäure)

# Exhibit 3-5

Epo new-  
Epo neu  
Epo new  
Epo neu



test layout (original-file: spectra.rar)

Data File name: C:\HPCHEM\1\DATA\MITTWOCH\HS000000->

Method name: C:\HPCHEM\1\METHODS\ISO1.M

Sample Name: So 1198/ 2. pos

Sample Info: HPLC\_MS ->

Injection Time: 10:39:32 AM

Sequence Name:

Report Style: screen1


data acquired by: Antje

vial: 7


on: [REDACTED]

10:39:32 AM


5




13.17 3093.6




13.94 11739.9



14.43 4049.8




15.61 9168.1




Data File name: C:\HPCHEM\1\DATA\MIL11WOCH\H500000-7

Method name: C:\HPCHEM\1\METHODS\ISO1.M


16.18 19416.5




17.11 2073.1




18.24 3374.3




18.59 3811.6




18.93 7105.5




19.24 10912.9



21.38 7113.0



23.48 47866.2




E10A

Data File name: C:\HPCHEM\1\DATA\MITTWOCH\HS000000->


Method name: C:\HPCHEM\1\METHODS\ISO1.M

24.29 21613.5




Epo 3

25.37 2684.5




29.25 4218.0




Epo 100  
(M7 = -16)

30.57 2971.2




Epo 100  
(M7 = -16)


33.39 2040.5




36.36 10485.2




39.41 5935.8




Plot of TIC from So 1198/2. pos




Plot of "Device A" from So 1198/2. pos




Plot of Spectrum from 23.38 min (8 scans) from So 1198/2. pos,




Plot of Spectrum from 24.18 min (6 scans) from So 1198/2. pos,




Plot of Spectrum from 29.15 min (9 scans) from So 1198/2. pos,



Plot of Spectrum from 30.49 min (9 scans) from So 1198/2. pos,



Plot of Spectrum from 36.19 min (12 scans) from So 1198/2. pos, subtracted



# Exhibit 3-6

6


EXHIBIT

DECLARATION OF TRANSLATOR

I, Ronald Richards, of the Technical Translation Agency, 2136  
Laa/Thaya, Austria,

do hereby avow and declare that I am conversant with the English and  
German languages and am a competent translator of German into English. I  
declare further that to the best of my knowledge and belief the following is a true  
and correct translation prepared and reviewed by me of the document in the  
German language attached hereto.

I hereby declare that all statements made herein of my own knowledge are  
true and that all statements made on information and belief are believed to be true;  
and further that these statements were made with the knowledge that willful false  
statements and the like so made are punishable by fine or imprisonment, or both,  
under Section 1001 of Title 18 of the United States Code and that such willful  
false statements may jeopardize the validity of U.S. Patent Application Serial No.  
09/313,524 or any patent issued thereon.

  
\_\_\_\_\_

Date: 5. August 2003

Confidential

6

Minutes no. 230 of meeting held at 09.00 on [REDACTED]

Present: Frau Kunze, Herr Augustiniak, Forche, Gerth, Höfle, Irschik, Jansen (not full time due to EDP Committee), Reichenbach, Sasse, Steinmetz, Washausen.

**Epothilone** (Reichenbach): During a visit to **Asta Medica**, first test results were presented. They showed that epothilone A and B had similar or better action than Taxol on four selected tumor cell lines. Epothilone B was up to ten times more potent. The LD<sub>50</sub> in mice was determined to be ca. 100 mg/kg. In the case of a rapidly growing leukemia, a therapeutic dose of 20 - 30 mg/kg was able to achieve approximately 40% prolongation of life, a value similar to that of the standard cyclophosphamide. There is great interest in testing epothilone further and developing it as a cytostatic agent. Preclinical trials would require 10 - 20 grams of substance, and ca. 100 g - 1 kg would be required up to clinical phase II. Authorisation could be applied for in three years in most favourable circumstances.

Testing of epothilones at **Boehringer Mannheim** has been delayed due to their tumor research being transferred to Italy. Nonetheless 100 mg of Epo A has been supplied, with which an in-vivo study in Mannheim will be run in parallel. **Bayer AG** has also expressed interest in epothilone and received 10 mg Epo A for initial trials. They are not primarily interested in using the natural product, but in derivatives for drug targeting. **Behringwerke** became aware of epothilone from a newspaper article and will be getting in touch with us shortly. Following a suggestion from Prof. Flohé, the US company **Sugen** has requested and already received 2 mg epothilone. **Bristol-Meyers-Squibb** has made an order for 100 mg epothilone. The state of testing at **Upjohn-Pharmacia** is not known. According to information from various sources, **Merck**, **Sharp and Dohme** stopped work on epothilone some time ago. **Ciba-Geigy** continues to express interest but has not yet made any definite orders.

**Ciba test results** (Reichenbach): Testing of Condramid is complete, the results do not require further work and the substance will shortly be released. The second attempt to carry out testing of Thuggacine against *Mykobakterium tuberculosis* in England again went wrong. Ciba will provide further material from its own supplies.

**Epothilone (Steinmetz):** Stocks of A/B mixture have decreased to 0.4 g after 100 mg were again used for derivatisation. The mixture needs further cleaning before using for test purposes. About 1-2 g epothilone mixture are available as raw extract.

**Epothilone (Gerth):** A 700 litre F-24 fermentor was run with So ce90 under sterile conditions but hardly produced anything. It was channelled using 0.5 mg/l epothilone and 6 - 7 mg Spirangien. The preculture for this fermentor had produced ca. 20 mg/l epothilone according to expectations, however. Another fermentor for producing epothilone using strain So ce660 is planned to run next week. This strain only produces Epo A and Spirangiens.

We have now been successful in plating strain So ce90. As already discovered for other strains, 10% of an old autoclaved culture has to be added to the strain. The usual procedures for optimising the strain can now be implemented.

It is very important that we should try to incorporate butyrate instead of acetate and propionate in the epoxide area, following the concept of mutasynthesis. The resulting ethyl analogue of epothilone could be more biologically active and would be patentable (ask Herr Boeters).

The following strains have been identified as new epothilone producers: So ce611, So ce498, So ce931, So ce618, So ce414, So ce320 and So ce1087. All these strains are less effective producers and also form Spirangins or Icumazole. An exception is strain So ce1198, which in addition to epothilones A and B also produces a small quantity of an unknown substance with two peaks in the uV spectrum that are more lipophilic than the epothilones. According to HPLC/MS studies carried out by Herr Steinmetz, the substance is composed of homologues ( $\Delta M 14^*$ ) possessing one oxygen less than epothilone A and B. There were ca. 1-2 mg of the new epothilones present in the shaken culture. They should be isolated individually or as a mixture. According to the NMR measurement, a biological test should be carried out. While optimising the media, the strains So ce90, So ce660 and So ce950 were cultivated in 10 different media. There were significant variations in growth and production between different strains.

\* (handwritten note)  $\Delta M 14$  is typing mistake, should be  $\Delta M 16$



6

1

Vertraulich

Protokoll Nr. 230 der Besprechung vom [REDACTED], 9.00 Uhr

Teilnehmer: Frau Kunze, die Herren Augustiniak, Forche, Gerth, Höfle, Irschik, Jansen (zeitweise wegen EDV-Kommission), Reichenbach, Sasse, Steinmetz, Washausen.

**Epothilon (Reichenbach):** Bei einem Besuch bei der Asta Medica wurden erste Versuchsergebnisse vorgestellt. Danach wirken Epothilon A und B bei vier ausgewählten Tumorzelllinien ähnlich oder besser als Taxol. Epothilon B erwies sich dabei bis zu 10x aktiver. Die LD<sub>50</sub> in der Maus wurde zu ca. 100 mg/kg bestimmt. Im Fall einer schnell wachsenden Leukämie konnte mit einer therapeutischen Dosis von 20 - 30 mg/kg eine ca. 40%ige Lebensverlängerung erzielt werden, ein Wert, der dem Standard Cyclophosphamid entspricht. Es besteht großes Interesse, Epothilon exklusiv weiter zu testen und als Cytostatikum zu entwickeln. Für Vorklinische-Versuche würden 10 - 20 Gramm Substanz, bis zur klinischen Phase II, ca. 100 g - 1kg benötigt. Eine Zulassung könnte im günstigsten Fall in drei Jahren beantragt werden.

Bei Boehringer Mannheim hat sich die Testung der Epothilone verzögert, da die Tumorforschung nach Italien verlegt worden ist. Es wurden jedoch 100 mg Nachsubstanz Epo A geliefert, mit denen parallel eine in-vivo Studie in Mannheim durchgeführt wird. Die Bayer AG hat ebenfalls Interesse an Epothilon bekundet und für erste Versuche 10 mg Epo A erhalten. Dort ist man nicht an der Anwendung des Naturstoffs primär interessiert, sondern an Derivaten im Sinne von Drug-Targeting. Die Behringwerke sind durch einen Zeitungsartikel auf Epothilon aufmerksam geworden und werden demnächst mit uns Kontakt aufnehmen. Auf einen Hinweis von Prof. Flohé hat die Firma Sugan (USA) um 2 mg Epothilon A gebeten und es bereits bekommen. Von Bristol-Meyers-Squibb liegt eine Bestellung für 100 mg Epothilon vor. Der Stand der Testung bei Upjohn-Pharmacia ist nicht bekannt. Nach Informationen aus verschiedenen Quellen hat Merck, Sharp and Dohme die Bearbeitung von Epothilon bereits seit längerer Zeit aufgegeben. Ciba-Geigy ist weiterhin interessiert, allerdings ist bis jetzt keine konkrete Substanzbestellung eingegangen.

**Testergebnisse Ciba (Reichenbach):** Die Testung von Condramid ist abgeschlossen, die Ergebnisse rechtfertigen keine weitere Bearbeitung, und die Substanz wird demnächst frei-

gegeben. Auch beim zweiten Anlauf ist die Testung des Thuggacins gegen *Mykobakterium tuberculosis* in England schiefgegangen. Die Ciba wird Nachsubstanz aus dem eigenen Vorrat bereitstellen.

**Epothilon (Steinmetz):** Der Vorrat an einem A/B-Gemisch ist auf 0,4 g geschrumpft nachdem für die Derivatisierung wieder 100 mg verbraucht worden sind. Für eine Abgabe zu Testzwecken muß das Gemisch noch weiter gereinigt werden. Als Rohextrakt liegen ca. 1-2 g Epothilongemisch vor.

**Epothilon (Gerth):** Ein 700 l Fermenter F-24 mit Soce 90 ist steril gelaufen, hat jedoch kaum produziert. Er wurde bei 0.5 mg/l Epothilon und 6 - 7 mg Spirangien kanalisiert. Die Vorkultur für diesen Fermenter hatte allerdings mit ca. 20 mg/l Epothilon erwartungsgemäß produziert. Nächste Woche soll ein weiterer Fermenter zur Herstellung von Epothilon mit Stamm Soce 660 laufen. Dieser Stamm produziert nur Epo A und Spirangiene.

Es ist jetzt gelungen den Stamm Soce 90 zu plattieren. Dazu muß, wie bereits früher bei anderen Stämmen gefunden, 10% einer alten, autoklavierten Kultur des Stammes zugegeben werden. Damit können nun die üblichen Verfahren zur Stammoptimierung eingesetzt werden.

Sehr wichtig ist es, u.a. zu versuchen, im Sinne einer Mutasyntese statt Acetat und Propionat im Bereich des Epoxids Butyrat einzubauen. Das resultierende Ethylanaloge Epothilon könnte biologisch aktiver sein und wäre patentierbar (bei Herrn Boeters nachfragen).

Als neue Epothilonproduzenten wurden identifiziert: Soce 611, Soce 498, Soce 931, Soce 618, Soce 414, Soce 320 und Soce 1087. Alle diese Stämme sind schlechtere Produzenten und bilden daneben Spriangiene oder Icumazole. Eine Ausnahme bildet der Stamm Soce 1198 der neben Epothilon A und B eine unbekannte Substanz und in geringer Menge zwei lipophilere Peaks mit dem UV-Spektrum der Epothilone. Nach HPLC/MS -Untersuchungen von Herrn Steinmetz handelt es sich dabei um Homologe ( $\Delta M 14$ ) die einen Sauerstoff weniger als Epothilon A und B besitzen. In der vorliegenden Schüttelkultur liegen ca. 1-2 mg der neuen Epothilone vor. Sie sollen einzeln oder als Gemisch isoliert werden. Nach der NMR-Messung soll ein biologischer Test versucht werden.

Bei einer Medienoptimierung wurden die Stämme Soce 90, Soce 660 und Soce 950 in 10 verschiedenen Medien kultiviert. Die Variation von Wachstum und Produktion waren groß und bei den einzelnen Stämmen unterschiedlich.

*AM14 ist  
schreibfehler  
muß AM16  
sein*

# Exhibit 3-7

EXHIBIT

DECLARATION OF TRANSLATOR

I, Ronald Richards, of the Technical Translation Agency, 2136  
Laa/Thaya, Austria,

do hereby avow and declare that I am conversant with the English and  
German languages and am a competent translator of German into English. I  
declare further that to the best of my knowledge and belief the following is a true  
and correct translation prepared and reviewed by me of the document in the  
German language attached hereto.

I hereby declare that all statements made herein of my own knowledge are  
true and that all statements made on information and belief are believed to be true;  
and further that these statements were made with the knowledge that willful false  
statements and the like so made are punishable by fine or imprisonment, or both,  
under Section 1001 of Title 18 of the United States Code and that such willful  
false statements may jeopardize the validity of U.S. Patent Application Serial No.  
09/313,524 or any patent issued thereon.



---

Date: S. August 2003

Soce 1198-45/30 / Screening / [REDACTED]

MeOH-extract: Weight: 198g *mg*

7

(Illegible)

(Illegible)

(HPLC Test Result)

(HPLC Test Result)

95 CH<sub>2</sub>Cl<sub>2</sub> / 5 MeOH

95 CH<sub>2</sub>Cl<sub>2</sub> / 5 MeOH

→ LH-20 - Separation

Column = LH-20, ≅ 70 cm long, diam. 1.5 cm

Solvent = MeOH, λ = 227 nm

Flow = 1.4 ml/min, Range = 0.1 -

Paper = 2mm/min, Fractionation time = 3 min

Fractionation

LH-1 = gl. 1-9 = —

HPLC

LH-2 = gl. 10-17 = Weight 82mg

→ RP separation [REDACTED]

HPLC

LH-3 = gl. 18-23 = —

HPLC

LH-4 = gl. 24-30 = —

discarded

HPLC

LH-5 = gl. 31-41 = —

HPLC

LH-6 = gl. 42-51 = —



HPLC

SoCC 1198-45130 / Screening

Sorangium cellulosum  
SoCC 1198-45130  
Screening!  
K. Fiedler, 02.05.96

MeOH-Extrakt: Gewicht: 198 mg

7



→ LH-20-Trennung

Säule = LH-20, ≈ 70 cm lang, Ø 1,5 cm

LM = MeOH,  $\lambda = 227 \text{ nm}$

Fluss = 1,4 ml/min, Range = 0.1 -

Papier = 2 mm/min, Fraktionierungszeit = 3 min

Fractionierung

~~LH-①~~ = 6.1-9 =  
HPCC

~~LH-②~~ = 6.10-17 = Gewicht: 82 mg → RP-Trennung

~~LH-③~~ = 6.18-23 =  
HPCC

~~LH-④~~ = 6.24-30 =  
HPCC


~~LH-⑤~~ = 6.31-41 =  
HPCC

~~LH-⑥~~ = 6.42-51 =  
HPCC

werfen

# Exhibit 3-8

8



97125 1

97125 1



# Exhibit 3-9

Sample: Sol198\_LH\_2

Report Method: Spectrum\_Index\_Plot


Printed: [REDACTED]

Page: 1 of 1

Project Name: Silke\_1  
Sample Name: Sol198\_LH\_2  
Date Acquired: 12.06.96 10:57:38  
Date Processed: 12.06.96 11:59:31  
SampleWeight: 1.00000  
Dilution: 1.00000  
Channel: 996 PDA 210.0 nm  
Acq Meth Set: Sora\_MS\_210nm  
Processing Method: Epothilon\_210\_PM

Sample Type: Unknown  
Vial: 3 Inj. 1  
Volume: 3.00  
Run Time: 20.0 min  
Laufmittel: 76MeOH/24H2O, NH4Ac

9




# Exhibit 3-10

EXHIBITDECLARATION OF TRANSLATOR

I, Ronald Richards, of the Technical Translation Agency, 2136  
Laa/Thaya, Austria,

do hereby avow and declare that I am conversant with the English and  
German languages and am a competent translator of German into English. I  
declare further that to the best of my knowledge and belief the following is a true  
and correct translation prepared and reviewed by me of the document in the  
German language attached hereto.

I hereby declare that all statements made herein of my own knowledge are  
true and that all statements made on information and belief are believed to be true;  
and further that these statements were made with the knowledge that willful false  
statements and the like so made are punishable by fine or imprisonment, or both,  
under Section 1001 of Title 18 of the United States Code and that such willful  
false statements may jeopardize the validity of U.S. Patent Application Serial No.  
09/313,524 or any patent issued thereon.



---

Date: 5. August 2003

RP separation of So1198 – LH-2

So1198 LH-2 Weight: 82mg separated in 2 runs

Column = Nucleosil 100 C18, 7 $\mu$ m, 20 x 250 mm

Solvent = 73 MeOH / 27 H<sub>2</sub>O +0.05M NH<sub>4</sub>Ac

$\lambda$  = 210 nm , Range = 016 - 128

Pump = 200 , Paper = 5 mm/ min

Fractions concentrated up to H<sub>2</sub>O phase, extracted 2 x with EE, EE phase washed with H<sub>2</sub>O and dried with Na<sub>2</sub>SO<sub>4</sub>.

(Test Result)

95 CH<sub>2</sub>Cl<sub>2</sub> / 5 MeOH

Sprayed with vanillin – H<sub>2</sub>SO<sub>4</sub>

Epo A/B = 2mg/ml

Fractionation

So 1198 – RP – 1 = Weight 0.7 mg, NMR 002549, .. 660 ..

→ 0.1 mg 2nd test 150 ng/ ml

→ prep. DC 19.6.96

So 1198 – RP – 2 = Weight: 1.0 mg, NMR 002550, .. 661 ..

COSY


→ 0.1mg 2nd test 100 ng/ml

## RP-Trennung von So.1198-LH-(2)

10

So.1198-LH-(2) = Gewicht: 82mg in 2 Zylinder getrennt  
Säule = Wulfsorb 100  $\mu$ m, 7  $\mu$ m, 20 x 250mm  
LM = 73 MeOH 127 H<sub>2</sub>O 7 + 0,05M NH<sub>4</sub>Ac  
 $\lambda$  = 210nm, Range: 0,16 - 1,28  
Pumpe = 200, Papier = Summa/min

Fractionen bis zu H<sub>2</sub>O-Phase eingefügt, 2x mit  
EE extrahiert, EE-Phase mit H<sub>2</sub>O gewaschen und mit  
Na<sub>2</sub>SO<sub>4</sub> getrocknet.



angesprüht mit Vanillin-H<sub>2</sub>SO<sub>4</sub>

Epo A/B = 2mg/ml

### Fractionierung:

~~So.1198-RP-(1)~~ = Gewicht: 0,7mg, NMR 002549, Wk. 660gr  
→ 0,1mg z. Test 150mg/ml

→ präp. DC (19.6.96)

So.1198-RP-(2) = " : 1,0mg, NMR 002550, Wk. 661gr  
cosy  
→ 0,1mg z. Test 100mg/ml

# Exhibit 3-11


Separation of  
B7 (freezing point)

11

50/148 CP-2

2100m 8

73/1401  
37 H<sub>2</sub>O J 4000 NH<sub>4</sub>Ac





# Exhibit 3-12

EXHIBIT

DECLARATION OF TRANSLATOR

I, Ronald Richards, of the Technical Translation Agency, 2136  
Laa/Thaya, Austria,

do hereby avow and declare that I am conversant with the English and  
German languages and am a competent translator of German into English. I  
declare further that to the best of my knowledge and belief the following is a true  
and correct translation prepared and reviewed by me of the document in the  
German language attached hereto.

I hereby declare that all statements made herein of my own knowledge are  
true and that all statements made on information and belief are believed to be true;  
and further that these statements were made with the knowledge that willful false  
statements and the like so made are punishable by fine or imprisonment, or both,  
under Section 1001 of Title 18 of the United States Code and that such willful  
false statements may jeopardize the validity of U.S. Patent Application Serial No.  
09/313,524 or any patent issued thereon.



Date: 5. August 2003

11  
12

**NMR REQUEST**  
GBF – Dept. of Molecular structure research

Date received: XXXXXXXXXX  
Spectrum no. 002550

Substance name: So 1198 – RP-2  
Substance producer: Pohlaus  
Dept.: NC (1.1-2) tel. 343  
Nuclear species:  $^1\text{H}_1$   
Amount of substance : 1.0 mg  
Suitable solvent:  $\text{CD}_3\text{OD}$   
Return substance? Yes

**General Information**

Store sample in fridge Y

Signal expected between  
 $\delta = 0$  and 9  
Requested: only spectra Y  
plus integral Y

**Type of experiment**

$^1\text{H}_1$  Standard spectrum Y

**Plot and Data manipulation**

$\delta = 8.9$  to  $-0.1$  (0.15 ppm/cm) Y

Special requests: COSY Y

Measured on AM-300 Y

Filed under no. SIPZ 2550110/ + COSY

Einlieferungsdatum:                     

Spektren-Nr.: 002550

# NMR-ANTRAG

GBF — Abt. Molekulare Strukturforschung

Substanz-Bez.: 61198-27-2

Strukturvorschlag:                     

Summenformel:                     

Substanzhersteller: Polysar

Abteilung: NC (1.1-7) Tel.: 343

Kernart ( $^1\text{H}$ ,  $^{13}\text{C}$ ,  $^{31}\text{P}$ , andere?)                     

Substanz-Menge: 1.0 mg, Molmasse:                     

geeignetes Lösungsmittel: CD<sub>3</sub>OD weitere Messung nach Zugabe von                     

Substanz zurück: ja ☒ nein ☐

Radioaktiv ☐ Toxisch ☐

## Allgemeine Angaben

Probe lagern im Kühlschrank ☒

im Tiefkühlfach ☐

im Dunkeln ☐

Probe auf Abruf beim Hersteller ☐

Signale erwartet zwischen  $\delta =$  0 und 9

Gewünscht: nur Spektrum ☒

plus Integral ☒

Interpretation ☐

Zahl der Akkumulationen (falls > 104):                     

## Art des Experiments

☒  $^1\text{H}$  Standardspektrum ☒

Entkopplung ☐ Differenz-NOE ☐

Differenz-Entkopplung ☐

Entkoppler-Frequenz(en):                     

☒  $^{13}\text{C}$   $^1\text{H}$ -Entkopplung:

Breitband ☐ selektiv ☐

DEPT ☐ ohne ☐

## Plot und Datenmanipulation

Gauss-Multiplikation ☐

Linienausdruck ☐

☒  $^1\text{H}$

$\delta =$  8.9 bis — 0.1 (0.15 ppm/cm) ☒

11.9 bis — 0.1 (0.2 ppm/cm) ☐

Drehungen:

10 Hz/cm ☐ von  $\delta =$                       bis                     

☒  $^{13}\text{C}$  normal ( $\delta = 220$  bis 0) ☐

anderes Format:                     

Sonderwünsche: COSY ☒

$^{13}\text{C}$  —  $^1\text{H}$  Korrel. Direkt ☐ Long-range ☐

gemessen auf ☒ AM-300  
☐ ARX-400  
☐ DMX-600

(Nicht vom Antragsteller auszufüllen)

gespeichert unter Nr. 51PE250140  
in copy

Bitte um Rücksprache ☐

Kommentar:                     

(Unterschrift)


See reports; DMK memo  
 She brings 1 sample for intemp.  
 He interprets

SIPZ2550 10 1 Pohlman

So 11/98

RP-2

Aug



Current Data Parameters  
 NAME SIPZ2550  
 EXPNO 10  
 PROCNO 1


F2 - Acquisition Parameters  
 Date\_ 4.55  
 Time 4.55  
 INSTRUM spect  
 PROBHD 5 mm QNP 1H  
 PULPROG zg30  
 TO 32768  
 SOLVENT MeOH  
 NS 246  
 DS 0  
 SMH 6172.839 Hz  
 FIDRES 0.188380 Hz  
 AQ 2.6542580 sec  
 RG 812.7  
 DM 81.000 usec  
 DE 4.50 usec  
 TE 300.0 K  
 D1 1.00000000 sec  
 P1 14.50 usec  
 DE 4.50 usec  
 SFO1 300.1318534 MHz  
 NUC1 1H  
 PL1 -3.00 dB

F2 - Processing parameters  
 S1 16384  
 SF 300.1300074 MHz  
 MDH no  
 SSB 0  
 LB 0.00 Hz  
 GB 0  
 PC 1.00


1D NMR plot parameters  
 CX 30.00 cm  
 F1P 4.400 ppm  
 F1 1320.57 Hz  
 F2P -0.100 ppm  
 F2 -30.01 Hz  
 PPHCH 0.15000 dps/cm  
 HZCH 45.01950 Hz/cm


12

SIPZ2550 10 1 Pohlman



22-2  
10/2








Epo D  
NMR Dekent 2012

SIPZ2550 10 1 Pohlen




Current Data Parameters  
NAME SIPZ2550  
EXPNO 10  
PROCNO 1

F2 - Acquisition Parameters  
Date\_   
Time 13.55  
INSTRUM spect  
PROBHD 5 mm QNP 1H  
PULPROG zg30  
TD 32768  
SOLVENT MeOH  
NS 246  
DS 0  
SWH 6172.839 Hz  
FIDRES 0.188380 Hz  
AQ 2.6542580 sec  
RG 812.7  
DM 81.000 usec  
DE 4.50 usec  
TE 300.0 K  
D1 1.00000000 sec  
P1 14.50 usec  
DE 4.50 usec  
SFO1 300.1318514 MHz  
NUC1 1H  
FL1 -3.00 dB


F2 - Processing parameters  
SI 12768  
SF 300.1306076 MHz  
WDW QSINE  
SSB 2  
LB 0.00 Hz  
GB 0  
PC 1.00

1D NMR plot parameters  
CX 30.00 cm  
FLP 4.400 ppm  
F1 1320.57 Hz  
F2 -0.100 ppm  
F2 -30.01 Hz  
P0PCX 0.15000 ppm/cm  
HZCX 45.01950 Hz/cm

SIPZ2550 10 1 Pohlen



Epo



Current Data Parameters

NAME S1721550  
EXPNO 11  
PROCNO 1

F2 - Acquisition Parameters

Date\_ 18.07  
Time 16.47  
INSTRUM spect  
PROBHD 5 mm QNP 1H  
PULPROG zgpg30  
TD 2048  
SOLVENT MeOH  
NS 72  
DS 16  
SWH 2495.010 Hz  
FIDRES 1.218267 Hz  
AQ 0.4184692 sec  
RG 9195.2  
IN 200.400 usec  
DE 4.75 usec  
TE 303.0 K  
0.0090030 sec  
d0 0.0000030 sec  
d1 1.00000000 sec  
F1 14.50 usec  
P1 1000.00 usec  
D16 0.00020000 sec  
DE 4.75 usec  
SFO1 300.1311813 MHz  
H1 IN  
PC1 -3.00 dB  
IN0 0.00940080 sec

F1 - Acquisition Parameters

WDW 1  
SSB 512  
F0 300.13113 MHz  
FIDRES 4.873087 Hz  
SN 8.313 PPM

F2 - Processing parameters

SI 2048  
SF 300.1300076 MHz  
WDW SINE  
SSB 0  
LB 0.00 Hz  
GB 0  
PC 1.00

F1 - Processing parameters

SI 1024  
AQ 1.024  
SF 300.1300076 MHz  
WDW SINE  
SSB 0  
LB 0.00 Hz  
GB 0

2D NS plot parameters

CY2 20.00 cm  
CX1 20.00 cm  
F2FLO 8.067 ppm  
F2FID 2421.16 Hz  
F2FHI -0.246 ppm  
F2FLO -73.85 Hz  
F2FHI 8.067 ppm  
F2FID 2421.17 Hz  
F2FHI -73.85 Hz  
F2FLO -0.246 ppm  
F2FHI 8.067 ppm  
F2FID 2421.17 Hz  
F2FHI -73.85 Hz  
F2FLO 0.41565 ppm/cm  
F2FHI 124.75050 Hz/cm  
F2FID 0.41556 ppm/cm  
F2FHI 124.75102 Hz/cm

# Exhibit 3-13

13

EXHIBIT

DECLARATION OF TRANSLATOR

I, Ronald Richards, of the Technical Translation Agency, 2136  
Laa/Thaya, Austria,

do hereby avow and declare that I am conversant with the English and  
German languages and am a competent translator of German into English. I  
declare further that to the best of my knowledge and belief the following is a true  
and correct translation prepared and reviewed by me of the document in the  
German language attached hereto.

I hereby declare that all statements made herein of my own knowledge are  
true and that all statements made on information and belief are believed to be true;  
and further that these statements were made with the knowledge that willful false  
statements and the like so made are punishable by fine or imprisonment, or both,  
under Section 1001 of Title 18 of the United States Code and that such willful  
false statements may jeopardize the validity of U.S. Patent Application Serial No.  
09/313,524 or any patent issued thereon.

  
\_\_\_\_\_

Date: \_\_\_\_\_

S. August 2003

Preparative DC of Sol1198 / RP-1

13

Sol1198/ RP-1 = Weight = 0.6 mg put in  $\text{CH}_2\text{Cl}_2$

K.G<sub>60</sub> F<sub>254nm</sub>, diam. 0.2 mm, Al foil, 7 x 7 cm

DC solvent = 95  $\text{CH}_2\text{Cl}_2$  / 5 MeOH

front

Bands cut out, extracted 3 x with MeOH in centrifuge tube, oil pump, absorbed in  $\text{CH}_2\text{Cl}_2$ , filtered through cotton wool, concentrated, oil pump

(Test result)

(254nm)

Sol1198-RP-1 / DC Weight 0.4 mg, NMR 002630  
→ 0.1 mg 2nd test

(Test result)

Sprayed with  
vanillin -  $\text{H}_2\text{SO}_4$


Präparatives DC von 501198 RP-①

13

501198 RP-① = Gewicht = 0,6 mg in  $\text{CH}_2\text{Cl}_2$  aufgetragen


K660 F25mm, Ø 0,2mm, Alufolie, 7 x 7 cm

DC-LM = 95  $\text{CH}_2\text{Cl}_2$  5 MeOH



Bande herausgeschnitten, mit MeOH in Fein-  
brüfugenglas 3x extrahiert, MeOH-Extrakt ein-  
geengt, Ölperle, in  $\text{CH}_2\text{Cl}_2$  aufgenommen,  
über Watte filtriert, eingeeengt, Ölperle.

501198-RP-①/DC = Gewicht: 0,4 mg, NMR 200 2630  
→ 0,1 mg 2. Teil



angereicht mit  
Vanillin-H<sub>2</sub>O<sub>2</sub>

# Exhibit 3-14



14


EXHIBIT

DECLARATION OF TRANSLATOR

I, Ronald Richards, of the Technical Translation Agency, 2136  
Laa/Thaya, Austria,

do hereby avow and declare that I am conversant with the English and  
German languages and am a competent translator of German into English. I  
declare further that to the best of my knowledge and belief the following is a true  
and correct translation prepared and reviewed by me of the document in the  
German language attached hereto.

I hereby declare that all statements made herein of my own knowledge are  
true and that all statements made on information and belief are believed to be true;  
and further that these statements were made with the knowledge that willful false  
statements and the like so made are punishable by fine or imprisonment, or both,  
under Section 1001 of Title 18 of the United States Code and that such willful  
false statements may jeopardize the validity of U.S. Patent Application Serial No.  
09/313,524 or any patent issued thereon.

  
\_\_\_\_\_

Date: S. August 2003

Soce 1198 - 45/30 / screening / [redacted]

from 5-10 stake flasks 14

Soce 1198 - 45/30 / screening / [redacted]

MeOH extract  
198 mg

from microbiology

LH-20 separation

LH-1  
discarded

LH-2  
82 mg  
Epo A,B + Epo un

LH-3, LH-4, LH-5, LH-6  
discarded

RP separation

See # 9

See # 10

RP-1  
0.7 mg  
NMR 002549  
Epo un

RP-2  
1.0 mg  
NMR 002550, COSY  
Epo un

Prep DC

selects RP-2

EPO D

not pure so further separation #13

Rp-1 / DC  
0.4 mg  
NMR 002630  
Epo un

→ EPOC

14

Soc 1198 - 45130 Screening

MeOH-Extrakt  
198 mg

CH-20-Trennung

CH-①  
Verworfen

CH-②  
82 mg  
Epo AD-Epo-mer

CH-③  
Verworfen

CH-⑥

RP-Trennung

RP-①  
0,7 mg  
NMR 002549  
Epo-mer

präp. DC

RP-① (DC)  
0,4 mg  
NMR 002630  
Epo-mer

RP-②  
1,0 mg  
NMR 002550, cosy  
Epo-mer

# Exhibit 3-15

15


EXHIBIT

DECLARATION OF TRANSLATOR

I, Ronald Richards, of the Technical Translation Agency, 2136  
Laa/Thaya, Austria,

do hereby avow and declare that I am conversant with the English and  
German languages and am a competent translator of German into English. I  
declare further that to the best of my knowledge and belief the following is a true  
and correct translation prepared and reviewed by me of the document in the  
German language attached hereto.

I hereby declare that all statements made herein of my own knowledge are  
true and that all statements made on information and belief are believed to be true;  
and further that these statements were made with the knowledge that willful false  
statements and the like so made are punishable by fine or imprisonment, or both,  
under Section 1001 of Title 18 of the United States Code and that such willful  
false statements may jeopardize the validity of U.S. Patent Application Serial No.  
09/313,524 or any patent issued thereon.



---

Date: S. August 2003

**NMR REQUEST**  
GBF – Dept. of Molecular structure research

Date received: [REDACTED]

Spectrum no. 002630

Substance name: So 1198 – RP-1 / DC

Substance producer: Pohlaus

Dept.: NC (1.1-2) tel. 343

Nuclear species:  $^1\text{H}_1$

Amount of substance : 0.4 mg

Suitable solvent:  $\text{CD}_3\text{OD}$

Return substance? Yes

**General Information**

Store sample in fridge Y

Signal expected between

$\delta = 0$  and 9

Requested: only spectra Y

plus integral Y

**Type of experiment**

$^1\text{H}_1$  Standard spectrum Y

**Plot and Data manipulation**

$\delta = 8.9$  to  $-0.1$  (0.15 ppm/cm) Y

Filed under no. SIPR 2640\ ??

Einlieferungsdatum:                     

Spektren-Nr.:                     

**002630**

# **NMR-ANTRAG**

GBF — Abt. Molekulare Strukturforschung

Substanz-Bez.: 501198 127-① DC

Strukturvorschlag:                     

Summenformel:                     

Substanzhersteller: Powles

Abteilung: MC (1.1.2) Tel.: 343

Kernart: (<sup>1</sup>H, <sup>13</sup>C, <sup>31</sup>P, andere?)                     

Substanz-Menge: 0.4 mg, Molmasse:                     

geeignetes Lösungsmittel: CD<sub>3</sub>OD weitere Messung nach Zugabe von                     

Substanz zurück: ja ☒ nein ☐

Radioaktiv ☐ Toxisch ☐

## **Allgemeine Angaben**

Probe lagern im Kühlschrank ☒  
im Tiefkühlfach ☐  
im Dunkeln ☐  
Probe auf Abruf beim Hersteller ☐

Signale erwartet zwischen  $\delta =$  0 und 9  
Gewünscht: nur Spektrum ☒  
plus Integral ☒  
Interpretation ☐  
Zahl der Akkumulationen (falls > 104):                     

## **Art des Experiments**

☒ Standardspektrum ☒  
Entkopplung ☐ Differenz-NOE ☐  
Differenz-Entkopplung ☐  
Entkoppler-Frequenz(en):                     

☒ <sup>13</sup>C <sup>1</sup>H-Entkopplung:  
Breitband ☐ selektiv ☐  
DEPT ☐ ohne ☐

## **Plot und Datenmanipulation**

Gauss-Multiplikation ☐

Linienausdruck ☐

☒ <sup>1</sup>H  
 $\delta =$  8.9 bis — 0.1 (0.15 ppm/cm) ☒  
11.9 bis — 0.1 (0.2 ppm/cm) ☐

Drehungen:  
10 Hz/cm ☐ von  $\delta =$                       bis                     

☐ <sup>13</sup>C normal ( $\delta =$  220 bis 0) ☐

anderes Format:                     

Sonderwünsche: COSY ☐

<sup>13</sup>C — <sup>1</sup>H Korrel. Direkt ☐ Long-range ☐

(Nicht vom Antragsteller auszufüllen)

gemessen auf ☐ AM-300  
☐ ARX-400  
☐ DMX-600

gespeichert unter Nr. SIRBAS/101  
                      
                    

Bitte um Rücksprache ☐

Kommentar:                     

(Unterschrift)

So 1198.1  
RP-1 DC

0.4mg

SIPR2630 10 1





15

Current Data Parameters  
NAME SIPR2630  
EXPNO 10  
PROCNO 1

F2 - Acquisition Parameters  
Date\_ 11.39  
Time 11.39  
INSTRUM 300K400  
PROBHD 5 mm QNP 1H  
PULPROG zg30  
TD 32768  
SOLVENT MeOH  
NS 280  
DS 2  
SMH 8333.333 Hz  
FIDRES 0.254313 Hz  
AQ 1.9661300 sec  
RG 4096  
DH 60.000 usec  
DE 85.71 usec  
TE 300.0 K  
D1 1.0000000 sec  
P1 14.00 usec  
DE 85.71 usec  
SFO1 400.1324710 MHz  
NUCLEUS 1H

F2 - Processing parameters  
SI 16384  
SF 400.1259912 MHz  
WDW no  
SSB 0  
LB 0.00 Hz  
GB 0  
PC 1.40

10 NMR plot parameters  
CX 30.00 cm  
F1P 4.400 ppm  
F1 1760.57 Hz  
F2P -40.100 ppm  
F2 -40.01 Hz  
PPMCH 0.15000 ppm/cm  
HZCH 60.01950 Hz/cm





Epo C  
NMR Patient -  
Anispride

20.6.96


SIPR2630 10 1

Current Data Parameters  
NAME SIPR2630  
EXPO 10  
PROCNO 1

F2 - Acquisition Parameters  
Date\_   
Time 11.39  
INSTRUM arx400  
PROBHD 5 mm QNP 1H  
PULPROG zg30  
TD 32768  
SOLVENT MeOH  
NS 280  
DS 2  
SWH 8331.333 Hz  
FIDRES 0.254313 Hz  
AQ 1.9661300 sec  
RG 4096  
DE 60.000 usec  
TE 300.0 K  
D1 1.60000000 sec  
P1 14.00 usec  
DE 85.71 usec  
SFO1 400.1324710 MHz  
NUCLEUS 1H

F2 - Processing parameters  
SI 32768  
SF 400.1299912 MHz  
WDW QSIING  
SSB 2  
LB 0.00 Hz  
GB 0  
PC 1.40


1D NMR plot parameters  
CK 30.00 cm  
F1P 4.400 ppa  
F1 1760.57 Hz  
F2P -0.100 ppa  
F2 -40.01 Hz  
FREQ 0.15000 ppa/cm  
HZCH 60.01950 Hz/cm



1H NMR spectrum of compound 11, showing peaks at 0.214, 0.513, 0.810, and 0.413 ppm. The spectrum is labeled '11, 11, 4' and 'ppm'.

dd, vv, y

SIPR2630 10 1



# Exhibit 3-16

16

EXHIBIT

DECLARATION OF TRANSLATOR

I, Ronald Richards, of the Technical Translation Agency, 2136  
Laa/Thaya, Austria,

do hereby avow and declare that I am conversant with the English and  
German languages and am a competent translator of German into English. I  
declare further that to the best of my knowledge and belief the following is a true  
and correct translation prepared and reviewed by me of the document in the  
German language attached hereto.

I hereby declare that all statements made herein of my own knowledge are  
true and that all statements made on information and belief are believed to be true;  
and further that these statements were made with the knowledge that willful false  
statements and the like so made are punishable by fine or imprisonment, or both,  
under Section 1001 of Title 18 of the United States Code and that such willful  
false statements may jeopardize the validity of U.S. Patent Application Serial No.  
09/313,524 or any patent issued thereon.

  
\_\_\_\_\_

Date: \_\_\_\_\_

S. August 2003

Confidential

16

Minutes no. 231 of meeting held at 09.00 on [REDACTED]

Present: Frau Herrmann, Herr Augustiniak, Forche, Gerth, Höfle, Irschik, Jansen, Reichenbach, Sasse, Steinmetz, Washausen

Information:

- Two written offers of contract for epothilone have now been received, three more are expected in the near future.
- ASTRA has expressed an interest in testing etnangien in its polymerase tests. A test sample should be sent after concluding a confidentiality agreement.
- The NBI department plans to make TA culture extracts from myxobacteria available on payment to each of several firms for them to screen. GBF retains the rights to the strains and will have an appropriate share in any success.
- Rhone-Poulenc has applied for two patents covering a substance from actinomycete A 9738 (CBS 162.94), which is identical to Cittilin from Mx x48. The compound is described as a neurotensin antagonist. GBF's patent law firm is being asked to check whether these patents can't be challenged on the basis of our two publications and if necessary overturned.
- Ciba-Geigy Pharma in a letter dated [REDACTED] has released the following substances: eliamid, etnangien, argyris A and B. A verbal indication suggested that chondramid might be released, since the substance has no in vivo effect.
- Myxothiazol must be given further medium-term fermentation (Kunze).

**Epothilone** (Gerth, Sasse, Steinmetz): In addition to the 15 known producers we already have, 9 were recently added from Herr Irschik's screening; eight of them formed Spirangien at the same time, one formed icumazol, one only formed epothilone A; productivity was not significantly high in any of the cases. A suitable production strain has to be selected from these known 24 producers. The strains So ce90 (the original producer), So ce660 (only forms

epothilone A), So ce950 (only forms icumazol), So cel198 (free of extra substances) as well as So cel275 and 1294 (both originate from the same sample, grow better, form no known extra substances) are at present being investigated in more detail (adaptation to homogenous growth, plating, clone selection). Tests on medium optimisation indicated that the different strains react in different ways. The addition of propionate with So ce90 caused increased formation of epothilone B; for the other strains this did not occur, synthesis in part even being totally blocked despite good growth (So cel198, So e1275). The addition of formate to So ce90 caused increased formation of epothilone B and a reduction in epothilone A as well as increased synthesis overall; succinate on the other hand had no effect. So cel198 and So ce1275 formed no epothilone at all with formate. The type of starch added also has dramatic effect on epothilone synthesis: for So ce90 the best results are achieved with Cerestar (100 %); the results with wheatmeal or ryemeal were considerably worse; using Ciba starch 37 % of the Cerestar yield was achieved, soluble starch achieved 74 %; with soluble starch the ratio of A to B shifted from 1.26 (Cerestar) to 0.83. For skimmed milk and yeast extract a quality comparison still needs to be carried out. Using complex substrates such as banana, plum, or mushroom flour resulted in good growth, but epothilone production was poor or even completely suppressed. So ce477 grew well in full-fat soya flour, So ce90 only grew in low-fat soya flour. While So ce90 did not produce anything on agar plates, So cel198 still appeared to form epothilone on certain types of agar, which would make strain selection much easier. None of the strains have grown homogeneously so far. Plating is possible with So cel148: 50 clones have recently been isolated, of which 10 produce epothilone and 40 do not. The formation of Spirangien can easily be detected during cloning, so Spirangien producers can quickly be eliminated.

**Fermentations of 9.8.-18.8.:** F25 (10 l), starter culture: good growth, clean; transferred to F26 (100 l): good growth, clean; F27 (1000 l) had meanwhile been prepared: it frothed over and lost 450 l medium; the reactor was filled again and autoclaved, but again lost another 80 l; it nonetheless remained sterile and was then inoculated from F27. An aliquot of 10 l was at the same time taken from F26 and inoculated into F28 (830 l): F27 and F28 were both infected with a bacillus after 1 day; it was discovered that the inoculation tube used for both inoculations had a hole. F29 (3000 l) was still planned however: after autoclaving the skimmed milk medium the reactor was unsterile after 1 day; it was then autoclaved again and was then still sterile after 4 days; the rest of the medium was then added: 2 days later the reactor was again unsterile and was autoclaved again, it was again unsterile shortly afterwards and was discarded. Since the Biotechnikum is totally closed from Week 37 to 40 due to the ITP, the next series of fermentations will only be possible from 20.9 to 10.10; it might be possible to

add a second cascade afterwards; it is planned to run with a total of 6340 l and 5400 l production volume. From Week 43 to 46 the brine plant is being repaired but this should not affect epothilone production. Due to frothing over of reactors and infections at early stages, 50 l of expensive XAD were lost.

At present there is about 600 mg epothilone A and 400 mg epothilone B available in very impure extracts and the material is being purified at great expense. The substance is urgently needed for test samples.

100 mg epothilone A and 100 mg epothilone B have recently been sold to Bristol-Myers, 150 mg A and 150 mg B to Boehringer.

A batch of about 1.5 g epothilone got lost during recovery. Following preparative HPLC the substance was still all right, it was rotated and left overnight; It was then chlorinated by adding HCl and was inactive.

The strains So cel198, So cel275 and So cel294 form two new epothilones as well as epothilone, but with the epoxide missing. They had considerably reduced action, but were not inactive: The  $IC_{50}$  for L929 cells was 150 ng/ml for RP1 (from So cel198), and 100 ng/ml for RP2. Noticeable effect on Tubulin could be detected in cell cultures. Perhaps patenting could be possible?

Epothilone had lost all activity in mouse serum after 4 h at 37°, and showed a similar result in rat serum; however, the substance was not inactivated in serum of humans, cattle, rabbits, goat or sheep (only serums that were not heat inactivated were used). In human serum the substance was completely stable for 2 days at 37° (HPLC analysis). Lyophilized mouse serums were inactive and hamster serum slightly active. Pig liver esterase opens the lactone ring.

**Ambruticin** (Gerth): Following feeding of  $^{14}C$  Ambruticin A, only VS3 and S could subsequently be detected; radioactivity could not be found anywhere else.

**From 2.9** Herr Dipl. Biol. Knauth will be working on the mechanism of action of ambruticin and jerangolid (NBI Dept.).



Vertraulich

16

Protokoll Nr. 231 der Besprechung am [REDACTED] 9.00 Uhr

Teilnehmer: Frau Herrmann, die Herren Augustiniak, Forche, Gerth, Höfle, Irschik, Jansen, Reichenbach, Sasse, Steinmetz, Washausen

Zur Information:

- Für Epothilon liegen inzwischen zwei schriftliche Vertragsangebote vor, drei weitere sind in naher Zukunft zu erwarten.
- Die Firma ASTRA hat Interesse angemeldet, Etnangien in ihren Polymerase-Tests zu prüfen. Nach Abschluß eines Vertraulichkeitsabkommens soll ein Prüfmuster versandt werden.
- Die Abteilung NBI plant, an mehrere Firmen gegen Bezahlung von je einer TA Kulturextrakte von Myxobakterien für deren Screening zur Verfügung zu stellen. Die GBF behält die Rechte auf die Stämme und wird am Erfolg angemessen beteiligt.
- Rhône-Poulenc hat zwei Patente für eine Substanz aus Actinomyces A 9738 (CBS 162.94) angemeldet, die mit Citalin aus Mx x48 identisch ist. Die Verbindung ist als Neurotensin-Antagonist beschrieben. Das Patentbüro der GBF wird gebeten zu prüfen, ob diese Patente nicht auf Basis unserer zwei Veröffentlichungen angegriffen werden können und diese dann ggf. auch zu kippen.
- Ciba-Geigy Pharma hat mit Schreiben vom [REDACTED] folgende Substanzen freigegeben: Eliamid, Etnangien, Argyrin A und B. Mündlich wurde eine Freigabe von Chondramid in Aussicht gestellt, da die Substanz in vivo nicht wirkt.
- Myxothiazol muß mittelfristig nachfermentiert werden (Kunze).

Epothilon (Gerth, Sasse, Steinmetz): Zu den 15 schon bekannten eigenen Produzenten kamen neuerdings 9 aus dem Screening von Herrn Irschik hinzu; von diesen bildeten 8 gleichzeitig Spirangien, einer Icumazol, einer nur Epothilon A; die Produktivität war in keinem Fall ungewöhnlich hoch. Aus den somit bekannten 24 Produzenten muß ein geeigneter

Produktionsstamm ausgewählt werden. Zur Zeit werden So ce90 (der ursprüngliche Produzent), So ce660 (bildet nur Epothilon A), So ce950 (bildet nur Icumazol), So ce1198 (frei von Begleitsubstanzen) sowie So ce1275 und 1294 (stammen beide aus derselben Bodenprobe, wachsen besser, bilden keine bekannten Begleitsubstanzen) näher charakterisiert (Adaption zu homogenem Wachstum, Plattierung, Klonselektion). Versuche zur Mediumsoptimierung zeigten, daß die einzelnen Stämme unterschiedlich reagieren. Zusatz von Propionat führt bei So ce90 zu einer verstärkten Bildung von Epothilon B; für die anderen Stämme gilt dies nicht, zum Teil wird die Synthese sogar trotz guten Wachstums total blockiert (So ce1198, So ce1275). Zusatz von Formiat zu So ce90 führt zur verstärkten Bildung von Epothilon B und einer Reduzierung von Epothilon A sowie insgesamt zu einer verstärkten Synthese; Succinat hat dagegen keinen Effekt. So ce1198 und So ce1275 bilden mit Formiat überhaupt kein Epothilon. Auch die Art der zugesetzten Stärke beeinflusst die Epothilonsynthese dramatisch: Bei So ce90 werden die besten Ergebnisse mit Cerestar erhalten (100 %); mit Weizen- oder Roggenmehl Ergebnisse erheblich schlechter; mit Ciba-Stärke werden 37 %, mit löslicher Stärke 74 % der Ausbeute mit Cerestar erreicht; mit löslicher Stärke verschiebt sich dabei das Verhältnis von A zu B von 1.26 (Cerestar) zu 0.83. Für Magermilch und Hefeextrakt muß erst noch ein Qualitätsvergleich durchgeführt werden. Komplexe Substrate wie Bananen-, Pflaumen- oder Pilzmehl erhält man gutes Wachstum, aber die Epothilon-Produktion ist schlecht oder ganz unterdrückt. So ce477 wächst gut in Sojamehl vollfett, So ce90 dagegen nur in entfettetem Sojamehl. Während So ce90 auf Agarplatten nichts produziert, scheint So ce1198 auf bestimmten Agarsorten noch Epothilon zu bilden, was die Stammselektion sehr erleichtern würde. Keiner der Stämme wächst bisher homogen. So ce1148 kann plattiert werden: Inzwischen sind 50 Klone isoliert, von denen 10 Epothilon produzieren, 40 dagegen nicht. Die Bildung von Spirangien läßt sich beim Klonieren leicht erkennen, so daß Spirangien-Produzenten schnell ausgeschieden werden können.

Fermentationen vom 9.8.-18.8.: F25 (10 l), Starterkultur: gut gewachsen, sauber; überführt in F26 (100 l): gutes Wachstum, sauber; inzwischen war F27 (1000 l) vorbereitet worden: dieser schäumte aber über und verlor 450 l Medium; der Reaktor wurde wieder aufgefüllt und autoklaviert, verlor aber anschließend nochmals 80 l; trotzdem blieb er steril und wurde dann aus F27 angeimpft. Ein Aliquot von 10 l wurde aus F26 parallel in F28 (830 l) überimpft: F27 sowie F28 waren beide nach 1 d mit einem Bacillus infiziert; wie sich herausstellte, hatte der für beide Impfvorgänge verwendete Impfschlauch ein Loch. Außerdem war noch F29 (3000 l) geplant: Nach Autoklavieren des Magermilch-Mediums war der Reaktor nach 1 d unsteril; er wurde danach nochmals autoklaviert und war anschließend für 4 d steril; danach wurde der Rest des Mediums zugesetzt: 2 d später war der Reaktor wieder unsteril und wurde erneut autoklaviert, kurz danach war er wieder unsteril und wurde verworfen. Da das Biotechnikum in

der 37. - 40. Woche durch den ITP total blockiert ist, ist die nächste Fermentationsserie erst vom 20.9.-10.10. möglich; vielleicht läßt sich danach eine zweite Kaskade anschließen; vorgesehen ist insgesamt 6340 l Arbeitsvolumen mit 5400 l Produktionsvolumen. Von der 43.-46. Woche wird die Soleanlage repariert, was aber keine Auswirkungen für die Epothilon-Produktion haben sollte. Durch das Überschäumen der Reaktoren und die Infektionen auf frühem Stadium gingen 50 l teures XAD verloren.

Derzeit liegen in +/- stark verunreinigten Extrakten rund 600 mg Epothilon A und 400 mg Epothilon B vor und werden unter großem Aufwand gereinigt. Die Substanz wird dringend für Prüfmuster benötigt.

Vor kurzem wurden an Bristol-Myers 100 mg Epothilon A und 100 mg Epothilon B verkauft, an Boehringer 150 mg A und 150 mg B.

Eine Charge von rund 1.5 g Epothilon gingen bei der Aufarbeitung verloren. Nach präparativer HPLC war die Substanz noch in Ordnung, sie wurde einrotiert und stand über Nacht. Danach war sie durch HCl Addition chloriert und unwirksam.

Die Stämme So ce1198, So ce1275 und So ce1294 bilden neben Epothilon auch zwei neue Epothilone, denen das Epoxid fehlt. Deren Wirksamkeit war stark reduziert, jedoch nicht aufgehoben: Die  $IC_{50}$  für L929-Zellen betrug für RP1 (aus So ce1198) 150 ng/ml, für RP2 100 ng/ml. In Zellkulturen war auch eine deutliche Wirkung auf Tubulin zu erkennen. Vielleicht wäre eine Patentierung möglich?

Epothilon in Serum der Maus hatte nach 4 h bei 37° alle Aktivität verloren, ebenso in Serum der Ratte; in Serum von Mensch, Rind, Kaninchen, Ziege und Schaf wurde die Substanz dagegen nicht inaktiviert (verwendet wurden ausschließlich Seren, die nicht hitzeinaktiviert waren). In Humanserum war die Substanz über 2 d bei 37° völlig stabil (HPLC-Analytik). Lyophilisierte Seren der Maus waren unwirksam, des Hamsters etwas wirksam. Schweineleberesterase öffnet den Lactonring.

Ambruticin (Gerth): Nach Verfütterung von  $^{14}C$ -Ambruticin A waren anschließend nur VS3 und S nachweisbar; nirgendwo sonst war Radioaktivität zu entdecken.

Ab 2.9. wird sich Herr Dipl.Biol. Knauth mit dem Wirkmechanismus von Ambruticin und Jerangolid beschäftigen (Abt. NBI).

- Soce M98 A5 (A37) in Flüssigkultur + auf Platte angeimpft 4-10  
 (konservieren 10x + als Stammkultur weiter laufen lassen  $\Rightarrow$   
 M98 A2 verwerfen (5x konservieren))

## Laborbuch Fischer

- F-Medium für 1L  
 Grundmedium

1,5g MgSO<sub>4</sub>  
 1g Pepton (Münch)  
 1g H<sub>2</sub>PO<sub>4</sub> } pH 7,2

- Stammesg. herstellen

- Flaschen !!!

- Medium fermenter !!! F 15/

- Screening-Stämme weiterimpfen (3x 250 ml): 1332, 1335, 1358, 1339

- E-Medium für 1,2L (Protokoll 89 + 90)

4,8g Magermilch  
 2,4g Yeast extrakt  
 12g Stärke  
 1,2g CaCl<sub>2</sub>  
 1,2g MgSO<sub>4</sub>  
 14,3g H<sub>2</sub>PO<sub>4</sub>  
 1,2ml Fe-EDTA  
 4,8g Sojamehl  
 6g Glycerin } pH 7,2 + 2ml XAD / Folben

- Protokoll 88 1. Vorbereitung P 89 + P 90 auch

- Protokoll 88, 89 + 90 animpfen

- Soce 90 A3 weitergeimpft für Fermentation (4x 250ml)

- Ernte von Protokoll 77 + DSLI-Stamm

- Stammesg. herstellen:

Protokoll 89: 10 l. ige Propionatlg. = 10g Propionat + 90 ml H<sub>2</sub>O

Protokoll 90: 10 l. ige Formiatlg. = 10g Formiat + 90 ml H<sub>2</sub>O

Mied. 7: Stammesg. 1: 20 l. KNO<sub>3</sub> = 20g KNO<sub>3</sub> + 80 ml H<sub>2</sub>O

1,25 l. K<sub>2</sub>HPO<sub>4</sub> = 1,25g K<sub>2</sub>HPO<sub>4</sub> + 8,75 ml H<sub>2</sub>O

Stammesg. 2: 10 l. CaCl<sub>2</sub> = 10g CaCl<sub>2</sub> + 89,92 ml H<sub>2</sub>O

800mg/l Fe/EDTA = 80mg Fe/EDTA

Stammesg. 3: 25 l. Glukose = 25g Glukose + 75 ml H<sub>2</sub>O

35 l. ige Glukose: 70g Glukose + 130 ml H<sub>2</sub>O

Ca<sub>2</sub>-Agar: Stammesg. 1: 15,0g KNO<sub>3</sub>, 15,0g K<sub>2</sub>HPO<sub>4</sub> + 180 ml H<sub>2</sub>O

Stammesg. 2: 20g MgSO<sub>4</sub> + 99,0 ml H<sub>2</sub>O

Stammesg. 3: 0,4g CaCl<sub>2</sub> + 0,8g FeCl<sub>3</sub> + 199,80 ml H<sub>2</sub>O

$\rightarrow$  PROPIONAT?  $\leftarrow$

- 10 Kolben 1498 A2 ernten: XAD zu Hr. Steinmetz
- Glukose 25 l + 35 l
- E-Medium für 5 l

4-11

- 20 g Magermilch
- 20 g Sojamehl
- 10 g Yeast extrakt
- 50 g Stärke
- 5 g  $CaCl_2$
- 5 g  $MgSO_4$
- 5 g H<sub>2</sub>O<sub>2</sub>
- 5 ml Fe-EDTA

pH 7.4

20 Kolben à 250 ml  
+ 10 ml XAD

- F-Medium für 2 l
- Grundmedium

- 5 g  $MgSO_4$
- 2 g Pepsin (Marcor)
- 23.8 g H<sub>2</sub>O<sub>2</sub>

pH 7.2

20 Kolben à 250 ml  
+ 2 ml XAD

→ wenn das Medium erkaltet ist, die Stammlog. + Glukose dazugeben ←

- Stammkulturen überimpfen
- 10 x 250 ml Soe 90 A3 → Fermenter
- Versuch 10 x 250 ml mit 1498 A2 angereicht
- Zugabe von Propionat + Formiat bei den Protokollen 88 + 90
- Ernte von Protokoll 86

Fr

- Stammkulturen überimpfen
- 12 von 100 Soe 90 A3 weiterimpfen
- Screening: Stämme weiterimpfen (1 x 250 ml) Soe 1241 + 1270
- Stämme Soe 1305 verworfen → wächst nicht!
- Ernte der Screening-Stämme Soe 1325 + 1328
- Protokoll 77 + 86 ansetzen + HPLC
- In Lager Medium abgeben (F-150 2x)
- Floue von Soe 90 ernten (148 Stk.)
- F-Medium mit Soe 90 A3, 1275 + 1498 A5 → Fr. 27.9
- (nach 1 Woche HPLC Überprüfung) ansetzen
- Konservierung -71°C Soe 90 A3 10x, Soe 1498 A2 10x
- 50 Kolben à 20 ml E-Medium kochen (AC)

Mu.

- 50 Kolben (100 ml) à 20 ml E-Medium + XAD (0.5 ml)
- (AC Medium)

Da

- restliche Medien für Fermentation abgeben + ins Fermenter bringen (Flaschen + XAD ebenfalls)
- F-15 (2x) ansetzen → parallel (10<sup>15</sup> Uhr)
- Ernte: Screening-Stämme Soe 1333

⇒ Biologisches Screening! ←

E-Medium für 1L

- 4g - Magermilch
- 4g - Sojamehl, entfettet
- 2g - Yeast extrakt
- 10g - Stärke
- 1g -  $\text{CaCl}_2$
- 1g -  $\text{MgSO}_4$
- 11,9g - HEPES
- 1ml - Fe-EDTA

50 Kolben à 20ml  
+ 0,5 ml XAD

pH 7,4

Ernte von Protokoll 88, 89, 90

E-Medium animpfen mit 90A3, 1198 A5 + 1275 je 5 Kolben

- Soce 1198 A5 10 x 2ml bei  $-71^\circ\text{C}$  konservieren
- E-Medium animpfen mit Soce 90A3 1275 + 1198 A5  
je 5 Kolben (Zusätze nicht vergessen beim Medium)  $\Rightarrow$   
nach 1 Woche HPLC-Überprüfung ✓ dr. g.H.
- Biologisches Screening (22 Stämme) ✓

Krank  $\rightarrow$  Stirnhöhlenvereiterung

1110, 07.10.96

- Soce 1198 A5 + 1300 10 x 2ml konserv.  $-71^\circ\text{C}$
- Ernte von Screening-Stämmen Soce 1334, 1336, 1332, 1335, 1338 + 1339
- Soce 1328 schon in 50 ml
- Screening-Stämme Soce 1241 + 1270 weiterimpfen in 3 x 250 ml
- Versuche mit 1198 A2 ein Kolben ernten + Analytik erstellen
- Screening-Stämme neu animpfen Soce 1308, 1342, 1349, 1352, 1364
- MIC-Bestimmung für Oliver Gronwald 6 Stck.
- 600ml E-Medium kochen für Protokoll 91 (Soce 90A3 aussetzen)

Konzentrationsbestimmung von Epothilon

$$\text{Fläche} \times 4 \times 0,0007072 \text{ ng}$$

$$= \frac{\text{Fläche} \times 0,0028288}{(\text{Area [mAU} \cdot \text{s]})}$$

4-13

FERMENTER: Blatt 1

Kostenstelle: 103320

Vers.Nr.: 96/145/02/02

Betreiber: K. Gule

Betreuer: M. Haezel

Organismus: Soe 90

Kulturführung: Aerob: ☐ Anaerob: ☐ Phototroph: ☐

Prozeßführung: Batch: ☐ Feed-Batch: ☐ Konti: ☐

Fermenteraufbau:

Fermenter Nr. 15 L 150 A	Verwendung: Fermentation <input checked="" type="radio"/> Vorlage <input type="radio"/>	Steriltest <input type="radio"/> für Protokoll-Nr. __/__/__
Sicherheitsmaßnahmen	Abluftfilter: Nein <input checked="" type="radio"/> Ja <input type="radio"/> Handschuhe tragen: Nein <input type="radio"/> Ja <input checked="" type="radio"/>	
Betrieb-Beginn	Datum: <del>24.08.02</del>	Uhrzeit: 8 <sup>00</sup>
Rührerart	3 x Schube	
Sondergeräte		
Pumpe für	Typ: <input checked="" type="checkbox"/>	Pumprate: <input checked="" type="checkbox"/> Durchm. Schlauch: <input checked="" type="checkbox"/>
Pumpe für	Typ: <input checked="" type="checkbox"/>	Pumprate: <input checked="" type="checkbox"/> Durchm. Schlauch: <input checked="" type="checkbox"/>
	Medien mit Uebervorteil	

Elektroden:

pH-Elektrode	Nr.: 8	Puffer 1: H <sub>2</sub> O Poti/ mV: <input checked="" type="checkbox"/>	Puffer 2: H <sub>2</sub> O Poti/ mV: <input checked="" type="checkbox"/>
pH-Elektrode	Nr.:	Puffer 1: Poti/ mV:	Puffer 2: Poti/ mV:
pO <sub>2</sub> -Elektrode	Nr.: 02013	Nr.:	

Reaktorgewicht:

Gesamtgewicht

Sollgewicht	82	[KG]
leer		[KG]
Wassermenge	14.52	[l] [KG]
Medium-Zugabe	Name: Soce-80	Herk.: Nutzer <input checked="" type="checkbox"/> SE: __/__/__ [KG]
	XAD Zugabe: <input checked="" type="checkbox"/> Ja <input type="checkbox"/> Nein	
Antischaum	Art: Geserp	Volumen: 30 [ml]
pH vor Sterilisation	Ist: 6.25	Soll: 7.0
pH eingestellt mit	Name: 12011	Konz.: 5% Menge: 13 ml

Sterilisation:

Steril. Gleitringdichtung	Datum:	Uhrzeit:	Dauer:	min
1. Sterilisation Fermenter	Datum: <del>24.08.02</del>	Uhrzeit: 8 <sup>00</sup>	Dauer: 60	min
2. Sterilisation Fermenter	Datum:	Uhrzeit:	Dauer:	min
pH nach Sterilisation	14.0	Reaktorgewicht nach St.	82	kg

4-14

Fermenter: Blatt 3

Kostenstelle: 103310

Vers.Nr.: 96/145702/02

Inokulation:

Inokulum 1	Herk.: Nutzer 0 Protokoll-Nr.: 1 Flasche Nr. 6510	Volumen: 1 [l]
Fermentation Beginn	Datum: [redacted] Uhrzeit: 10:54	
Inokulum 2	Herk.: Nutzer 0 Protokoll-Nr.: 1 Flasche Nr. [redacted]	Volumen: [redacted] [l]
Inokulum 2 Zeitpunkt	Datum: [redacted] Uhrzeit: [redacted]	
Fermentergewicht nach	Inokulierung 1: 10 [kg]	Inokulierung 2: [redacted] [kg]

Fermentation-Ende:

Fermentation-Ende	Datum: [redacted] Uhrzeit: 14:25
Fermenter-Gewicht	92
Korrekturmittel: Volumen nach der Fermentation	Säure 1: [redacted] Säure 2: [redacted] Säure 3: [redacted] Lauge 1: [redacted] Lauge 2: [redacted] Lauge 3: [redacted] Antischaum 1: [redacted] Antischaum 2: [redacted]
Volumen nach Ferm. von	Zufütterung 1: [redacted] Zufütterung 2: [redacted]
wie geplant	<input checked="" type="checkbox"/>
Kontamination	<input type="checkbox"/> Zeitpunkt: Vor den Animpfen <input type="checkbox"/> Vorkultur <input type="checkbox"/> <input type="checkbox"/> Nach dem Animpfen <input type="checkbox"/>
Defekt	<input type="checkbox"/> Art: [redacted]
übergeschäumt	<input type="checkbox"/> Zeitpunkt: Sterilisation: <input type="checkbox"/> Kultivierung: <input type="checkbox"/> <input type="checkbox"/> Aufheizphase <input type="checkbox"/> Vor den Animpfen <input type="checkbox"/> <input type="checkbox"/> Haltephase <input type="checkbox"/> während der Kultivierung <input type="checkbox"/> <input type="checkbox"/> Abkühlphase <input type="checkbox"/> am Ende der Kultivierung <input type="checkbox"/>
sonstiges	[redacted]

Weiterverarbeitung:

Transferleitung Sterilisat.	Datum: [redacted] Uhrzeit: [redacted] Dauer: [redacted]
Ablassleitung Sterilisat.	Datum: [redacted] Uhrzeit: 15:00 Dauer: 120
Nächster Schritt der Weiterverarbeitung:	Aufarbeitung <input type="checkbox"/> Übergeimpft auf einen Fermenter <input checked="" type="checkbox"/> An Nutzer übergeben <input type="checkbox"/> Übergeimpft auf mehrere Fermenter <input type="checkbox"/>
Volumen [l]	92 [redacted] [redacted] [redacted] [redacted] [redacted] [redacted] [redacted]
Protokoll-Nr. der nächsten Schritte	02105 1 1 1 1 1 1 1

Entsorgung:

Sterilisation Abluftfilter	Datum: [redacted] Zeit: [redacted] Dauer: [redacted]
Inaktivierung Fermenter	gesamter Inhalt <input type="checkbox"/> restl. Inhalt <input type="checkbox"/> Vol: [redacted] [l] Überstand <input type="checkbox"/> Datum: 30.08.86 Uhrzeit: 8:00 Dauer: 5 Temp: 80°C
Besonderheiten	Reinigung des Fermenters + Desinfektion
Betriebs-Ende	Datum: 30.08.86 Uhrzeit: 11:00



Anmerkungen/Besonderheiten zur Fermentation: Medium mit Ultratrac suspendieren

30 ml Tegospin

**Aufarbeitung**

Zielsetzung: .....

Feststoffabtrennung: ☐ Zentrifugation ☐ Mikrofiltration ☐ Dead-end-Filtration  
☐ Adsorbentharzabtrennung

benötigt werden ---> ☐ Filtrat/Überstand ☐ Feststoff

☐ Lyophilisation ☐ Ultrafiltration

→ Verdampfung gewünschtes Endvolumen :.....(L/al)  
 max.Temp. :.....(°C)

Extraktion: ☐ Kulturbühe ☐ Überstand ☐ Feststoff

Verteilungskoeffizient: .....

Lösungsmittel/Zusätze: .....

Phasenverhältnis: ..... Stufenzahl: .....

Zusatzprotokolle: .....

Produktspezifische Besonderheiten/weitergehende Aufarbeitungsschritte/Analytik:

.....  
 .....  
 .....  
 .....  
 .....

Toxische Eigenschaften/Sicherheitsmaßnahmen: .....

Besonderheiten der Entsorgung/Dekontamination von Mikroorganismen bzw. toxischen Produkten:

.....  
 .....

**ACHTUNG!!** Lagerzeiten von Kühlgut max.3 Arbeitstagen, von Gefriergut max.3 Monate! Nach Terminüberschreitung erfolgt Entsorgung!!

Datum/Unterschrift: [Signature]

Wägeprotokoll

## Versuchsnummer:

20120/520

**Mikroorganismus:**

5000 40

Code für Medium:

Ein-Medien

Behälterbeschriftung:

1	
2	
3	
4	

Datum:

Unterschrift:

C. T. S.

linweise:

Die Versuchsnummer wird bei Abgabe des Anmeldeprotokolls vergeben / Mediumscode und Bezeichnung des Mikroorganismus wie Im:Anmeldeprotokoll / das Protokoll kann für mehrere Behälter benutzt werden, die Nr. in der ersten Spalte ist die Behälternummer (siehe 'Behälterbeschriftung') / bei Lösungen unter Substanz z.B. Wasser aufführen und Menge in ml oder l angeben / werden Stammlösungen verwendet, ist deren Zusammensetzung beizufügen

[illegible]

4-16

Organismus: \_\_\_\_\_

Betreiber: \_\_\_\_\_

[illegible]

# Fermentations- und Arbeitsprotokoll

Technikumleitung : Dr. A. Roß Tel. 130 priv: (0 53 1) 34 35 49  
 Fermentation : H. Schürer Tel. 131 priv: (0 53 71) 79 48  
 Aufarbeitung : R. Krüdel Tel. 137 priv: (0 53 02) 45 62  
 Technik : 70 33 40 Tel. priv: 0

Antragnummer: 0049-0000 Reaktions-Nr.: 0049-0 Anmeldeort: 00-00-00

Name: Fischer Übersicht: NBI

Dienststelle: 465 Privatstelle:

Stamm/Medium: Soe 90 o Stammzucht liegt vor / ☒ liegt bei

Ziel: Vorfermenter

Prozessbeginn am: 10.15 Uhr Prozessende am: 10.00 Uhr

Startvolumen: 10.2 Volumen der Vorkultur + sonst. Zugabe: 1.2

Vorkultur: Schüttelkultur ☒ Anforderer o Biotechnikum o Reaktion Exp.Nr.

Medium: Nr. E ☒ trocken o gelöst a. kg/l ☒ AS

o wird vom Biotechnikum angesetzt o wird geliefert von am

Vorlage	Art	Inhalt	/ Konz.	/ Menge(l oder kg)	Pumpentyp	Rate (max)
1						
2						
3						
4						
5						

Waagen o 1 o 2 o 3 o 4 o 5 Timer o 1 o 2 o 3 o 4 o 5  
 o pH-Einstellung vor Sterilisation auf 7.0 mit a. ml KOH

Sterilisation bei 121°C für min / Inkubation h / Sterilität h

Startwerte für die Kultivierung:

Temperatur: 30 °C Belüftung: 0.1 vvm Norm<sup>3</sup>/h

Drehzahl: 150 rpm pH-Wert: > 7.0 <

Druck: mbar

pO<sub>2</sub>-Messung o nein ☒ ja

Abgemessung o nein o ja, Kanal

pO<sub>2</sub>-Regelung o nein o ja, Sollwert % Sättigung

Druckregelung o nein o ja, Sollwert mbar

Rechnererfassung o nein o ja, Exp.Nr.

Parameter ändern

nach h -> nach h -> nach h ->  
 nach h -> nach h -> nach h ->  
 nach h -> nach h -> nach h ->

Weitere Angaben siehe versiegelt ->

GVO o ja

Das Fermentationsprotokoll sollte in der der Fermentation vorangehenden Kalenderwoche fröhlich vorliegen, spätestens jedoch zwei Tage vor Beginn der Fermentation. Mündliche Fermentationsvereinbarungen werden nur bis zu diesem Zeitpunkt berücksichtigt. Der Nutzer verpflichtet sich, nicht abgesprochene Manipulationen an Geräten zu unterlassen und im Technikum die Sicherheitsvorschriften (z.B. UVV 102) einzuhalten.

Substrat und Hilfsmittel Zugabe:(nach Sterilisation)

Art	Herkunft	Vol. [ml]	Datum	Zeit	Gew. <sub>gesamt</sub> (KG)
	Flasche Nr. _____				
	Flasche Nr. _____				
	Flasche Nr. _____				
	Flasche Nr. _____				
	Flasche Nr. _____				
	Flasche Nr. _____				
	Flasche Nr. _____				

Vorlagen und Korrekturmittel:

Lauge 1:	Vol. <sub>Anfang</sub> :	Dat./Zeit:	Herk.: Flasche Nr. _____
Lauge 2:	Vol. <sub>Anfang</sub> :	Dat./Zeit:	Herk.: Flasche Nr. _____
Lauge 3:	Vol. <sub>Anfang</sub> :	Dat./Zeit:	Herk.: Flasche Nr. _____
Säure 1:	Vol. <sub>Anfang</sub> :	Dat./Zeit:	Herk.: Flasche Nr. _____
Säure 2:	Vol. <sub>Anfang</sub> :	Dat./Zeit:	Herk.: Flasche Nr. _____
Säure 3:	Vol. <sub>Anfang</sub> :	Dat./Zeit:	Herk.: Flasche Nr. _____
Antischaum 1:	Vol. <sub>Anfang</sub> :	Dat./Zeit:	Herk.: Flasche Nr. _____
Antischaum 2:	Vol. <sub>Anfang</sub> :	Dat./Zeit:	Herk.: Flasche Nr. _____
Zufütterung 1 Art:	Vol. <sub>Anfang</sub> :	Dat./Zeit:	Herk.: Flasche Nr. _____
Zufütterung 2 Art:	Vol. <sub>Anfang</sub> :	Dat./Zeit:	Herk.: Flasche Nr. _____

Regelung und Fermentationsstrategie, Startwerte:

pH-Sollwert: <u>7.0</u>	eingestellt mit : _____		pH-Regelung von ..... bis .....
pO <sub>2</sub>	Messung nein o ja <input checked="" type="checkbox"/>	Regelung nein <input checked="" type="checkbox"/> ja o	
pO <sub>2</sub> Sollwert:	Strategie: Drehzahl o Zuluft o sonstige o _____		
Temperatur <u>30</u> [°C]	Druck _____ [mbar]	Drehzahl <u>150</u> [rpm]	
Parameter:	Sollwert: [.....]	Strategie:	
Parameter:	Sollwert: [.....]	Strategie:	
Parameter:	Sollwert: [.....]	Strategie:	
Begasung: Luft	_____ l/min	<u>0.1 vvm = 60.0 l/min</u> vvm	
andere: <u>Distillat</u>	_____ l/min	_____ vvm	
	_____ l/min	_____ vvm	
	_____ l/min	_____ vvm	
Überlagerung GLRD	Luft o Dampf o		
Abgasmessung	nein <input checked="" type="checkbox"/> ja o	Kanal: <u>K1/K2</u>	
Rechnererfassung nein o	Exp.: <u>26.14.2</u>	Start-Datum: _____	Zeit: _____

# Fermentations- und Anlagentechnikprotokoll

Technikumleitung : Dr. A. Roll Td. 130 priv: (0 53 1) 34 35 49  
 Fermentation : H. Schuler Td. 131 priv: (0 53 71) 79 48  
 Anlagentechnik : R. Kriegl Td. 137 priv: (0 53 02) 45 62  
 Technik : 103 210 Td. priv: 0

4-20

Anlagennummer: 0048:0000 Reaktor-Nr.: 0019.2 Anmeldeformular: 00.00.00

Name: Gerth / Fischer Überprüfer: NBL

Dienstleider: 433 / 465 Privatleider:

Stamm/Medium: Soce 90 o Stammzucht liegt vor ☒ liegt bei

Ziel: Vorfermenter o Produktion von Epothilon (wenn anderer Stoff)

Prozessbeginn am: 10.12.10 um 10:00 Uhr Prozessende am: 10.12.10 um 10:00 Uhr

Startvolumen: 100 l Volumen der Vorkultur + sonst. Zugabe: 10 l

Vorkultur: Schüttelkultur o Anforderer o Biotechnikum o Reaktor/Exp.Nr.

Medium: Nr. E o trocken o gelöst ca. kg/l ☒ 15

o wird vom Biotechnikum angesetzt o wird geliefert von am

Vorgabe	Art	Inhalt	/ Konz.	/ Menge (l oder kg)	Pumptyp	Rate (max.)
1	Lauge	KOH	10 l			
2	A.S.	Agar, 15, 100 ml				
3						
4						
5						

Waagen o 1 o 2 o 3 o 4 o 5 Timer o 1 o 2 o 3 o 4 o 5  
 o pH-Einstellung vor Sterilisation auf 7,6 mit ca. ml KOH

Sterilisation bei 121°C für min / fraktioniert h / Sterilität h

Startwerte für die Kultivierung:

Temperatur: 30 °C Belüftung: 0,1 von Nm³/h  
 Drehzahl: 200 rpm pH-Wert: 7,0 konstant  
 Druck: mbar

pO2-Messung o nein ☒ ja  
 Abgasmessung o nein o ja Kanal  
 pO2-Regelung o nein o ja Sollwert % Sättigung  
 Druckregelung o nein o ja Sollwert mbar  
 Rechneransteuerung o nein o ja Exp.Nr.

GVO o ja

Parameter ändern  
 nach h -> nach h -> nach h ->  
 nach h -> nach h -> nach h ->  
 nach h -> nach h -> nach h ->

Weitere Angaben siehe unten

Das Fermentationsprotokoll sollte in der der Fermentation vorausgehenden Kalenderwoche fertig vorliegen, spätestens jedoch zwei Tage vor Beginn der Fermentation. Mündliche Fermentationsvereinbarungen werden nur bis zu diesem Zeitpunkt berücksichtigt. Der Nutzer verpflichtet sich, nicht abgesprochene Manipulationen an Geräten zu unterlassen und im Technikum die Sicherheitsvorschriften (z.B. UVV 102) einzuhalten.

Anmerkungen/Besonderheiten zur Fermentation: Medium mit Ultrakimax suspendieren

70 ml Tryptone i. NAD-Zugabe

Anarbeitung

Zielsetzung: evk

Feststoffabtrennung:

- ☐ Zentrifugation
- ☐ Mikrofiltration
- ☐ Dead-end-Filtration
- ☐ Adsorberharzabtrennung

benötigt werden --->

- ☐ Filtrat/Überstand
- ☐ Feststoff

☐ Lyophilisation

- ☐ Ultrafiltration

Verdampfung gewünschtes Endvolumen :.....(L/ml)

max.Temp. :.....(°C)

Extraktion:

- ☐ Kulturbrühe
- ☐ Überstand
- ☐ Feststoff

Verteilungskoeffizient:.....

Lösungsmittel/Zusätze:.....

Phasenverhältnis:..... Stufenzahl:.....

Zusatzprotokolle:.....

Produktspezifische Besonderheiten/weitergehende Aufarbeitungsschritte/Analytik:

Toxische Eigenschaften/Sicherheitsmaßnahmen:.....

Besonderheiten der Entsorgung/Dekontamination von Mikroorganismen bzw. toxischen Produkten:

ACHTUNG!! Lagerzeiten von Kühlgut max. 3 Arbeitstagen, von Gefriergut max. 3 Monate! Nach Terminüberschreitung erfolgt Entsorgung!!

Datum/Unterschrift: C. Fischer

Wägeprotokoll 7 150/2

145102/03

Score 40

E-lection

\_\_\_\_\_

4	
---	--



1852

Die Versuchsnummer wird bei Abgabe des Anmeldeprotokolls vergeben / Mediumscode und Bezeichnung des Mikroorganismus wie im Anmeldeprotokoll / das Protokoll kann für mehrere Behälter benutzt werden, die Nr. in der ersten Spalte ist die Behälternummer (siehe 'Behälterbeschriftung') / bei Lösungen unter Substanz z.B. Wasser aufführen und Menge in ml oder l angeben / werden Stammlösungen verwendet, ist deren Zusammensetzung beizufügen

[illegible]



FERMENTER: Blatt 1

Kostenstelle: 103310

Vers.Nr.: 96/145702/03

Betreiber: IL. Jute Betreuer: StelinskiOrganismus: Succ 90Kulturführung: Aerob: ☐ Anaerob: ☐ Phototroph: ☐Prozeßführung: Batch: ☐ Feed-Batch: ☐ Konti: ☐

Fermenteraufbau:

Fermenter Nr. <u>100.2</u>	Verwendung: Fermentation <input checked="" type="checkbox"/> Vorlage <input type="checkbox"/>	Steriltest <input type="checkbox"/> für Protokoll-Nr. <u>1</u>
Sicherheitsmaßnahmen	Abluftfilter: Nein <input type="checkbox"/> Ja <input type="checkbox"/> Handschuhe tragen: Nein <input type="checkbox"/> Ja <input checked="" type="checkbox"/>	
Betrieb-Beginn	Datum: <u>          </u> Uhrzeit: <u>15:00</u>	
Rührerart	<u>3x Schenk</u>	
Sondergeräte		
Pumpe für <u>Lang</u>	Typ: <u>Flocon 500A</u> Pumprate: <u>          </u> Durchm. Schlauch: <u>          </u>	
Pumpe für <u>          </u>	Typ: <u>          </u> Pumprate: <u>          </u> Durchm. Schlauch: <u>          </u>	
<u>Medien mit Ultraschall</u>		

Elektroden:

pH-Elektrode	Nr.: <u>200.7</u> Puffer 1: <u>          </u> Puffer 2: <u>          </u> Poti/ mV: <u>5.16</u> Poti/ mV: <u>280</u>
pH-Elektrode	Nr.: <u>          </u> Puffer 1: <u>          </u> Puffer 2: <u>          </u> Poti/ mV: <u>          </u> Poti/ mV: <u>          </u>
pO <sub>2</sub> -Elektrode	Nr.: <u>500.7</u> Nr.: <u>520.7</u>

Reaktorgewicht:

Gesamtgewicht

Sollgewicht	<u>90.2</u>	[KG]
leer		<u>-5</u> [KG]
Wassermenge	<u>68</u> [l]	<u>63</u> [KG]
Medium-Zugabe	Name: <u>          </u> Herk.: Nutzer 0 SE: <u>03/02</u>	<u>84</u> [KG]
	XAD Zugabe: <input checked="" type="checkbox"/> Ja <input type="checkbox"/> Nein	<u>82</u>
Antischaum	Art: <u>Desuper</u> Volumen: <u>20</u> [ml]	
pH vor Sterilisation	Ist: <u>6.82</u> Soll: <u>2.6</u>	
pH eingestellt mit	Name: <u>120M</u> Konz.: <u>5N</u> Menge: <u>50</u> ml	

Sterilisation:

Steril. Gleitringdichtung	Datum: <u>          </u> Uhrzeit: <u>9:35</u> Dauer: <u>45</u> min
1. Sterilisation Fermenter	Datum: <u>          </u> Uhrzeit: <u>12:20</u> Dauer: <u>60</u> min
2. Sterilisation Fermenter	Datum: <u>          </u> Uhrzeit: <u>          </u> Dauer: <u>          </u> min
pH nach Sterilisation	<u>7.08</u> Reaktorgewicht nach St. <u>82</u> kg

Substrat und Hilfsmittel Zugabe:(nach Sterilisation)

Art	Herkunft	Vol. [ml]	Datum	Zeit	Gew. <sub>gesamt</sub> (KG)
	Flasche Nr. _____				
	Flasche Nr. _____				
	Flasche Nr. _____				
	Flasche Nr. _____				
	Flasche Nr. _____				
	Flasche Nr. _____				
	Flasche Nr. _____				

Vorlagen und Korrekturmittel:

Lauge 1: <i>120M 10%</i>	Vol. <sub>Anfang</sub> : <i>1500</i>	Dat./Zeit: <i>27.9/10<sup>10</sup></i>	Herk.: Flasche Nr. <i>533</i>
Lauge 2:	Vol. <sub>Anfang</sub> :	Dat./Zeit:	Herk.: Flasche Nr. _____
Lauge 3:	Vol. <sub>Anfang</sub> :	Dat./Zeit:	Herk.: Flasche Nr. _____
Säure 1:	Vol. <sub>Anfang</sub> :	Dat./Zeit:	Herk.: Flasche Nr. _____
Säure 2:	Vol. <sub>Anfang</sub> :	Dat./Zeit:	Herk.: Flasche Nr. _____
Säure 3:	Vol. <sub>Anfang</sub> :	Dat./Zeit:	Herk.: Flasche Nr. _____
Antischaum 1: <i>Tegospin</i>	Vol. <sub>Anfang</sub> : <i>500</i>	Dat./Zeit: <i>27.9/10<sup>10</sup></i>	Herk.: Flasche Nr. <i>554</i>
Antischaum 2:	Vol. <sub>Anfang</sub> :	Dat./Zeit:	Herk.: Flasche Nr. _____
Zufütterung 1 Art:	Vol. <sub>Anfang</sub> :	Dat./Zeit:	Herk.: Flasche Nr. _____
Zufütterung 2 Art:	Vol. <sub>Anfang</sub> :	Dat./Zeit:	Herk.: Flasche Nr. _____

Regelung und Fermentationsstrategie, Startwerte:

pH-Sollwert:	eingestellt mit : pH-Regelung von <i>7.0</i> bis .....	
pO <sub>2</sub>	Messung nein <input type="radio"/> ja <input checked="" type="radio"/>	Regelung nein <input checked="" type="radio"/> ja <input type="radio"/>
pO <sub>2</sub> Sollwert:	Strategie: Drehzahl <input type="radio"/> Zuluft <input type="radio"/> sonstige <input type="radio"/> _____	
Temperatur <i>30</i> [°C]	Druck _____ [mbar]	Drehzahl <i>200</i> [rpm]
Parameter:	Sollwert: _____ [.....]	Strategie:
Parameter:	Sollwert: _____ [.....]	Strategie:
Parameter:	Sollwert: _____ [.....]	Strategie:
Begasung: Luft	<i>10</i> l/min	<i>0,1</i> vvm
andere: _____	_____ l/min	_____ vvm
_____	_____ l/min	_____ vvm
_____	_____ l/min	_____ vvm
Überlagerung GLRD	Luft <input checked="" type="radio"/> Dampf <input type="radio"/>	
Abgasmessung	nein <input type="radio"/> ja <input checked="" type="radio"/>	Kanal: <i>2</i>
Rechnererfassung nein <input type="radio"/>	Exp.: <i>A61453</i>	Start-Datum: <i>[redacted]</i> Zeit: <i>10<sup>00</sup></i>


4-25

Fermenter: Blatt 3


Kostenstelle: 103310

Vers.Nr.: 96/145702/03


Inokulation:

Inokulum 1	Herk.: Nutzer 0 Protokoll-Nr.: 02104 Flasche Nr. 592	Volumen: 9 [l]
Fermentation Beginn	Datum:  Uhrzeit: 1025	
Inokulum 2	Herk.: Nutzer 0 Protokoll-Nr.: _/ _/ _ Flasche Nr. _	Volumen: [l]
Inokulum 2 Zeitpunkt	Datum: _/ _/ _ Uhrzeit: _	
Fermentergewicht nach	Inokulierung 1: 9.1 [kg]	Inokulierung 2: [kg]


Fermentation-Ende:

Fermentation-Ende	Datum:  Uhrzeit: 1310
Fermenter-Gewicht	90 kg
Korrekturmittel: Volumen nach der Fermentation	Säure 1: _____ Säure 2: _____ Säure 3: _____ Lauge 1: 1400 Lauge 2: _____ Lauge 3: _____ Antischaum 1: _____ Antischaum 2: _____
Volumen nach Ferm. von	Zufütterung 1: _____ Zufütterung 2: _____
wie geplant	<input checked="" type="checkbox"/>
Kontamination	<input type="checkbox"/> Zeitpunkt: Vor den Animpfen <input type="checkbox"/> Vorkultur <input type="checkbox"/> <input type="checkbox"/> Nach dem Animpfen <input type="checkbox"/>
Defekt	<input type="checkbox"/> Art: _____
übergeschäumt	<input type="checkbox"/> Zeitpunkt: Sterilisation: _____ Kultivierung: _____ Aufheizphase <input type="checkbox"/> Vor den Animpfen <input type="checkbox"/> Haltephase <input type="checkbox"/> während der Kultivierung <input type="checkbox"/> Abkühlphase <input type="checkbox"/> am Ende der Kultivierung <input type="checkbox"/>
sonstiges	

Weiterverarbeitung:

Transferleitung Sterilisat.	Datum:  Uhrzeit: 1100 Dauer: 120 min
Ablasseleitung Sterilisat.	Datum: _____ Uhrzeit: _____ Dauer: _____
Nächster Schritt der Weiterverarbeitung:	Aufarbeitung <input type="checkbox"/> Übergeimpft auf einen Fermenter <input checked="" type="checkbox"/> An Nutzer übergeben <input type="checkbox"/> Übergeimpft auf mehrere Fermenter <input type="checkbox"/>
Volumen [l]	100
Protokoll-Nr. der nächsten Schritte	02106 _/ _/ _ _/ _/ _ _/ _/ _

Entsorgung:

Sterilisation Abluftfilter	Datum: _____ Zeit: _____ Dauer: _____
Inaktivierung Fermenter	gesamter Inhalt <input type="checkbox"/> restl. Inhalt <input type="checkbox"/> Vol: _____ [l] Überstand <input type="checkbox"/> Datum: _____ Uhrzeit: _____ Dauer: _____ Temp.: _____
Besonderheiten	
Betriebs-Ende	Datum:  Uhrzeit: 1500



# Fermentations- und Anlagungsprotokoll

4-27

Technikumleitung : Dr. A. Rob. Td. 130. priv: (0 53 1) 34 35 49  
 Fermentation : H. Schuler. Td. 131. priv: (0 53 71) 79 48  
 Aufarbeitung : R. Krütsch. Td. 137. priv: (0 53 07) 45 62  
 Technik : 10 53 10. Td. . priv: 0

Auftragsnummer: 0045.0000 Reaktor-Nr.: 0003.7. Anmeldezeit: 00.00.00

Name: Fischer. Bericht/Titel: NBI

Dienststelle: 465 Privatstelle:

Stamm/Medium: Soe 90 o Stammzucht liegt vor / ☒ liegt bei

Ziel: Produktion von Epothilon

Prozessbeginn am: 10.15 Uhr Prozessende am: 10.00 Uhr

Startvolumen: 750 l Volumen der Vorkultur + sonst. Zugabe: 70 l

Vorkultur: Schüttelkultur o Anforderer o Biotechnikum o Reaktor/Exp.Nr.

Medium: Nr. E ☒ trocken o gelöst ca. kg/l ☒ AS

o wird vom Biotechnikum angesetzt o wird geliefert von am

Vorgang	Art	Inhalt	Konz.	Menge (l oder kg)	Pumpentyp	Rate (max)
1	Lauge	KOH	10%			
2	A.S.	Tegospin				
3						
4						
5						

Waagen o 1 o 2 o 3 o 4 o 5 Timer o 1 o 2 o 3 o 4 o 5  
 o pH-Einstellung vor Sterilisation auf 7.6 mit ca. ml KOH

Sterilisation bei 121°C für min / fraktioniert h / Sterilzeit h

Startwerte für die Kultivierung

Temperatur: 30 °C Belüftung: 0.1 von Nm³/h

Drehzahl: 200 rpm pH-Wert: > 7.0 konstant

Druck: mbar

pO2-Messung o nein ☒ ja

Abgemessung o nein o ja, Kanal

pO2-Regelung o nein o ja, Sollwert % Sättigung

Druckregelung o nein o ja, Sollwert mbar

Rechnererfassung o nein o ja, Exp.Nr.

Parameter ändern  
 nach h -> nach h ->  
 nach h -> nach h ->  
 nach h -> nach h ->

weitere Angaben siehe unten

GVO o ja

Das Fermentationsprotokoll sollte in der der Fermentation vorangehenden Kalenderwoche fertig vorliegen, spätestens jedoch zwei Tage vor Beginn der Fermentation. Mündliche Fermentationsreservierungen werden nur bis zu diesem Zeitpunkt berücksichtigt. Der Nutzer verpflichtet sich, nicht abgesprochene Manipulationen an Geräten zu unterlassen und im Technikum die Sicherheitsvorrichtungen (z.B. UVV 102) einzuhalten.

Anmerkungen/Besonderheiten zur Fermentation: Mertiline mit Ultrafiltrax suspendieren  
200 ml Tegosipon, KAD-Eigabe

### Aufarbeitung

Zielsetzung: Nach Rücksprache mit Hr. Steinwetz

#### Feststoffabtrennung:

- ☐ Zentrifugation      ☐ Mikrofiltration      ☐ Dead-end-Filtration  
☐ Adsorberharzabtrennung

#### benötigt werden -->

- ☐ Filtrat/Überstand      ☐ Feststoff

#### o Lyophilisation

- ☐ Ultrafiltration

Verdampfung      gewünschtes Endvolumen      :.....(L/ml)

max.Temp.      :.....(°C)

#### Extraktion:

- ☐ Kulturbrothe      ☐ Überstand      ☐ Feststoff

Verteilungskoeffizient:.....

Lösungsmittel/Zusätze:.....

Phasenverhältnis:..... Stufenzahl:.....

Zusatzprotokolle:.....

Produktspezifische Besonderheiten/weitergehende Aufarbeitungsschritte/Analytik:

Toxische Eigenschaften/Sicherheitsmaßnahmen:.....

Besonderheiten der Entsorgung/Dekontamination von Mikroorganismen bzw. toxischen Produkten:

**ACHTUNG!!** Lagerzeiten von Kühlgut max. 3 Arbeitstagen, von Gefriergut max. 3 Monate! Nach Terminüberschreitung erfolgt Entsorgung!!

Datum/Unterschrift:                     

C. FISCH

10

Soc 40

== Medication

1	
2	
3	
4	



C. 1558

4-29

[illegible]

FERMENTER: Blatt 1

Kostenstelle: 103310

Vers.Nr.: 96/115/02/06

Betreiber: K. Jule

Betreuer: S. Leisner

Organismus: S. cerevisiae

Kulturführung: Aerob: ☐ Anaerob: ☐ Phototroph: ☐Prozeßführung: Batch: ☐ Feed-Batch: ☐ Konti: ☐

Fermenteraufbau:

Fermenter Nr. 901	Verwendung: Fermentation <input checked="" type="radio"/> Vorlage <input type="radio"/>	Sterilttest <input type="radio"/> für Protokoll-Nr. ___/___
Sicherheitsmaßnahmen	Abluftfilter: Nein <input type="radio"/> Ja <input type="radio"/> Handschuhe tragen: Nein <input type="radio"/> Ja <input checked="" type="radio"/>	
Betrieb-Beginn	Datum: [redacted] Uhrzeit: 10:00	
Rührerart	2x 3-Blatt	
Sondergeräte		
Pumpe für Lauge	Typ: Pumprate: Durchm. Schlauch:	
Pumpe für	Typ: Pumprate: Durchm. Schlauch:	
	Medium mit Urturbesatz	

Elektroden:

pH-Elektrode	Nr.: 200.6 Puffer 1: 7 Poti/ mV: -5,1 Puffer 2: 4 Poti/ mV: 55,7
pH-Elektrode	Nr.: Puffer 1: Poti/ mV: Puffer 2: Poti/ mV:
pO <sub>2</sub> -Elektrode	Nr.: 150.16 Nr.: 150.15

Reaktorgewicht:

Gesamtgewicht

Sollgewicht	680	[KG]
leer		[KG]
Wassermenge	580 [l]	[KG]
Medium-Zugabe	Name: Herk.: Nutzer 0 SE: 03/04 XAD Zugabe: <input checked="" type="radio"/> Ja <input type="radio"/> Nein	680 [KG]
Antischaum	Art: Osmix Volumen: 200 [ml]	
pH vor Sterilisation	Ist: 6,18 Soll: 7,0	
pH eingestellt mit	Name: 1404 Konz.: 5N Menge: 360 ml	

Sterilisation:

Steril. Gleitringdichtung	Datum: [redacted] Uhrzeit: 9:40 Dauer: 50 min
1. Sterilisation Fermenter	Datum: [redacted] Uhrzeit: 11:10 Dauer: 60 min
2. Sterilisation Fermenter	Datum: Uhrzeit: Dauer: min
pH nach Sterilisation	6,79 Reaktorgewicht nach St. kg



4-31

Substrat und Hilfsmittel Zugabe:(nach Sterilisation)

Art	Herkunft	Vol. [ml]	Datum	Zeit	Gew. <sub>gesamt</sub> (KG)
AS	Flasche Nr. 96/0537	300		7:05	
	Flasche Nr. _____				
	Flasche Nr. _____				
	Flasche Nr. _____				
	Flasche Nr. _____				
	Flasche Nr. _____				
	Flasche Nr. _____				

Vorlagen und Korrekturmittel:

Lauge 1: KOH 10%	Vol. <sub>Anfang</sub> : 4700	Dat./Zeit: 13:45	Herk.: Flasche Nr. 539
Lauge 2:	Vol. <sub>Anfang</sub> :	Dat./Zeit:	Herk.: Flasche Nr. _____
Lauge 3:	Vol. <sub>Anfang</sub> :	Dat./Zeit:	Herk.: Flasche Nr. _____
Säure 1:	Vol. <sub>Anfang</sub> :	Dat./Zeit:	Herk.: Flasche Nr. _____
Säure 2:	Vol. <sub>Anfang</sub> :	Dat./Zeit:	Herk.: Flasche Nr. _____
Säure 3:	Vol. <sub>Anfang</sub> :	Dat./Zeit:	Herk.: Flasche Nr. _____
Antischaum 1: Tego 100	Vol. <sub>Anfang</sub> : 1500	Dat./Zeit:	Herk.: Flasche Nr. 537
Antischaum 2:	Vol. <sub>Anfang</sub> :	Dat./Zeit:	Herk.: Flasche Nr. _____
Zufütterung 1 Art:	Vol. <sub>Anfang</sub> :	Dat./Zeit:	Herk.: Flasche Nr. _____
Zufütterung 2 Art:	Vol. <sub>Anfang</sub> :	Dat./Zeit:	Herk.: Flasche Nr. _____

Regelung und Fermentationsstrategie, Startwerte:

pH-Sollwert:	eingestellt mit :		pH-Regelung von ..... bis 7.0	
pO <sub>2</sub>	Messung	nein o ja <input checked="" type="checkbox"/>	Regelung	nein <input checked="" type="checkbox"/> ja o
pO <sub>2</sub> Sollwert: 30%	Strategie: Drehzahl <input checked="" type="checkbox"/> Zuluft o sonstige o _____			
Temperatur 30 [°C]	Druck 300	[mbar]	Drehzahl 100	[rpm]
Parameter:	Sollwert:	[.....]	Strategie:	
Parameter:	Sollwert:	[.....]	Strategie:	
Parameter:	Sollwert:	[.....]	Strategie:	
Begasung: Luft	70	l/min	0.4?	vvm
andere:		l/min		vvm
		l/min		vvm
		l/min		vvm
Überlagerung GLRD	Luft <input checked="" type="checkbox"/>	Dampf o		
Abgasmessung	nein o	ja <input checked="" type="checkbox"/>	Kanal: 2	
Rechnererfassung nein o	Exp.: A61456	Start-Datum:	Zeit: 13:10	

Fermenter: Blatt 3

Kostenstelle: 402210

Vers.Nr.: 96/145/02/06

## Inokulation:

Inokulum 1	Herk.: Nutzer 0 Protokoll-Nr.: 02103 Flasche Nr. _____	Volumen: 100 [l]
Fermentation Beginn	Datum: _____ Uhrzeit: 13:10	
Inokulum 2	Herk.: Nutzer 0 Protokoll-Nr.: ____/____ Flasche Nr. _____	Volumen: [l]
Inokulum 2 Zeitpunkt	Datum: _____ Uhrzeit: _____	
Fermentergewicht nach	Inokulierung 1: 780 kg	Inokulierung 2: [kg]

## Fermentation-Ende:

Fermentation-Ende	Datum: _____ Uhrzeit: 08:5
Fermenter-Gewicht	
Korrekturmittel: Volumen nach der Fermentation	Säure 1: _____ Säure 2: _____ Säure 3: _____ Lauge 1: 3000 Lauge 2: _____ Lauge 3: _____ Antischaum 1: 100 Antischaum 2: _____
Volumen nach Ferm. von	Zufütterung 1: _____ Zufütterung 2: _____
wie geplant	<input checked="" type="checkbox"/>
Kontamination	<input type="checkbox"/> Zeitpunkt: Vor den Animpfen <input type="checkbox"/> Vorkultur <input type="checkbox"/> <input type="checkbox"/> Nach dem Animpfen <input type="checkbox"/>
Defekt	<input type="checkbox"/> Art: _____
übergeschäumt	<input type="checkbox"/> Zeitpunkt: Sterilisation: _____ Kultivierung: _____ Aufheizphase <input type="checkbox"/> Vor den Animpfen <input type="checkbox"/> Haltephase <input type="checkbox"/> während der Kultivierung <input type="checkbox"/> Abkühlphase <input type="checkbox"/> am Ende der Kultivierung <input type="checkbox"/>
sonstiges	

## Weiterverarbeitung:

Transferleitung Sterilisat.	Datum: _____ Uhrzeit: 10:55	Dauer: _____
Ablassleitung Sterilisat.	Datum: _____ Uhrzeit: _____	Dauer: _____
Nächster Schritt der Weiterverarbeitung:	Aufarbeitung <input checked="" type="checkbox"/> Übergeimpft auf einen Fermenter <input type="checkbox"/> An Nutzer übergeben <input type="checkbox"/> Übergeimpft auf mehrere Fermenter <input type="checkbox"/>	
Volumen [l]	780	
Protokoll-Nr. der nächsten Schritte	031/____/____/____/____/____/____/____/____	

## Entsorgung:

Sterilisation Abluftfilter	Datum: _____ Zeit: _____ Dauer: _____
Inaktivierung Fermenter	gesamter Inhalt <input type="checkbox"/> restl. Inhalt <input type="checkbox"/> Vol: _____ [l] Überstand <input type="checkbox"/> Datum: _____ Uhrzeit: _____ Dauer: _____ Temp.: _____
Besonderheiten	
Betriebs-Ende	Datum: _____ Uhrzeit: 15:00

Verlaufs-Protokoll: Fern-Nr.: -----

Kostenstelle: -----

Vers.Nr.: 96/ ----- /02/

Organismus: -----

Betreiber: -----

Datum	Zeit	Probe Vol.					Zugabe	Para- meter	Wert alt	Wert neu	Bemerkungen
	8 <sup>20</sup>	200 ml	UBC								
	15 <sup>35</sup>							Druck	500 mm	300 mm	Druckwert 300 mm = 0,6 bar eingestellt
	7 <sup>20</sup>	200 ml	agall								Wird direkt beim 300 mm Umschicken
"	"							Druck	200 mm	0 mm	
	2 <sup>50</sup>	200 ml	UBC								
	8 <sup>50</sup>										Polierleistung ein 30% mit 100 mm mit 300 mm
	10 <sup>00</sup>										Polierleistung ein mit 100 mm mit 300 mm
	7 <sup>45</sup>	200 ml	agall								Polierleistung ein mit 100 mm mit 300 mm
	8 <sup>45</sup>	200 ml	UBC								
	14 <sup>30</sup>						AS				ca. 20 mm
	17 <sup>30</sup>	<del>200 ml</del>					AS				2300 - 1 mm eingestellt
	12 <sup>35</sup>	200 ml									
	7 <sup>50</sup>	200 ml									
	8 <sup>50</sup>	200 ml									
	14 <sup>00</sup>	300									

8<sup>00</sup> 300  
8<sup>10</sup> 300

~~7-20-28~~

Stammzahlform

Socle 90 A4 27.11  
 Socle 90 A3 27.11  
 Socle 90 A2 27.11  
 Socle 90 W11 27.11  
 Socle 1275 27.11  
 Socle 1798 AS 27.11  
 Socle 1300 27.11

Bestellungen

Flugzeugteile

Wohn 1-50 Socle 90 11.11

~~7-20-28~~

P76 Erbeidig

Probe F900 17.11

Wohn 1-31

9.70

P78 Erneite V

800 F900 Probe V

Aufbau in Flug

F-111C

Dead-End-Filtration

Kostenstelle: 103310

Vers.Nr.: 96/0745 / 03 / 06

Bearbeiter: Prozeß: RSS.th Analysen: .....

Analysen-Protokoll: \_ \_ / \_ \_

Stamm / Medium: E.96 m Girth Mx

Besondere Sicherheitsmaßnahmen: Handseinh.

Zielsetzung: XAD Gewinn

Grunddaten

Anlagentyp	Prozeßfilter Wiga EFT 60/180		
Modifikationen	Process filter		
Betriebsbeginn	Datum: <u>          </u>	Uhrzeit: <u>08:00</u>	
Prozeßbeginn	Datum: <u>          </u>	Uhrzeit: <u>08:20</u>	

Bearbeitetes Material

Art:	<u>Firmeninteraktion</u>		
Vol.: <u>750</u>	[l]	Temp. <u>          </u>	[°C]
pH <u>          </u>	[-]	oD <sub>Medium</sub> <u>          </u>	[-]
		Herk.: Nutzer o	Protokoll-Nr: <u>02/06</u>
		Feststoffanteil	[g/l]

Filterhilfsmittel

Art:	eingewogen <u>          </u>		[kg]	gelöst in <u>          </u>	[l]
Konz <sub>Filterhilfsmittel</sub> <u>          </u>	[g/l]	konti.Dosierung <u>          </u>	[l/h]	Pumpentyp: <u>          </u>	

Geräteparameter

Eingesetzte Filtermedien/Siebgewebe: <u>250 µm</u>
--

Produkte

Gesamtlaufzeit <u>5,5</u>	[h]	mittlerer Flux <u>120</u>	[l/h]	oD <sub>Klarlauf</sub> <u>          </u>	[-]
Feststoff <sub>Klarlauf</sub> <u>          </u>	[mg/l]	Endvol. <sub>Klarlauf</sub> <u>          </u>	[l]	ProdKonz <sub>Klarlauf</sub> <u>          </u>	[.../ml]
Feststoff <sub>Konzentrat</sub> <u>          </u>	[mg/l]	Endvol. <sub>Konzentrat</sub> <u>          </u>	[kg/l]	ProdKonz <sub>Konzentrat</sub> <u>          </u>	[.../ml]
Konz.grad <u>          </u>	[%]	Konz.faktor <u>          </u>	[-]		

Verbleib der Produkte

	Klarlauf	Konzentrat <u>XAD</u>
Weiterverarbeitung	Protokoll-Nr. <u>          </u>	Protokoll-Nr. <u>03/</u>
An/Für Nutzer	übergeben o	eingelagert o
Inaktivierung (Art)		übergeben o
Entsorgung (Art)	<u>X</u>	eingelagert o

Versuchsprotokoll:

Bemerkungen.

Feststoffextraktion / Desorption

Kostenstelle: 703510

Vers.Nr.: 96/0747 / 03 / C 7

Bearbeiter: Prozeß: Roth Analysen: .....

Analysen-Protokoll: \_ \_ / \_ \_

Stamm / Medium: f. 940 Gerth MxT Epithelion

Besondere Sicherheitsmaßnahmen: .....

Zielsetzung: .....

**Grunddaten**

Anlagentyp	Sitz Prozeßf. für EFT 60/180		
Modifikationen			
Betriebsbeginn	Datum: <del>.....</del>	Uhrzeit: 14:30	
Prozeßbeginn	Datum: <del>.....</del>	Uhrzeit: 15:00	

**Bearbeitetes Material**

Art:	12g Produkt erwart.		
Vol.: [l]	Menge: [g]	Herk.: Nutzer o	Protokoll-Nr: /

**Extraktions-/ Desorptionsverlauf**

Probe Nr.	Menge [kg] ; [l]	Lösungsmittel		Kontaktzeit [min] ; [h]	Produkt [mg] ; [g]	Phasentrennung Q <sub>m</sub> [l/h]
		Art	Menge [l]			
Räuhung	ca. 15	Methanol 1)	45	über Nacht		
1. El.	15 kg	Methanol	15	3h		
2. El.	15 kg	Methanol	15	über das Wochenende		
3. El.	15 kg	Methanol	15	3h		
4. El.	15 kg	Methanol	30L	3h		
5. El.	15 kg	Methanol	30L	ü. Nacht		
6. El.	15 kg	Methanol	30L	3h		
7.	15 kg	Methanol	15L	2h		

1) Methanol + H<sub>2</sub>O Gemische 1 zu 2 !

## Verbleib der Produkte:

	Wässrige Phase	Organische Phase
Weiterverarbeitung	Protokoll-Nr.: /	Protokoll-Nr.: 03 / 08
An/Für Nutzer	übergeben o                      eingelagert o	übergeben o                      eingelagert o
Inaktivierung (Art)		
Entsorgung (Art)	o	

## Versuchsende

Prozeßende	Datum: [REDACTED]	Uhrzeit: 17:00
Betriebsende	Datum: [REDACTED]	Uhrzeit: 17:30

Bemerkungen:





4-40

### Batch-Reaktor

Vakuum	Konzentrattemp. [°C]
Destillatleistung	Heizmedium [°C]
Direkt-	Dampfdruck x 10 <sup>3</sup> [hPa]
bedampfung	Destillatleistung [l/h] Konzentrattemp. [°C]

### Produkte

End-Werte	Endvolumen [l]	Konz <sub>Produkt</sub> [ ]	Bemerkungen
Konzentrat	26 L		
Sumpf	130 L		

Ausbeute: Produktmenge Konzentrat/Menge Rohprodukt x 100 =	[%]
--	-----

### Verbleib der Produkte

	Konzentrat	Sumpf
Weiterverarbeitung	Protokoll-Nr.: 031	Protokoll-Nr.: 1
An/Für Nutzer	übergeben o eingelagert o	übergeben o eingelagert o
Inaktivierung (Art)		
Entsorgung (Art)		angetrocknet

### Versuchsende

Prozeßende	Datum: [redacted]	Uhrzeit 14:10
Betriebsende	Datum: [redacted]	Uhrzeit 16:00

Bemerkungen: Konzentrat Düs.w. Verdampfer in Rot. u. f. führt, in Rot. mithin. L. abzutrennen.  
Wasserbad Temp. 28°C, Max. Vakuum 11 m. Bar.

Flüssig-(Beß-)Flüssig - Extraktion

Kostenstelle: 1033.40

Vers.Nr.: 96/914.5/03/09

Bearbeiter: Prozeß: Perf. th ..... Analysen: .....

Analysen-Protokoll: \_\_ / \_\_

Stamm / Medium: Konzentrat f. Gae. North Epiphylla. Anisomorph. Kieselalgen

Besondere Sicherheitsmaßnahmen:.....

Zielsetzung: GFR extraktiv zur Wirtsephase.

## Grunddaten

Extraktortyp	Gieseler 4.5 qvf / SAN Separator	
Modifikationen		
Betriebsbeginn	Datum: <del>15.08.2006</del>	Uhrzeit: 14:30
Prozeßbeginn	Datum: <del>15.08.2006</del>	Uhrzeit: 15:00

## Bearbeitetes Material

Art: <i>Konzentrat</i>		<i>Dünnschichtverdampfer - Retort.-s Verdampfer</i>	
Vol.: <i>20</i>	[l]	Herkunft: Nutzer <i>o</i>	Protokoll-Nr.: <i>/</i>
Temp.: <i></i>	[°C]	pH <i></i>	[ <i>-</i> ] Konz. <small>Produkt</small> <i></i> [mg / <i>...</i> ]

## Zusätze

Zusatzstoffe	Art: Ammoniak 12.2. PH >
pH-Einstellungen	mit: 12.4 Essigsäure 20ml Lauf: 6,95
Lösungsmittel	Art: Ethylacetat

## Betriebsdaten

Durchfluß Phase <sub>Wassrige</sub>	[l/h]	Durchfluß Phase <sub>Organ</sub>	[l/h]	Phasen- verhältnis	1 zu 1 [-]
--	-------	-------------------------------------	-------	-----------------------	------------

### Extraktionsverlauf

Probe Nr.	Datum Uhrzeit	Stufe	Konz. [mg/l]		Vol. [l]		Bemer- kungen
			Phase <sub>Wässrig</sub>	Phase <sub>Organ</sub>	Phase <sub>Wässrig</sub>	Phase <sub>Organ</sub>	
1 Gehrst		1.			20	20	sch. Nahrung
2 Gehrst		2.					
		3.					
Start	842	4.			ca. 202 Crustaceen		
		5.					
		6.					

## Phasentrennung

Dekantierung	o	Filtration (Koaleszenz)	o	Zentrifugation
--------------	---	-------------------------	---	----------------

4-42

### Volumen der Produkte

Phase <small>Wässrige</small>	SA1 Endvol: 43 [l]	Phase <small>Organisch</small>	Endvol SA1 46 L [l]
-------------------------------	--------------------	--------------------------------	---------------------

### Konditionierung der Produkte

Phasen	Zusatzstoffe (z.B. Na <sub>2</sub> SO <sub>4</sub> )		pH-Einstellung	
	Art	Menge [g];[l]	mit	auf
Wässrige				
Organische				

### Verbleib der Produkte:

	Wässrige Phase	Organische Phase
Weiterverarbeitung	Protokoll-Nr.: /	Protokoll-Nr.: 03 / 10
An/Für Nutzer	übergeben o eingelagert o	übergeben o eingelagert o
Inaktivierung (Art)		
Entsorgung (Art)	×	

### Versuchsende

Prozeßende	Datum: [redacted]	Uhrzeit: 13:45
Betriebsende	Datum: [redacted]	Uhrzeit: 15:50

Bemerkungen: Auf Grund ungünstiger Trennungsvverhältnisse werden weiter

25 L Ethylacetat + 25 L Hexan hinzugefügt (16.10.96). 80 L Gesamtvolumen.

Dekantierung nicht erfolgreich.

Erneutes pH eingestellt, Trennung über Separator SA1.



4-44

### Batch-Reaktor

Vakuum	Konzentrattemp. [°C]
Destillatleistung	Heizmedium [°C]
Direkt-	Dampfdruck
bedampfung	Destillatleistung [l/h] x 10 <sup>3</sup> [hPa]
	Konzentrattemp. [°C]

### Produkte

End-Werte	Endvolumen [l]	Konz <sub>Produkt</sub> [ ]	Bemerkungen
Konzentrat	55 L	ca. 4 L	
Sumpf	42 L		

Ausbeute: Produktmenge Konzentrat/Menge Rohprodukt x 100 = [%]

### Verbleib der Produkte

	Konzentrat	Sumpf
Weiterverarbeitung	Protokoll-Nr.: /	Protokoll-Nr.: /
An/Für Nutzer	übergeben <input checked="" type="checkbox"/> eingelagert o	übergeben o eingelagert <input checked="" type="checkbox"/>
Inaktivierung (Art)		
Entsorgung (Art)		

### Versuchsende

Prozeßende	Datum: [redacted] Uhrzeit: 09:30
Betriebsende	Datum: [redacted] Uhrzeit: 11:00

Bemerkungen:

# EPOTHILON - Aufarbeitung 900L (750L AV)

XAD: ca 15L

4-45

Analytik Integral 3500 \$ 4,1 kg EPOB


XAD-Clution

EPOTHILON  
AB-Gehalt

		1: 2 (ca 30L)		(g)
1	MeOH / Wasser			
2	MeOH	15L	2,60	1,40
3	"	15L	2,10	1,23
4	"	15L	2,95	1,80
5	"	30L	2,63	1,58
6	"	30L	1,37	0,80
7	"	30L	0,50	0,32
8	"	15L	0,22	0,14
			12,37	7,27

Gesamtelutionsvolumen: 150L  $\approx$  10 BV

Gehalt  
EPOA  
[g]



Data File name: C:\HPCHEM\1\DATA\ANTJE\EPO00014.D

Method name: C:\HPCHEM\1\METHODS\SCREEN1.M

Sample Name: ~~no 9001~~ pos

Sample Info: HPLC\_MS ->

4-46

Injection Time: 10:36:36 AM

Sequence Name:


Report Style: screen1

data acquired by: Antje


on: ~~15-10-97~~

10:36:36 AM


vial:1




23.04 17886.2




23.79 9602.7



25.38 439.0




26.10 1218.6







26.92 324.4



28.64 519.7



29.87 367.6



4-47

4-4p


Integration Results

Signal 1: DAD1 A, Sig=210,4 Ref=550,100


5 $\mu$ L

Peak#	Time [min]	Type	Area [mAU*s]	Height [mAU]	Width [min]	Start [min]	End [min]
1	23.043	PV	17886.182	1359.741	0.201	22.614	23.457
2	23.791	VV	9602.701	705.784	0.206	23.457	24.265
3	25.383	PV	438.981	17.453	0.335	24.965	25.738
4	26.101	VV	1218.557	71.913	0.252	25.738	26.433
5	26.919	VV	324.419	11.990	0.364	26.720	27.240
6	28.637	PV	519.701	26.924	0.279	28.245	29.056
7	29.872	VV	367.567	23.116	0.235	29.647	30.264


Plot of TIC from Epo 900L pos




Plot of "Device A" from Epo 900L pos




Plot of Spectrum from 22.95 min (7 scans) from Epo 900L pos, subtracted




Plot of Spectrum from 23.71 min (8 scans) from Epo 900L pos, subtracted



Plot of Spectrum from 28.55 min (7 scans) from Epo 900L pos, subtracted



Plot of Spectrum from 29.82 min (7 scans) from Epo 900L pos, subtracted



XAD- Eluate 1-8 über Ethanol dampf konzentriert  
und anschließend mit EE verteilt

Wasserphase 42 l

ca 280 mg EPO in H<sub>2</sub>O Phase

EE-Phase 46 l

- Wasserphase wurde nur 1 mal extrahiert
- EE-Phase wurde im Roti konzentriert und von NCH übernommen.
- EE-Extrakt war sauer pH 5.0 und wurde mit 1M  $\text{NH}_4\text{COOCH}_3$  gepuffert. Schwierigkeiten bei der PufferEinstellung.
- Nachteil bei EE-konzentrierung dampft  $\text{NH}_3$ -Gas ab.  
Das Eluat wird sauer.
- Korrekte Pufferung mit 0,5M  $\text{K}_2\text{PO}_4$ -Puffer
- Extraktion ist neutral!

- EE-gesamt - Extrakt 246 g 407.

Gehalt EPO A = 13,82 g A

" " B = 8,40 g B

- von den 246 g Rohextrakt wurden 80 g entnommen  
und einmal mit n-Heptan ausgeschüttelt!

- n-Heptanphase : 20 g

MeOH-Phase : 60 g 408.

Essigester Extrakt 407

(2/3 davon)  
(1/3 s. 408)

4-51

wurden mit Heptan - verteilt

413.

Auswagel:

Heptanphase: 20.4 g

EPO A

30 mg

verworfen

4 B


19 mg

"


~~900.2~~ LH20 - Chromatographie von So 90. 408

900.2 4-52

Teil 1 60g Rohextrakt



8 - Spirangien  
spots for epo-  
thilone faded



Fractionierung

1-3

4+5

6-12

13- Rest (19)

Spirangene

Spirang. + EPO

~~409.~~ 409. 409

410.

411.

412.

Data File name: C:\HPCHEM\1\DATA\MITTWOCH\A000000->

Method name: C:\HPCHEM\1\METHODS\SCREEN1.M

Sample Name: Fr.6-12

Sample Info:

Injection Time: 4:02:31 PM

Sequence Name:

Report Style: screen1

data acquired by: Antje


on: [REDACTED]

4:02:31 PM


vial:1

LH20


4-53




19.63 15304.8




20.53 7617.5



23.25 894.8




23.65 648.1



Data File name: C:\HPCHEM\1\DATA\MITTWOCH\A000000->


Method name: C:\HPCHEM\1\METHODS\SCREEN1.M

26.40 1092.9




4-54


27.89 956.6




28.91 706.5




30.42 1476.6




31.59 640.9




32.60 4815.9



33.04 2642.5




33.41 1459.5






Data File name: C:\HPCHEM\1\DATA\MITTWOCH\A000000->  
Method name: C:\HPCHEM\1\METHODS\SCREEN1.M  
34.85 3627.0


4-55



36.17 2105.6



38.56 4021.5



Integration Results

4-56

Signal 1: DAD1 A, Sig=210,4 Ref=550,100


Peak#	Time [min]	Type	Area [mAU*s]	Height [mAU]	Width [min]	Start [min]	End [min]
1	19.635	VV	15304.782	955.700	0.241	18.503	20.154
2	20.529	VV	7617.517	457.590	0.249	20.154	21.403
3	22.436	VV	533.199	14.543	0.483	21.932	22.681
4	23.254	VV	894.821	30.803	0.393	22.681	23.520
5	23.651	VV	648.111	19.969	0.421	23.520	24.170
6	26.398	VV	1092.884	30.130	0.485	25.898	26.916
7	27.892	VV	956.623	28.627	0.448	27.530	28.423
8	28.909	VV	706.516	16.989	0.574	28.423	29.178
9	30.425	VV	1476.579	27.183	0.694	29.676	31.044
10	31.586	VV	640.901	16.239	0.506	31.044	31.743
11	32.604	VV	4815.902	130.967	0.544	31.743	32.884
12	33.035	VV	2642.511	127.065	0.294	32.884	33.321
13	33.407	VV	1459.489	68.039	0.298	33.321	33.810
14	34.846	VV	3626.960	57.609	0.901	34.234	35.738
15	36.168	VV	2105.642	36.447	0.792	35.738	37.210
16	38.556	VV	4021.517	61.308	0.907	38.050	40.871

A  
 B  
 C  
 D  
 400 mg  
 400 mg


5800 - 2 µg  
 1000  
 4.1 mg/ml  
 2.0 µg  
 5ul  

$$\frac{2000 + 200 + 5500}{5800 + 1000} = 400 \text{ mg}$$


Plot of TIC from Fr.6-12 pos




Plot of "Device A" from Fr.6-12 pos




Plot of Spectrum from 19.55 min (7 scans) from Fr.6-12 pos, subtracted




Plot of Spectrum from 20.48 min (9 scans) from Fr.6-12 pos, subtracted




Plot of Spectrum from 23.31 min (10 scans) from Fr.6-12 pos, subtracted



Plot of Spectrum from 26.38 min (12 scans) from Fr.6-12 pos, subtracted



Plot of Spectrum from 27.82 min (9 scans) from Fr.6-12 pos, subtracted



4-57

Prop. RP-18 Chromatographie (Merck Prop-bar)

4-58

Soce 90 411 ca 30 g

Probe wurde nicht getrocknet, da sie sich in MeOH schlecht löst.

LM M/W 6:4

Fraktion 1-30 mit Gradientensprung auf M/W=67/33

1	verworfen		
2		150 mg	414.
3			
4	verworfen		
5			
6		0,521 g	415.
7		1,892 g	<del>416.</del>
8		1,850 g	417.
9			
10		1,847 g	418.
11		1,486 g	419.
12			
13		ca 940 g	420.
14		98 mg	421.
15			
16	verworfen		
17			
18		334 mg	422.
19		Gewicht= 232 mg	Epo E 423.
20			
21		500 mg	424.
22			

(250 mg Bristol  
200 mg Boehringer  
Fraktion verbraucht) 1,3 g Silke

Fraktion:

4-5<sup>th</sup>

~~22~~

23

172 mg ~~524 mg~~ 425.

24

Gewicht = 153 mg, Epo < 426.

25

verworfen

26 }

27 }

257 mg 427.

28

Gewicht = 150 mg, Epo D 428.

29

Gewicht = 103 mg 429.

30

verworfen

# RP 18 Chromatography

So 90. 411 ~ 1. Teil cc 30 g Rohextrakt

R. 1.28

l: M/W 6:4

→ ~ 160 ml/min  
180

R. 256.

A

7

8

9

10

11

12

13

14

15

16

17

18

19

20

21

22

23

24

25

26

27

28

29

30

31

32

33

34

35

36

37

38

39

40

41

42

43

44

45

46

47

48

49

50

51

52

53

54

55

56

57

58

59

60

61

62

63

64

65

66

67

68

69

70

71

72

73

74

75

76

77

78

79

80

81

82

83

84

85

86

87

88

89

90

91

92

93

94

95

96

97

98

99

100

101

102

103

104

105

106

107

108

109

110

111

112

113

114

115

116

117

118

119

120

121

122

123

124

125

126

127

128

129

130

131

132

133

134

135

136

137

138

139

140

141

142

143

144

145

146

147

148

149

150

151

152

153

154

155

156

157

158

159

160

161

162

163

164

165

166

167

168

169

170

171

172

173

174

175

176

177

178

179

180

181

182

183

184

185

186

187

188

189

190

191

192

193

194

195

196

197

198

199

200

201

202

203

204

205

206

207

208

209

210

211

212

213

214

215

216

217

218

219

220

221

222

223

224

225

226

227

228

229

230

231

232

233

234

235

236

237

238

239

240

241

242

243

244

245

246

247

248

249

250

251

252

253

254

255

256

257

258

259

260

261

262

263

264

265

266

267

268

269

270

271

272

273

274

275

276

277

278

279

280

281

282

283

284

285

286

287

288

289

290

291

292

293

294

295

296

297

298

299

300

301

302

303

304

305

306


307

308

4-61

LC A 254,4 450,80 of AK1.D  
DATA:AK1.D

Peak#	Ret Time	Type	Width	Area	Start Time	End Time
1	3.572	BV	0.234	13828	2.544	3.977
2	4.066	VV	0.255	1568	3.977	5.115




1,454gA

0,164gB

4-62

LC A 254,4 450,80 of AK1.D  
DATA:AK1.D

Peak#	Ret Time	Type	Width	Area	Start Time	End Time
1	3.585	BV	0.227	5375	3.079	3.892
2	4.103	VB	0.274	4495	3.892	5.447



0,566g A


0,474g B




4-63

LC A 254,4 450,80 of AK1.D  
DATA:AK1.D


Peak#	Ret Time	Type	Width	Area	Start Time	End Time
1	4.101	VV	0.275	3461	3.795	5.254




(Original-File: Spectra.FRP)  
 Data File name: C:\HPCHEM\1\DATA\ANTJE\EPO000004.D  
 Method name: C:\HPCHEM\1\METHODS\SCREEN1.M  
 Sample Name: Fr.24  
 Injection Time: 2:40:26 PM  
 Sequence Name:  
 Report Style: screen1  
 data acquired by: Antje  
 vial:1  
 Sample Info: HPLC\_MS ->  
 Screening gradient  
 5090.426  
 on: 2:40:26 PM




23.30 132.3




27.27 66.4




28.00 83.1




28.54 2561.2




Plot of TIC from Fr.24 pos




Plot of "Device A" from Fr.24 pos



Plot of Spectrum from 29.22 min (7 scans) from Fr.24 pos, subtracted



Plot of Spectrum from 28.42 min (9 scans) from Fr.24 pos, subtracted



4-65

Data File name: C:\HPCHEM\1\DATA\ANTJE\EPO00005.D

Method name: C:\HPCHEM\1\METHODS\SCREEN1.M

4-66

Sample Name: Fr.28

Sample Info: HPLC\_MS ->

Injection Time: 3:27:54 PM

So 90.428

Sequence Name:


Report Style: screen1

data acquired by: Antje


on: [REDACTED]

3:27:54 PM


vial:1




23.35 607.2




24.08 573.3




29.80 3346.7




31.94 906.1




Plot of TIC from Fr.28 pos



Plot of "Device A" from Fr.28 pos




Plot of Spectrum from 23.38 min (20 scans) from Fr.28 pos




Plot of Spectrum from 27.58 min (14 scans) from Fr.28 pos



Plot of Spectrum from 29.79 min (10 scans) from Fr.28 pos



Plot of Spectrum from 31.99 min (8 scans) from Fr.28 pos





RP-Trennung von RP-Fl. 24 (70mg) von 153mg

in 2 Läufen

4-6,9

Säule = Wundwoll C18, 7µm, 20x250mm

LM = 70 MeOH  
30 Na<sub>2</sub>HPO<sub>4</sub>, 0,01M, pH 7,0 → umgestellt auf 73 MeOH  
27 Na<sub>2</sub>HPO<sub>4</sub>, 0,01M, pH 7,0

Pumpe = 200

Papier = Smaragdweiß

λ = 254nm

Range = 0.66 - ∞

Fraktion bis in H<sub>2</sub>O-Phase injiziert, 2x mit EE extrahiert, EE-Phase mit H<sub>2</sub>O gewaschen und mit MgSO<sub>4</sub> getrocknet.

Fraktionierung:


RP-Fl. 24/RP-① - Peak: Gewicht: 8mg → Ansatz Epo (117)  
426-11

RP-Fl. 24/RP-② - Peak: " = 26,7mg, NMR 004033, Wt. 132g  
" = C<sup>13</sup>  
→ V6-Filtration (5.11.96)

↓  
Epo (117)

DC-95 CH<sub>2</sub>Cl<sub>2</sub> / 5 MeOH

Yergem. 24 - RP-① = 19,3964g  
RP-② = 22,7225g





Silica gel

KG-Filtration von So 90.411 IRP-Fl. 24 IRP-②  $\hat{=}$  "C"

4-71


So 90.411 IRP-Fl. 24 IRP-② = 26,7 mg in  $\text{CH}_2\text{Cl}_2$  auflösen

offene Säule = Material = KG.60, 0,063 mm - 0,300 mm

18 cm Höhe,  $\varnothing$  0,9 cm

→ Kopf gepackt in 97  $\text{CH}_2\text{Cl}_2$  / 3  $\text{MeOH}$ , eluiert mit

97  $\text{CH}_2\text{Cl}_2$  / 3  $\text{MeOH}$



Epo C

angereichert mit Vanillin /  $\text{HNO}_3$

### Fractionierung:

So 90.411 IRP-Fl. 24 IRP-② / KG = Gewicht: 426.2

Duio-d<sub>6</sub>

NMR 004174, M. 273 g/L

NMR 004069,

M. 167 g/L

9,4 mg Takt. Masse

4,2 mg  
Epo ④

10 mg  
Epo ⑤

2 mg  
Epo ⑥

Ans ⑧

Epo C

MTC = Songline

Leergewicht: 18,4406 g

# RP-Trennung von RP-Fl. 28 / 90mg von 159mg

4-72

So 90.411 RP-Fl. 28 = 90mg von 159mg in 3 Läufen getrennt

Säule = Wundschl. (18, 25mm), 20 x 250mm

LM = 75  $\text{NaOH}$  / 25  $\text{Na}_2\text{HPO}_4$  / 0,01M, pH 7.0

Flussrate = 200

, Papier = Sinterstein

k = 25mm

, Range = 0.64 - 0


Fractionen bis zur H<sub>2</sub>O-Phase angelangt, 3x mit  
EE extrahiert, EE-Phase mit H<sub>2</sub>O gewaschen und mit  
MgSO<sub>4</sub> getrocknet.

## Fraktionierung:


So 90.411 RP-Fl. 28 / RP-① = 31,2mg, NMR 004045, Nr. 145 gelb  
13C KG = Filtration (S. 11.96)

So 90.411 RP-Fl. 28 / RP-② = 14,6mg, NMR 004046, Nr. 146 gelb  
? 428.1

So 90.411 RP-Fl. 28 / RP-③ = 17,8mg  
428.2 Reduktion





Seagewichte = RP-① = 22, 1893g  
RP-② = 17, 9169g  
RP-③ = 10, 9302g




D

RP18 chromatogr.



4-75

97. 23 июн



12

Ans. (37)

$E_{500} = MIC = 15 \mu g/ml$

England - 17, 74, 98 40-0  
4, 78, 42 40-0

Einlieferungsdatum:                     

Spektren-Nr.: 004033

# NMR-ANTRAG

GBF — Abt. Molekulare Strukturforschung

132.

4-78

Substanz-Bez.: <u>S 90.411 IRP-F124 IRP-(2) 1<sup>4</sup>C<sup>4</sup></u> Summenformel: <u>                    </u> Substanzhersteller: <u>Powlan</u> Abteilung: <u>NC (1.1.2)</u> Tel.: <u>343</u> Kernart: <u>(<sup>1</sup>H, <sup>13</sup>C, <sup>31</sup>P, andere?)</u> Substanz-Menge: <u>20</u> mg, Molmasse: <u>                    </u> geeignetes Lösungsmittel: <u>CD<sub>3</sub>OD</u> weitere Messung nach Zugabe von <u>                    </u> Substanz zurück: ja <input checked="" type="checkbox"/> nein <input type="checkbox"/>	Strukturvorschlag: <u>                    </u>          Radioaktiv <input type="checkbox"/> Toxisch <input type="checkbox"/>
--	--

## Allgemeine Angaben

Probe lagern im Kühlschrank ☒  
 im Tiefkühlfach ☐  
 im Dunkeln ☐  
 Probe auf Abruf beim Hersteller ☐

Signale erwartet zwischen  $\delta =$  0 und 9  
 Gewünscht: nur Spektrum ☒  
 plus Integral ☒  
 Interpretation ☐  
 Zahl der Akkumulationen (falls > 104):                     

## Art des Experiments

☒ <sup>1</sup>H Standardspektrum ☒  
 Entkopplung ☐ Differenz-NOE ☐  
 Differenz-Entkopplung ☐  
 Entkoppler-Frequenz(en):                     

☒ <sup>13</sup>C <sup>1</sup>H-Entkopplung:  
 Breitband ☒ selektiv ☐  
 DEPT ☒ ohne ☐

## Plot und Datenmanipulation

Gauss-Multiplikation ☐

Linienausdruck ☐

☒ <sup>1</sup>H  
 $\delta =$  8.9 bis -0.1 (0.15 ppm/cm) ☒  
 11.9 bis -0.1 (0.2 ppm/cm) ☐

Drehungen:  
 10 Hz/cm ☐ von  $\delta =$                       bis                     

☒ <sup>13</sup>C normal ( $\delta = 220$  bis 0) ☒

anderes Format:                     

Sonderwünsche: COSY ☐

<sup>13</sup>C — <sup>1</sup>H Korrel. Direkt ☐ Long-range ☐

(Nicht vom Antragsteller auszufüllen)

gemessen auf ☒ AM-300  
☐ ARX-400  
☐ DMX-600

gespeichert unter Nr. S1P24033110  
                    , 121  
                    , 120

Bitte um Rücksprache ☐

Kommentar:                     

(Unterschrift)


Epikilon

So 90.411 127-FL-24 (27-2)

426

20mg -> 26.7mg

SIPZ4033 10 1 Pohlan



Current Data Parameters  
NAME SIPZ4033  
EXPNO 10  
PROCNO 1


F2 - Acquisition Parameters  
Date\_ 13.33  
Time 13.33  
INSTRUM spect  
PROBHD 5 mm QNP 4H  
PULPROG zg30  
TD 32768  
SOLVENT MeOH  
NS 88  
DS 0  
SWH 6172.839 Hz  
FIDRES 0.188380 Hz  
AQ 2.6542580 sec  
RG 128  
DK 81.000 usec  
DE 4.50 usec  
TE 300.0 K  
D1 1.00000000 sec  
P1 15.00 usec  
DE 4.50 usec  
SF01 300.1318534 MHz  
NUC1 1H  
PL1 -4.00 dB

F2 - Processing parameters  
SI 16384  
SF 300.1299924 MHz  
WDW no  
SSB 0  
LB 0.00 Hz  
GB 0  
PC 1.00

1D NMR plot parameters  
CX 30.00 cm  
F1P 4.400 ppm  
F1 1320.57 Hz  
F2P -0.100 ppm  
F2 -30.01 Hz  
PPHCH 0.15000 ppm/cm  
HZCK 45.01950 Hz/cm

4-79

SIPZ4033 10 1 Pohlman



4-80



Epilobium

So 90.411 (RT - 71.24) RT. 2  
428

C

4-81

SIPZ4033 20 1 Pohlan

Current Data Parameters  
NAME SIPZ4033  
EXPNO 20  
PROCNO 1

F2 - Acquisition Parameters  
Date\_ 18.21  
Time 18.21


INSTRUM Spect  
PROBHD 5 mm QNP 1H  
PULPROG zgpg30  
TD 32768  
SOLVENT MeOH  
NS 8000  
DS 2

SWH 23809.523 Hz  
FIDRES 0.726609 Hz  
AQ 0.6881780 sec  
RG 2298.8  
DM 21.000 usec  
DE 5.50 usec  
TE 300.0 K

D12 0.00002800 sec  
PL13 25.00 dB  
D1 1.00000000 sec  
CPDPRG2 waltz16  
PCPD2 97.00 usec  
SF02 300.1312005 MHz  
NUC2 1H  
PL2 -3.00 dB  
PL12 12.20 dB  
P1 11.00 usec  
DE 5.50 usec  
SF01 75.4772501 MHz  
NUC1 13C  
PL1 -3.00 dB  
D11 0.03000000 sec

F2 - Processing parameters  
SI 32768  
SF 75.4676419 MHz  
WDW EM  
SSB 0  
LB 2.00 Hz  
GB 0  
PC 1.40

1D NMR plot parameters  
CX 30.00 cm  
FLP 230.000 ppm  
F1 17357.56 Hz  
F2 0.000 ppm  
F2 0.00 Hz  
PPMCM 7.66667 ppm/cm  
HZCM 578.58527 Hz/cm



DU=u, USER=chk, NAME=SIPZ4033, EXPNO=20, PROCNO=1  
 F1=230.000ppm, F2=0.000ppm, MI=0.00cm, MAXI=10000.00cm, PC=1.400

4-82

#	ADDRESS	FREQUENCY	INTENSITY
		[Hz]	[PPM]
1	6711.5	16634.855	220.4237 c 5
2	11729.2	12988.931	172.1126 c 1
3	12268.5	12597.046	166.9198 - 20
4	13677.1	11573.543	153.3577 - 18
5	15059.4	10569.200	140.0494 - 16
6	15615.0	10165.478	134.6998
7	15640.6	10146.885	134.4534
8	16421.5	9579.455	126.9346
9	16539.5	9493.721	125.7986
10	17145.8	9053.167	119.9609 - 17
11	17174.2	9032.517	119.6873
12	17398.9	8869.258	117.5240
13	17411.1	8860.417	117.4068 - 19
14	21268.5	6057.571	80.2671
15	21435.0	5936.585	78.6640
16	21589.0	5824.736	77.1819 - 15
17	22074.7	5471.765	72.5048
18	22109.2	5446.690	72.1725
19	23214.8	4643.369	61.5279
20	23911.4	4137.193	54.8208 - 4
21	24427.1	3762.479	49.8555
22	24456.7	3741.021	49.5712
23	24486.1	3719.640	49.2879
24	24515.7	3698.156	49.0032
25	24545.1	3676.791	48.7201
26	24574.6	3655.367	48.4362
27	24604.0	3633.971	48.1527
28	24859.0	3448.694	45.6977
29	24909.2	3412.198	45.2141 - 6
30	25330.7	3105.970	41.1563
31	25423.8	3038.285	40.2594 - 2
32	25580.7	2924.324	38.7494 - 8
33	25632.7	2886.490	38.2481
34	25754.1	2798.305	37.0795
35	26126.9	2527.458	33.4906
36	26224.1	2456.771	32.5540
37	26298.1	2403.007	31.8416
38	26550.6	2219.588	29.4111
39	26578.0	2199.672	29.1472
40	26607.5	2178.218	28.8629
41	26699.7	2111.245	27.9755
42	27183.1	1760.011	23.3214 - 23
43	27335.5	1649.258	21.8538
44	27438.5	1574.374	20.8616
45	27497.0	1531.911	20.2989
46	27551.4	1492.387	19.7752 - 22
47	27665.8	1409.260	18.6737 - 21
48	27735.3	1358.757	18.0045
49	27819.7	1297.446	17.1921 - 24
50	27855.0	1271.734	16.8514
51	27928.3	1218.500	16.1460
52	27951.2	1201.869	15.9256
53	27965.0	1191.823	15.7925
54	28102.7	1091.776	14.4668

4-83

Einlieferungsdatum:                     

Spektren-Nr.: 004069

**167.**

# NMR-ANTRAG

GBF — Abt. Molekulare Strukturforschung

+ - 84

Substanz-Bez.: Jo 90.4 / 1.1 RP - Fl. 24 / RP-2 / KG = C

Summenformel:                     

Substanzhersteller: Polwan

Abteilung: ML (1.1.2)

Tel.: 343

Kernart (<sup>1</sup>H, <sup>13</sup>C, <sup>31</sup>P, andere?)                     

Substanz-Menge: 17,6 mg, Molmasse:                     

geeignetes Lösungsmittel: CD<sub>3</sub>OD weitere Messung nach Zugabe von                     

Substanz zurück: ja ☒  
nein ☐

Strukturvorschlag:                     

Radioaktiv ☐

Toxisch ☐

## Allgemeine Angaben

Probe lagern im Kühlschrank ☒

im Tiefkühlfach ☐

im Dunkeln ☐

Probe auf Abruf beim Hersteller ☐

Signale erwartet zwischen  $\delta =$  0 und 9

Gewünscht: nur Spektrum ☒

plus Integral ☒

Interpretation ☐

Zahl der Akkumulationen (falls > 104):                     

## Art des Experiments

☒ <sup>1</sup>H Standardspektrum ☒

Entkopplung ☐ Differenz-NOE ☐

Differenz-Entkopplung ☐

Entkoppler-Frequenz(en):                     

☐ <sup>13</sup>C <sup>1</sup>H-Entkopplung:

Breitband ☐ selektiv ☐

DEPT ☐ ohne ☐

## Plot und Datenmanipulation

Gauss-Multiplikation ☐

☒ <sup>1</sup>H

Linienausdruck ☐

$\delta =$  8.9 bis -0.1 (0.15 ppm/cm) ☒

11.9 bis -0.1 (0.2 ppm/cm) ☐

Drehungen:

10 Hz/cm ☐ von  $\delta =$                       bis                     

☒ <sup>13</sup>C normal ( $\delta =$  220 bis 0) ☐

anderes Format:                     

Sonderwünsche: COSY ☐

<sup>13</sup>C — <sup>1</sup>H Korrel.

Direkt ☐

Long-range ☐

(Nicht vom Antragsteller auszufüllen)

gemessen auf ☒ AM-300

☐ ARX-400

☐ DMX-600

gespeichert unter Nr. SLP 24.069/10


Bitte um Rücksprache ☐

Kommentar:                     

(Unterschrift)

So 90.411 RR-FL-24 (RP-Q) 146  
 426.2  
 E110 C

SIPZ4069 10 1 Pohlan



Current Data Parameters  
 NAME SIPZ4069  
 EXPNO 10  
 PROCNO 1

F2 - Acquisition Parameters  
 Date\_ 9.18  
 Time 17.16  
 INSTRUM spect  
 PROBHD 5 mm QNP 1H  
 PULPROG zg30  
 TD 32768  
 SOLVENT MeOH  
 NS 64  
 DS 0  
 SWH 6172.839 Hz  
 FIDRES 0.188360 Hz  
 AQ 2.6542560 sec  
 RG 181  
 DW 81.000 usec  
 DE 4.50 usec  
 TE 300.0 K  
 D1 1.00000000 sec  
 P1 15.00 usec  
 DE 4.50 usec  
 SF01 300.1318534 MHz  
 NUC1 1H  
 PL1 -4.00 dB


F2 - Processing parameters  
 SI 16384  
 SF 300.1299927 MHz  
 WDW no  
 SSB 0  
 LB 0.00 Hz  
 GB 0  
 PC 1.00

1D NMR plot parameters  
 CX 30.00 cm  
 F1P 4.400 ppm  
 F1 1320.57 Hz  
 F2P -0.100 ppm  
 F2 -30.01 Hz  
 PPMCH 0.15000 ppm/cm  
 HZCH 45.01950 Hz/cm

4-85

4-86

SIPZ4069 10 1 Pohlen



Einlieferungsdatum:                     

Spektren-Nr.: 004174

**273.**

# **NMR-ANTRAG**

GBF — Abt. Molekulare Strukturforschung

4-87

Substanz-Bez.: Epo C 15090.426.2

Summenformel:                     

Substanzhersteller: Poldan

Abteilung: NC (1.1.2) Tel.: 343

Kernart: (<sup>1</sup>H, <sup>13</sup>C, <sup>31</sup>P, andere?)                     

Substanz-Menge: 5 mg, Molmasse:                     

geeignetes Lösungsmittel: DMSO-d<sub>6</sub> weitere Messung nach Zugabe von                     

Substanz zurück: ja ☒ nein ☐

Strukturvorschlag:                     

Radioaktiv ☐ Toxisch ☐

## **Allgemeine Angaben**

Probe lagern im Kühlschrank ☒  
im Tiefkühlfach ☐  
im Dunkeln ☐  
Probe auf Abruf beim Hersteller ☐

Signale erwartet zwischen  $\delta =$  0 und 9  
Gewünscht: nur Spektrum ☒  
plus Integral ☒  
Interpretation ☐  
Zahl der Akkumulationen (falls > 104):                     

## **Art des Experiments**

☒ <sup>1</sup>H Standardspektrum ☒  
Entkopplung ☐ Differenz-NOE ☐  
Differenz-Entkopplung ☐  
Entkoppler-Frequenz(en):                     

☒ <sup>13</sup>C <sup>1</sup>H-Entkopplung:  
Breitband ☒ selektiv ☐  
DEPT ☒ ohne ☐

## **Plot und Datenmanipulation**

Gauss-Multiplikation ☐

☒ <sup>1</sup>H

Linienausdruck ☐

$\delta =$  8.9 bis -0.1 (0.15 ppm/cm) ☒  
11.9 bis -0.1 (0.2 ppm/cm) ☐

Drehungen:  
10 Hz/cm ☐ von  $\delta =$                       bis                     

☒ <sup>13</sup>C normal ( $\delta = 220$  bis 0) ☒

anderes Format:                     

Sonderwünsche: COSY ☒

<sup>13</sup>C—<sup>1</sup>H Korrel. Direkt ☒ Long-range ☐

(Nicht vom Antragsteller auszufüllen)

gemessen auf ☐ AM-300  
☐ ARX-400  
☐ DMX-600

Bitte um Rücksprache ☐

Kommentar:                     


gespeichert unter Nr. 8182.4129/61  
                    , 1.11.87 SGL  
                    , 1.30.87 1H  
                    , 1.31.87 direkt  
                    , 1.32.87 13C  
                    , 1.30.87 Dept

(Unterschrift)

S030.426.2  
 Epithalon C  
 5-g  
 DMSO

4-88

SIPZ4174 10 1






4-89

DMSO

SIPZ4174 10 1



4-90

EPO, C

in DMSO


SIPZ4174 32 1

Current Data Parameters  
NAME SIPZ4174  
EXPNO 32  
PROCNO 1

F2 - Acquisition Parameters  
Date\_   
Time 3.47  
INSTRUM spect  
PROBHD 5 mm QNP 1H  
PULPROG zgpg30  
TD 32768  
SOLVENT DMSO  
NS 16000  
DS 2  
SWH 23809.523 Hz  
FIDRES 0.726609 Hz  
AQ 0.6881780 sec  
RG 1024  
DM 21.000 usec  
DE 5.50 usec  
TE 300.0 K  
D12 0.0002000 sec  
PL13 25.00 dB  
D1 1.0000000 sec  
CPDPRG2 waltz16  
PCPD2 97.00 usec  
NUC2 1312005 MHz  
LH 13  
PL2 -3.00 dB  
PL12 12.20 dB  
P1 11.00 usec  
DE 5.50 usec  
SFO1 75.4772501 MHz  
NUC1 13C  
PL1 -3.00 dB  
D11 0.0300000 sec

F2 - Processing parameters  
SI 32768  
SF 75.4677834 MHz  
WDW EM  
SSB 0  
LB 2.00 Hz  
GB 0  
PC 1.40

1D NMR plot parameters  
CX 30.00 cm  
FLP 220.000 ppm  
F1 16602.91 Hz  
F2P 0.000 ppm  
F2 0.00 Hz  
PPMCH 7.33333 ppm/cm  
HZCH 553.43042 Hz/cm



/u/data/chk/nmr/SIPZ4174/32/pdata/1/screen

Fri 11:53:32

DU=u, USER=chk, NAME=SIPZ4174, EXPNO=32, PROCNO=1  
 F1=220.000ppm, F2=0.000ppm, MI=0.00cm, MAXI=10000.00cm, PC=1.400

#	ADDRESS	FREQUENCY [Hz]	INTENSITY [PPM]
1	6816.6	16416.988	217.5364°
2	11741.3	12838.661	170.1211°
3	12351.2	12395.488	164.2487°
4	13595.5	11491.410	152.2691°
5	15141.4	10368.112	137.3846°
6	15584.0	10046.493	133.1229°
7	16473.1	9400.500	124.5631°
8	17077.2	8961.576	118.7470°
9	17205.5	8868.332	117.5115°
10	21259.5	5922.664	78.4794°
11	21602.4	5673.510	75.1779°
12	22047.1	5350.393	70.8964°
13	23883.2	4016.285	53.2185°
14	24807.3	3344.793	44.3208°
15	25221.4	3043.879	40.3335°
16	25250.1	3023.054	40.0575°
17	25278.9	3002.137	39.7804°
18	25307.6	2981.238	39.5034°
19	25336.4	2960.312	39.2262°
20	25365.2	2939.395	38.9490°
21	25394.1	2918.407	38.6709°
22	25613.2	2759.223	36.5616°
23	26182.3	2345.693	31.0820°
24	26292.8	2265.400	30.0181°
25	26542.5	2083.948	27.6137°
26	27070.1	1700.635	22.5346°
27	27306.2	1529.078	20.2613°
28	27448.7	1425.549	18.8895°
29	27601.3	1314.677	17.4204°
30	27743.0	1211.706	16.0559°
31	27885.2	1108.357	14.6865°

78.5  
75.2  
70.9

9 + 11

EPO C

in DMSO

Einlieferungsdatum:                     

Spektren-Nr.:                     

004045

4-92

145.

# NMR-ANTRAG

GBF — Abt. Molekulare Strukturforschung

Substanz-Bez.: So 90. 411 RP-Fl. 2P 1RP-①  
 Summenformel:                     (?)  
 Substanzhersteller: Roldan  
 Abteilung: NC (1.1-2) Tel.: 343  
 Kernart: <sup>1</sup>H <sup>13</sup>C <sup>31</sup>P, andere?  
 Substanz-Menge: 10.5 mg, Molmasse:                       
 geeignetes Lösungsmittel: CD<sub>3</sub>OD weitere Messung nach Zugabe von                       
 Substanz zurück: ja ☒ nein ☐

Strukturvorschlag:                     

Radioaktiv ☐

Toxisch ☐

## Allgemeine Angaben

Probe lagern ☒ im Kühlschrank  
☐ im Tiefkühlfach  
☐ im Dunkeln  
 Probe auf Abruf beim Hersteller ☐

Signale erwartet zwischen

$\delta$  = 0 und 9

Gewünscht: nur Spektrum ☒

plus Integral ☒

Interpretation ☐

Zahl der Akkumulationen (falls > 104):                     

## Art des Experiments

☒ <sup>1</sup>H Standardspektrum ☒  
 Entkopplung ☐ Differenz-NOE ☐  
 Differenz-Entkopplung ☐  
 Entkoppler-Frequenz(en):                     

☒ <sup>13</sup>C <sup>1</sup>H-Entkopplung:

Breitband ☒ selektiv ☐

DEPT ☒ ohne ☐

## Plot und Datenmanipulation

Gauss-Multiplikation ☐

☒ <sup>1</sup>H

Linienausdruck ☐

$\delta$  = 8.9 bis -0.1 (0.15 ppm/cm) ☒

11.9 bis -0.1 (0.2 ppm/cm) ☐

Drehungen:

10 Hz/cm ☐ von  $\delta$  =                      bis                     

☒ <sup>13</sup>C normal ( $\delta$  = 220 bis 0) ☒

anderes Format:                     

Sonderwünsche: COSY ☐

<sup>13</sup>C — <sup>1</sup>H Korrel.

Direkt ☐

Long-range ☐

(Nicht vom Antragsteller auszufüllen)

gemessen auf ☒ AM-300  
☐ ARX-400  
☐ DMX-600

gespeichert unter Nr. 51P24045/10

120  
121

Bitte um Rücksprache ☐

Kommentar:                     

(Unterschrift)


Epothelone

5050411 RP-FL-28 RP-1

428

195mg → 31.2g

SIPZ4045 10 1 Pohlman



Current Data Parameters  
NAME SIPZ4045  
EXPNO 10  
PROCNO 1

F2 - Acquisition Parameters  
Date\_ 10.13

Time 10.13

INSTRUM spect

PROBHD 5 mm QNP 1H

PULPROG zg30

TD 32768

SOLVENT MeOH

NS 88

DS 0

SMH 6172.839 Hz

FIDRES 0.188380 Hz

AQ 2.6542580 sec

RG 128

DK 81.000 usec

TE 4.50 usec

300.0 K

1.00000000 sec

P1 15.00 usec

DE 4.50 usec

SFO1 300.1318534 MHz

NUC1 1H

PL1 -4.00 dB

F2 - Processing parameters

SI 16394

SF 300.1299824 MHz

WDW no

SSB 0

LB 0.00 Hz

GB 0

PC 1.00

1D NMR plot parameters

CX 30.00 cm

F1P 4.400 ppm

F1 1320.57 Hz

F2P -0.100 ppm

F2 -30.01 Hz


PPMCK 0.15000 ppm/cm

HZCK 45.01950 Hz/cm

4-93

4-94

SIPZ4045 10 1 Pohlman



5090.411 RP-F28/RP-11

428

1

SIPZ4045 20 1 Pohlman

19.5

4-95

Current Data Parameters  
NAME SIPZ4045  
EXPNO 20  
PROCNO 1

F2 - Acquisition Parameters


Date 11.15  
Time 11.15  
INSTRUM spect  
PROBHD 5 mm QNP 1H  
PULPROG zgpg30  
TO 32768  
SOLVENT MeOH  
NS 3760  
DS 2  
SHH 23609.523 Hz  
FIDRES 0.726609 Hz  
AQ 0.6881780 sec  
RG 2048  
DK 21.000 usec  
DE 4.50 usec  
TE 300.0 K  
PL13 0.00002000 sec  
PL12 25.00 dB  
D1 0.80000001 sec  
PCPD2 waitz16  
SF02 300.1312005 MHz  
NUC2 1H  
PL2 -3.00 dB  
PL12 12.20 dB  
P1 11.00 usec  
DE 4.50 usec  
SF01 75.4772501 MHz  
NUC1 13C  
PL1 -3.00 dB  
D11 0.03000000 sec

F2 - Processing parameters

SI 32768  
SF 75.4675434 MHz  
WDW EM  
SSB 0  
LB 2.00 Hz  
GB 0  
PC 1.00

1D NMR plot parameters

CX 30.00 cm  
F1 230.000 ppm  
F2 17357.56 Hz  
F2 0.000 ppm  
F2 0.00 Hz  
PP4CM 7.6667 ppm/cm  
HZCM 578.56627 Hz/cm



DU=u, USER=chk, NAME=SIPZ4045, EXPNO=20, PROCNO=1  
 F1=285.042ppm, F2=-30.451ppm, MI=0.00cm, MAXI=10000.00cm, PC=1.000

4-96

#	ADDRESS	FREQUENCY	INTENSITY
		[Hz]	[PPM]
1	6727.1	16623.523	220.2735
2	11639.8	13053.866	172.9730
3	11728.8	12989.188	172.1160
4	12268.2	12597.304	166.9232
5	13675.7	11574.554	153.3711
6	15038.5	10584.390	140.2507
7	15046.4	10578.604	140.1740
8	15482.8	10261.529	135.9726
9	15633.7	10151.841	134.5191
10	15860.8	9986.855	132.3329
11	16164.2	9766.376	129.4114
12	16979.5	9174.034	121.5625
13	17026.3	9139.974	121.1112
14	17142.8	9055.376	119.9902
15	17413.0	8859.035	117.3885
16	19273.8	7506.939	99.4723
17	21223.3	6090.396	80.7021
18	21426.0	5943.124	78.7506
19	21563.9	5842.914	77.4228
20	21672.9	5763.763	76.3740
21	21702.0	5742.604	76.0936
22	21782.6	5684.032	75.3175
23	22079.0	5468.641	72.4634
24	22226.1	5361.791	71.0475
25	22380.5	5249.614	69.5611
26	23214.8	4643.351	61.5277
27	23532.2	4412.754	58.4721
28	23849.3	4182.324	55.4188
29	23899.7	4145.718	54.9337
30	24427.1	3762.483	49.8556
31	24456.7	3741.021	49.5712
32	24486.1	3719.649	49.2880
33	24515.7	3698.168	49.0034
34	24545.1	3676.783	48.7200
35	24574.6	3655.310	48.4355
36	24604.1	3633.926	48.1521
37	24831.9	3468.352	45.9581
38	24886.7	3428.535	45.4305
39	25410.1	3048.276	40.3918
40	25435.5	3029.783	40.1468
41	25569.3	2932.605	38.8591
42	25601.9	2908.881	38.5447
43	25693.4	2842.427	37.6642
44	25729.8	2815.970	37.3136
45	25834.2	2740.097	36.3082
46	25910.7	2684.506	35.5716
47	25932.5	2668.660	35.3616
48	26078.1	2562.849	33.9596
49	26144.2	2514.872	33.3238
50	26188.5	2482.672	32.8972
51	26319.9	2387.194	31.6320
52	26484.0	2267.948	30.0519
53	26750.9	2074.054	27.4827
54	26942.3	1934.968	25.6397
55	27169.2	1770.080	23.4548
56	27172.5	1767.651	23.4226
57	27303.7	1672.336	22.1596
58	27438.6	1574.322	20.8609



59	27582.0	1470.167	19.4808	6.62
60	27666.0	1409.089	18.6714	5.77
61	27735.5	1358.635	18.0029	0.80
62	27784.2	1323.231	17.5338	6.44
63	27929.8	1217.439	16.1319	6.08
64	27959.4	1195.924	15.8468	2.67
65	27973.0	1186.054	15.7161	7.35
66	28015.2	1155.338	15.3090	1.07
67	28070.7	1115.056	14.7753	1.22
68	28102.7	1091.755	14.4665	4.71
69	28199.8	1021.210	13.5318	1.65
70	28503.1	800.822	10.6115	0.73
71	28571.0	751.523	9.9582	1.60
72	28776.3	602.322	7.9812	0.91
73	29178.2	310.304	4.1118	1.55

4-97

Einlieferungsdatum:                     

Spektren-Nr.: 004070

**168.**

**NMR-ANTRAG**  
GBF — Abt. Molekulare Strukturforschung

4-98

<p>Substanz-Bez.: <u>Jo 90.411/PP-Fl.28PP-1/UG-1</u></p> <p>Summenformel: <u>= 1</u></p> <p>Substanzhersteller: <u>Pohlan</u></p> <p>Abteilung: <u>NC (1.1.2)</u> Tel.: <u>343</u></p> <p>Kernart: (<sup>1</sup>H, <sup>13</sup>C, <sup>31</sup>P, andere?) <u>                    </u></p> <p>Substanz-Menge: <u>11,2</u> mg, Molmasse: <u>                    </u></p> <p>geeignetes Lösungsmittel: <u>CD<sub>3</sub>O</u> weitere Messung nach Zugabe von <u>                    </u></p> <p>Substanz zurück: ja <input checked="" type="checkbox"/> nein <input type="checkbox"/></p>	<p>Strukturvorschlag: <u>                    </u></p>     <p>Radioaktiv <input type="checkbox"/> Toxisch <input type="checkbox"/></p>
---	--

**Allgemeine Angaben**

Probe lagern ☒ im Kühlschrank ☐ im Tiefkühlfach ☐ im Dunkeln ☐

Probe auf Abruf beim Hersteller ☐

Signale erwartet zwischen  $\delta =$  0 und 9

Gewünscht: nur Spektrum ☒ plus Integral ☒ Interpretation ☐

Zahl der Akkumulationen (falls > 104):                     

**Art des Experiments**

☒ <sup>1</sup>H Standardspektrum ☐ Differenz-NOE ☐

Entkopplung ☐ Differenz-Entkopplung ☐

Entkoppler-Frequenz(en):                     

☐ <sup>13</sup>C <sup>1</sup>H-Entkopplung:

Breitband ☐ selektiv ☐

DEPT ☐ ohne ☐

**Plot und Datenmanipulation**

Gauss-Multiplikation ☐

☒ <sup>1</sup>H  $\delta =$  8.9 bis -0.1 (0.15 ppm/cm) ☒ 11.9 bis -0.1 (0.2 ppm/cm) ☐

Linienausdruck ☐

Drehungen: 10 Hz/cm ☐ von  $\delta =$                       bis                     

☒ <sup>13</sup>C normal ( $\delta =$  220 bis 0) ☐ anderes Format:                     

Sonderwünsche: COSY ☐ <sup>13</sup>C — <sup>1</sup>H Korrel. Direkt ☐ Long-range ☐

(Nicht vom Antragsteller auszufüllen)

gemessen auf ☒ AM-300 ☐ ARX-400 ☐ DMX-600

gespeichert unter Nr. S182407010

Bitte um Rücksprache ☐


Kommentar:                     

(Unterschrift)

SIPZ4070 10 1 Pohlen

So 90.411 RR 77.281 RR 10 (K6)  
428.3  
11,2 mg

Epo D




Current Data Parameters  
NAME SIPZ4070  
EXPNO 10  
PROCNO 1

F2 - Acquisition Parameters  
Date\_ 9.26  
Time 9.26  
INSTRUM spect  
PROBHD 5 mm GNP 1H  
PULPROG zg30  
TD 32768  
SOLVENT MeOH  
NS 64  
DS 0  
SWH 6172.839 Hz  
FIDRES 0.186360 Hz  
AQ 2.6542560 sec  
RG 161  
DN 81.000 usec  
DE 4.50 usec  
TE 300.0 K  
D1 1.00000000 sec  
P1 15.00 usec  
DE 4.50 usec  
SF01 300.1318534 MHz  
NUC1 1H  
PL1 -4.00 dB

F2 - Processing parameters  
SI 16384  
SF 300.1299927 MHz  
WDW no  
SSB 0  
LB 0.00 Hz  
GB 0  
PC 1.00

1D NMR plot parameters  
CX 30.00 cm  
F1P 4.400 ppm  
F1 1320.57 Hz  
F2P -0.100 ppm  
F2 -30.01 Hz  
PPMCM 0.15000 ppm/cm  
HZCM 45.01950 Hz/cm

SIPZ4070 10 1 Pohlman



4 - 100

Einlieferungsdatum: 22.05.2023

Spektren-Nr.: 004175

# NMR-ANTRAG

GBF — Abt. Molekulare Strukturforschung

274.

4 - 107

Substanz-Bez.: Epo.D1 So 90. 42p.3  
 Summenformel: \_\_\_\_\_  
 Substanzhersteller: Molvan  
 Abteilung: NC (1.1.2) Tel.: 363  
 Kernart (<sup>1</sup>H, <sup>13</sup>C, <sup>31</sup>P, andere?) \_\_\_\_\_  
 Substanz-Menge: 11,2 mg, Molmasse: \_\_\_\_\_  
 geeignetes Lösungsmittel: DMSO-d<sub>6</sub> weitere Messung nach Zugabe von \_\_\_\_\_  
 Substanz zurück: ja ☒ nein ☐  
 Radioaktiv ☐ Toxisch ☐

**Allgemeine Angaben**  
 Probe lagern im Kühlschrank ☒  
 im Tiefkühlfach ☐  
 im Dunkeln ☐  
 Probe auf Abruf beim Hersteller ☐  
 Signale erwartet zwischen  $\delta =$  0 und 9  
 Gewünscht: nur Spektrum ☒  
 plus Integral ☒  
 Interpretation ☐  
 Zahl der Akkumulationen (falls > 104): \_\_\_\_\_

**Art des Experiments**  
☒ <sup>1</sup>H Standardspektrum ☒  
 Entkopplung ☐ Differenz-NOE ☐  
 Differenz-Entkopplung ☐  
 Entkoppler-Frequenz(en): \_\_\_\_\_  
☒ <sup>13</sup>C <sup>1</sup>H-Entkopplung:  
 Breitband ☒ selektiv ☐  
 DEPT ☒ ohne ☐

**Plot und Datenmanipulation**  
 Gauss-Multiplikation ☐  
☒ <sup>1</sup>H Linienausdruck ☐  
 $\delta =$  8.9 bis -0.1 (0.15 ppm/cm) ☒  
 11.9 bis -0.1 (0.2 ppm/cm) ☐  
 Drehungen: 10 Hz/cm ☐ von  $\delta =$  \_\_\_\_\_ bis \_\_\_\_\_  
☒ <sup>13</sup>C normal ( $\delta = 220$  bis 0) ☒ anderes Format: \_\_\_\_\_

Sonderwünsche: COSY ☒ <sup>13</sup>C-<sup>1</sup>H Korrel. Direkt ☒ Long-range ☐

gemessen auf ☐ AM-300  
☐ ARX-400  
☐ DMX-600

Bitte um Rücksprache ☐

Kommentar:

(Nicht vom Antragsteller auszufüllen)

gespeichert unter Nr.


SIR 2, 475/6 ✓  
SIR 2, 475/6 ✓  
—, — ✓  
— 12 ✓  
— 20 ✓  
— 21 ✓  
— 130 ✓  
  
(Unterschrift) 21 die

(Unterschrift)



4-103

M



EPO 2

4-104


SIPR4175 20 1

Current Data Parameters  
NAME SIPR4175  
EXPNO 20  
PROCNO 1

F2 - Acquisition Parameters  
Date\_   
Time 10.47  
INSTRUM art400  
PROBHD 5 mm QNP 1H  
PULPROG zgpg30  
TD 32768  
SOLVENT DMSO  
NS 2800  
DS 2  
SWH 33333.332 Hz  
FIDRES 1.017253 Hz  
AQ 0.4915700 sec  
RG 16384  
DM 15.000 usec  
DE 21.43 usec  
TE 300.2 K  
DL2 0.0002000 sec  
DL6 22.00 dB  
DI 1.00000000 sec  
CPDPRG waltz16  
F31 75.00 usec  
DL1 0.03000000 sec  
DL5 19.00 dB  
P1 10.00 usec  
DE 21.43 usec  
SFOL 100.6254358 MHz  
NUCLEUS 13C

F2 - Processing parameters  
SI 32768  
SF 100.6128210 MHz  
WDW EM  
SSB 0  
LB 2.00 Hz  
GB 0  
PC 1.00

1D NMR plot parameters  
CX 30.00 cm  
F1P 220.000 ppm  
F1 22134.82 Hz  
F2P 0.000 ppm  
F2 0.00 Hz  
FPCMH 7.33333 ppm/cm  
B2CH 737.82739 Hz/cm





EPO "D"

/u/data/beate/nmr/SIPR4175/20/pdata/1/screen

Thu 11:00:25

DU=u, USER=beate, NAME=SIPR4175, EXPNO=20, PROCNO=1  
 F1=220.000ppm, F2=0.000ppm, MI=0.00cm, MAXI=10000.00cm, PC=1.400

4-105

#	ADDRESS	FREQUENCY [Hz]	INTENSITY [PPM]
1	7281.0	21874.760	217.4152 - 5
2	7869.0	21276.707	211.4711
3	7890.0	21255.285	211.2582
4	7910.4	21234.568	211.0523
5	8761.3	20368.973	202.4491
6	11965.3	17109.715	170.0550 - 1
7	12546.3	16518.658	164.1805 - 20
8	13727.1	15317.464	152.2417 - 18
9	15100.4	13920.527	138.3574 - 13
10	15178.4	13841.204	137.5690 - 16
11	16887.6	12102.428	120.2871 - 13
12	17037.1	11950.432	118.7764 - 17
13	17170.7	11814.513	117.4255 - 19
14	20978.3	7941.236	78.9287
15	21325.2	7588.323	75.4210 - 7
16	21783.3	7122.308	70.7893 - 3
17	23516.1	5359.594	53.2695 - 4
18	23983.1	4884.557	48.5481
19	24390.8	4469.850	44.4262 - 6
20	24816.1	4037.134	40.1254
21	24836.8	4016.094	39.9163
22	24857.5	3995.087	39.7075
23	24878.1	3974.093	39.4989
24	24898.8	3953.056	39.2898
25	24919.4	3932.073	39.0812 - 2
26	24939.8	3911.324	38.8750
27	25192.9	3653.922	36.3167 - 8
28	25633.0	3206.134	31.8661 - 11
29	25662.2	3176.465	31.5712 - 14
30	25827.4	3008.435	29.9011 - 9?
31	26223.3	2605.679	25.8981 - 10?
32	26515.7	2308.261	22.9420 - 26
33	26556.8	2266.417	22.5261 - 23
34	26819.3	1999.370	19.8719 - 22
35	26921.3	1895.664	18.8412 - 21
36	26945.8	1870.694	18.5930
37	27036.9	1778.058	17.6723 - 25
38	27182.0	1630.464	16.2053 - 24
39	27348.0	1461.631	14.5273 - 27

\* können vertauscht sein

**Reply to the**  
**Opposition Statement against EP-B-1186606**  
concerning the identity of epothilones C and D

by Gerhard Höfle  
GBF, September 8, 2005

Contributions by Dr. K. Gerth, H. Steinmetz (GBF),  
Prof. D. Schinzer (University of Magdeburg)

1. Introduction.....	p. 2
2. NMR spectroscopy.....	p. 3
3. Production with strain DSM 6773.....	p. 4
4. Statement of Prof. Schinzer.....	p. 6
5. References.....	p. 6
6. Tables.....	p. 7
7. Appendix.....	pdf-file

## Introduction

In 1990/1991 epothilones A and B have been produced for the first time on the gram-scale with *Sorangium cellulosum* So ce90 wild strain. A patent was filed for epothilones A and B November 19, 1991, and the strain was deposited at the German Strain Collection DSMZ under the code DSM 6773. At that time during large scale isolation work more lipophilic epothilones were observed during RP chromatography of epothilones A and B<sup>1)</sup> which however were not isolated because of the small amounts present and, later lack of interest in epothilones.

After the tubulin activity was published by Bollag et al. in 1995 work was resumed, and a number of other *Sorangium cellulosum* strains were identified to produce epothilones A and B. From these strain So cel198 was selected for further work because of low abundance of other unwanted metabolites which otherwise interfere with isolation. As side products from this strain epothilones C and D were isolated in pure state and the structures elucidated in June 1996 as documented already. In September/October a 350 L fermenter with strain So cel198 was run for the production of epothilones A and B. From this as side products several hundred milligrams of epothilones C and D were isolated. With this material a complete set of NMR spectra in DMSO-D<sub>6</sub> was obtained, and the chromatographic behavior determined as basis for the patent application November 18, 1996 which represents the first description of epothilones C and D. Later epothilones C and D were obtained by total synthesis by the Danishefsky, Nicolaou and Schinzer groups. Interestingly, in papers by Danishefsky et al. they were named desoxyepothilones A and B<sup>2)</sup>. From natural sources epothilones C and D were re-isolated using So ce90/B2, a mutant with improved epothilone A and B production (Hardt et al.<sup>3)</sup>), and P450 knock-out mutants obtained by UV-irradiation (Gerth et al.<sup>4)</sup>) or genetic engineering (Lau et al.<sup>5)</sup>). To demonstrate feasibility of epothilone A and B total synthesis via C and D, epoxidation using dimethyl dioxirane and *m*-chloroperbenzoic acid were performed in June 1996 and in more detail in November/ December 1996 (Höfle et al.<sup>6)</sup>).

In the present Opposition Statement of Sloan Kettering Institute for Cancer Research it is claimed that the compounds described in the GBF patent EP-B-1186606<sup>(2)</sup> are not epothilones C and D but rather other epothilones of non-defined structure. This conclusion is based on two observations:

1. Certain signals in the proton and carbon NMR spectra taken from MSKCC epothilones C and D differ significantly from those given in EP-B-1186606 (Prof. J. D. Roberts),
2. Attempts to obtain epothilones C and D by cultivation of *Sorangium cellulosum* strain So ce90 obtained as DSM6773 from the DSMZ failed (Dr. P. J. Licari, KOSAN).

In the following it is clearly proven that the compounds isolated in 1996 had the structures claimed in EP-B-1186606, today known as epothilones C and D.

It is further demonstrated that by re-fermentation of strain DSM6773 and isolation as described indeed epothilones C and D are obtained.

Prof. Schinzer confirms that epothilones C and D obtained in 1996 from GBF were identical with his synthetic compounds.

1a) EP 2 07 129 322 7

Wild So ce 90 = DSM 6773

## NMR spectroscopy of epothilones C and D

From the first proton NMR spectra recorded in June 1996 (as documented before) and biosynthetic considerations the structures of epothilones C and D were unequivocally derived. When larger amounts of the compounds became available in Oct./Nov. complete sets of 1D and 2D spectra were recorded in DMSO- $D_6$ . The particular solvent was chosen to allow observation of hydroxy proton signals and couplings, and to facilitate comparison with the published data for epothilones A and B (Höfle et al.<sup>7)</sup>). In the appendix NMR Request Forms, shift records and 1D proton and carbon spectra of epothilones C and D are given. In Tab. 1 and 2 carbon shifts, in Tab. 3 and 4 proton shifts determined November 14 and 15, 1996,<sup>1)</sup> are compared with those from recent samples measured May 18, 2005 at GBF, and KOSAN (Opposition Statement, p. 54-60).

### <sup>13</sup>C shifts for epothilone C (Tab. 1)

The shifts for all carbons except C1 and C2 are found within  $\pm 0.1$  ppm. C2 deviates by 0.2 ppm due to partial overlap with solvent signals, whereas carbonyl C1 is by 0.5 ppm too high in the KOSAN spectrum. This may be attributed to a solvent induced shift, which is common with carbonyl carbons.

### <sup>13</sup>C shifts for epothilone D (Tab. 2)

The shifts for all carbons except C2 are in excellent agreement within  $\pm 0.1$  ppm. C2 deviates by 0.3-0.4 ppm due to overlap with solvent signals. The values reported by KOSAN are consistently too high by 0.2 ppm which is attributed to an offset of the reference.

### <sup>1</sup>H shifts for epothilone C (Tab. 3)

In general most proton NMR signals of complex natural products are complex because of multiple couplings and signal overlap. Under these circumstances the shift differences of  $\pm 0.03$  ppm between GBF, Nov. 96 and KOSAN measurements indicate excellent Übereinstimmung.

### <sup>1</sup>H shifts for epothilone D (Tab. 4)

The majority of shifts are identical for GBF, Nov. 96 and KOSAN measurements, and only few deviate up to  $\pm 0.3$  ppm.

**The above comparison of chemical shift data unequivocally prove that the epothilones isolated in 1996 were indeed epothilones C and D.**

How can the deviating values in the table of EP-B-1186606 be explained although they were extracted from the spectra measured November 14/15, 1996 which contain the correct ones ?

As basis for writing the table for the patent application Mr. Steinmetz used an existing table with the data for epothilones A and B and replaced the values atom for atom with the corresponding values from the epothilone C and D spectra. Obviously, he started with the signals around the 12,13 double bond and epoxide, respectively, and others differing significantly in the olefin and epoxide series. These values are marked in blue in Tab. 1-4. Later, he apparently forgot to adjust also the slightly deviating values. They are marked in red like those for epothilones A and B in Tab. 1-4. When I checked the table fabricated by Mr. Steinmetz for plausibility before submission of the patent application I had no chance to discover the small differences.

(ii) after PLOS DE 195 42 986 17.11.1995

Coming back to the claim on p. 39 of the Opposition Statement that the compounds described in the patent were actually isomers of epothilones C and D with the same molecular mass, viz.  $m/z$  477 and 491, respectively. This can be ruled out by comparison of the  $^{13}\text{C}$  shifts for e.g. the known epothilone C isomers, 12*E*-epothilone C, epothilones D1 and D2. As shown in Tab. 5 shift differences of 2 up to 6.6 ppm are observed which is by far above the (slightly wrong) values in the patent.

Even though the values given in EP-B-1186606 deviate slightly more than is generally observed as experimental error, they are not misleading in a structural assignment. To demonstrate the variability of chemical shifts for complex natural products published data for epothilone C in  $\text{CDCl}_3$  from different authors are summarized in Tab. 6.

### Production of epothilones C and D with *Sorangium cellulosum* DSM 6773

Epothilones C and D are the primary products of epothilone biosynthesis. After release from the polyketide synthase complex they are modified by so-called "decorating enzymes" to the epoxides, epothilone A and B<sup>4,8,9)</sup>, and then to the 21-hydroxy derivatives, epothilones E and F<sup>10)</sup>.

Thus any *Sorangium* strain capable of epothilone A and/or B synthesis has to produce as intermediates epothilones C and/or D. (Molnar et al.<sup>8)</sup> and Tang et al.<sup>9)</sup>). Whether these intermediates can be observed and isolated from a culture depends on a variety of preconditions which are not well defined and mostly unknown. Certainly the harvest time, media composition, export activity of the organism and presence of XAD adsorber resin are essential factors. It is not surprising for the expert that in a single run these compounds may be missed as they are minor side products with wild strains or mutants generated for production of epothilones A and B. Only P450 knock-out mutants produce reliably high amounts of epothilones C and D (Gerth et al.<sup>4)</sup>, Lau et al.<sup>5)</sup>).

KOSAN ordered strain DSM6673 from DSMZ three times, November 26, 1999, March 9, 2000 and May 18, 2004. From this and data in the Opposition Statement (Appendix 4, p.1-2) it follows that the second shipment of March 9, 2000 was used for the attempt to reproduce the production of epothilones C and D. No information is given whether and how the strain was preserved or kept in culture for several years until the experiments were performed in August-November 2004. It is well known that myxobacteria like other microorganisms change their properties during extended cultivation due to clonal selection or unfavourable conditions for preservation. Thus without an analytical check for epothilone production on the shake flask level it was high-risk to start with a 70 L batch. Even though the procedure in the patent could be reproduced yielding 167 g of crude extract (180 g in the patent). This material was separated on LH-20 under supposedly the same conditions as described in the patent. The fraction eluting between 240-300 minutes was collected without checking the presence of epothilones by TLC or HPLC. The fraction contained only 35 mg instead of 72 g in the patent. From the fact that epothilones A, B, C and D co-elute from Sephadex LH-20 it must be concluded at this point that the right fraction was missed. This is however not too surprising, as it is common experience that retention times are very sensitive to a number of parameters which, particularly in large scale chromatography, cannot be controlled. The expert in the field in such case certainly would have checked adjoining fractions for the presence of epothilones before discarding them and not continued tedious work on a tiny fraction (2,250-times less than expected). In addition, this fraction was not analysed for the presence of epothilones C and D but stupidly processed further without a result. Obviously, this was the actual purpose of the exercise.

It is important to notice at this point that KOSAN reported the isolation of epothilone C and D from a not specified *Sorangium cellulosum* (wild) strain: "They are secreted as minor products during the fermentation process with a combined yield of about 0.4 mg/L" (Lau et al.<sup>5</sup>).

To demonstrate that wild strain DSM 6673 indeed produces epothilones C and D it was newly ordered from DSMZ, and obtained May 24, 2005. The culture (now coded as So ce90wild DSM 6773) on slant agar was propagated on agar plates and taken into liquid culture as described by Gerth et al.<sup>11</sup> and in epothilone A/B patents.<sup>12/13</sup>

In detail,

1. agar plates with probion medium<sup>13</sup> were inoculated on May 24, and propagated,
2. H medium<sup>12</sup> plus 1.2% HEPES buffer (500 mL) was inoculated on June 16, and the culture propagated
3. 22 shaking flasks with H medium<sup>12</sup> plus 1.2% HEPES buffer (550 mL each) were inoculated on July 26,
4. a 150 L fermentor with H medium<sup>12</sup> (100 L) and 2 Kg of wet XAD-16 adsorber was inoculated with 10 l of the above culture on August 1 (pH adjustment with 10% aq. acetic acid, and 10% aq. KOH, 32°C, 30% oxygen saturation, see also Figure 1),
5. the adsorber resin was harvested by sieving on August 15 and immediately processed further as described below.

When the production of epothilones C and D was determined on the shake flask level a constantly high proportion of spirangiens was observed and only very little of epothilones A-D. This unfavourable production profile may be due to the short time of adaption of the strain to the liquid medium. It was later reproduced with the production fermentor containing 2.4 g of spirangiens A and C, and only 3.1 mg of epothilone A, 1.8 mg of epothilone B, 1.4 mg of epothilone C, and 0.5 mg of epothilone D (Figure 3 - 7). To facilitate the isolation of such small amounts of epothilones in presence of co-eluting spirangiens an additional extraction step with sodium carbonate solution was introduced which removed most of the spirangiens as carboxylic acid salts.

The entire isolation process from wet XAD adsorber resin to pure epothilones C and D is given in Figures 2a and 2b. It should be noted that the presence of epothilones in LH20 and RP-silica gel chromatography fractions was monitored by HPLC/MS. Thus no loss of material occurred, and the expected amounts of 1.4 mg of epothilone C and 0.5 mg of epothilone D were obtained in pure state. From physical data, in particular proton and carbon NMR spectra, the identity of the compounds is equivocally proven (Table 7).

Thus, So ce90wild DSM 6773 (patent strain of DE 4138042) is indeed producing Epothilones C and D.

(1) DE 4138042

### Statement from Prof. Schinzer

In 1996 Dieter Schinzer was Professor for Organic Chemistry at the University of Braunschweig and a colleague of mine. Like other synthetic chemists he obtained the absolute configuration of epothilone A and B around November 1995. He developed plans for a total synthesis and discussed certain crucial steps with me. In summer 1996 I mentioned to him the isolation of epothilones A and B and my preliminary experiments on the epoxidation to give preferably the desired stereoisomer. In October he received samples of ca. 5 mg each for comparison purposes. Both were found to be identical with his compounds from total synthesis. This was acknowledged for epothilone C in a paper on epothilone A total synthesis.<sup>14</sup>

From my recent contacts with Prof. Schinzer I know that he is willing to witness this.

### References

- 1) H. Steinmetz, unpublished.
- 2) D.-S. Su et al. *Angew. Chem. Int. Ed.*, 36, 757, 1997.
- 3) I. H. Hardt et al., *J. Nat. Prod.* 2001, 64, 847.
- 4) Gerth, K. et al., *J. Antibiot.*, 54, 144, 2001.
- 5) J. Lau et al., *Biotech & Bioengin.* 78, 280, 2002.
- 6) G. Höfle et al. *Pure Appl. Chem.* 71, 41, 2002.
- 7) Höfle, G. et al *Angew. Chem. Int. Ed.*, 35, 1567, 1996.
- 8) I. Molnar et al. *Chemistry & Biology* 7, 97, 2000.
- 9) L. Tang et al., *Science*, 2000, 287, 640.
- 10) Gerth, K. et al., *J. Antibiot.*, 55, 41, 2002.
- 11) Gerth, K. et al., *J. Antibiot.*, 49, 560, 1996.
- 12) H medium is the production medium used in DE 4138042 (Nov. 19,1991).
- 13) Pradella et al. *Arch Microbiol*, 178, 484, 2002.
- 14) D. Schinzer et al., *Chem. Eur. J.* 5, 2483, 1999

Tab. 1  $^{13}\text{C}$ -NMR chemical shifts of epothilone C in DMSO- $\text{D}_6$ 

C-Atom	Epo A <sup>1</sup>	EP-B-1186606 18.11.96	GBF <sup>3</sup> 15.11.96	GBF <sup>4</sup> 18.5.05	Kosan <sup>5</sup>
1	170.3	170.3	170.1	170.0	170.6
2	38.4	38.4	38.7	38.9	38.8
3	71.2	71.2	70.9	70.8	70.8
4	53.1	53.1	53.2	53.2	53.2
5	217.1	217.1	217.5	217.5	217.5
6	45.4	45.4	44.3	44.2	44.3
7	75.9	75.9	75.2	75.1	75.1
8	35.4	35.4	36.6	36.6	36.6
9	29.6	27.6	27.6	27.5	27.6
10	23.6	30.0	30.0	30.0	30.0
11	27.2	27.6	27.6	27.6	27.6
12	56.6	133.1 <sup>2</sup>	133.1	133.0	133.1
13	54.4	124.6 <sup>2</sup>	124.6	124.6	124.5
14	32.1	31.1	31.1	31.1	31.1
15	76.3	76.3	78.5	78.4	78.4
16	137.3	137.3	137.4	137.4	137.3
17	119.1	119.1	118.7	118.7	118.7
18	152.1	152.1	152.3	152.2	152.3
19	117.7	117.7	117.5	117.4	117.5
20	164.2	164.2	164.2	164.2	164.2
21	18.8	18.8	18.9	18.8	18.9
22	20.8	20.8	20.3	20.1	20.2
23	22.6	22.6	22.5	22.4	22.5
24	16.7	16.7	16.1	15.9	16.0
25	18.4	18.4	17.4	17.3	17.5
26	-	-	-	-	-
27	14.2	14.2	14.7	14.7	14.7

#### Refences and comments:

1. G. Höfle et al. Angew. Chem. Int. Ed. Engl. 1996, 35, 1567-1569.
2. Misalignment corrected.
3. Spectrum taken Nov. 15, 1969 from a sample of epothilone C isolated Oct./Nov. 1996.
4. Recent sample of epothilone C.
5. Opposition Statement, p. 55-56.

#### Conclusion:

Chemical shifts for C1-C8 and C15-C27 in EP-B-1186606 are identical with epothilone A (red), those for C9-C14 are identical with epothilone C (blue).



Tab. 2  $^{13}\text{C}$ -NMR chemical shifts of epothilone D in DMSO- $\text{D}_6$ 

C-Atom	Epo B <sup>1</sup>	EP-B-1186606 18.11.96	GBF <sup>2</sup> 14.11.96	GBF <sup>3</sup> 17.5.05	Kosan <sup>4</sup>
1	170.1	170.1	170.1	170.1	170.3
2	38.2	39.0	39.0	38.7	39.1
3	70.0	70.8	70.8	70.8	71.0
4	53.2	53.2	53.3	53.3	53.5
5	217.4	217.4	217.4	217.5	217.6
6	44.9	44.4	44.4	44.4	44.7
7	75.5	75.5	75.4	75.4	75.6
8	35.6	36.3	36.3	36.3	36.5
9	29.6	29.9	29.9	29.9	30.1
10	23.0	25.9	25.9	25.9	26.1
11	32.1	31.6	31.6	31.6	31.8
12	61.0	138.3	138.4	138.4	138.6
13	61.5	120.3	120.3	120.3	120.5
14	33.0	31.9	31.9	31.9	32.1
15	76.6	76.6	78.9	79.0	79.1
16	137.2	137.2	137.6	137.6	137.8
17	119.2	119.2	118.8	118.8	119.0
18	152.1	152.1	152.2	152.3	152.5
19	117.7	117.7	117.4	117.4	117.7
20	164.3	164.3	164.2	164.2	164.4
21	18.9	18.9	18.8	18.9	19.1
22	19.7	19.7	19.9	19.9	20.1
23	22.5	22.5	22.5	22.6	22.7
24	16.4	16.4	16.1	16.2	16.4
25	18.4	18.4	17.7	17.7	17.9
26	22.1	22.9	22.9	23.0	23.2
27	14.1	14.1	14.5	14.6	14.7

**Refences and comments:**

1. G. Höfle et al. Angew. Chem. Int. Ed. Engl. 1996, 35, 1567-1569.
2. Spectrum taken Nov. 14, 1969 from a sample of epothilone D isolated Oct./Nov. 1996.
3. Recent sample of epothilone D.
4. Opposition Statement, p. 59-60.

**Conclusion:**

Chemical shifts for C7, C15-C25, and C27 in EP-B-1186606 are identical with epothilone B(red), those for C1-C6, C8-C14, and C26 are identical with epothilone D (blue).

Tab. 3  $^1\text{H}$ -NMR chemical shifts of epothilone C in DMSO- $\text{D}_6$ 

H-Atoms	Epo A <sup>1</sup>	EP-B-1186606 18.11.96	GBF <sup>2</sup> 15.11.96	Kosan <sup>3</sup>
2a	2.38	2.38	2.35	2.35
2b	2.50	2.50	2.41	2.43
3	3.97	3.97	4.11	4.14
3OH	5.12	5.12	5.10	-
6	3.07	3.07	3.08	3.10
7	3.49	3.49	3.48	3.51
7OH	4.46	4.46	3.18	-
8	1.34	1.34	1.35	1.38
9a	1.15	1.15	1.03	1.05
9b	1.40	1.40	1.55	1.56
10a	1.15	1.15	1.15	1.19
10b	1.46	1.35	1.35	1.37
11a	1.35	1.90	1.88	1.90
11b	1.66	2.18	2.21	2.22
12	2.84	5.38	5.44	5.48
13	3.06	5.44	5.39	5.40
14a	1.76	2.35	2.15	2.14
14b	2.10	2.70	2.70	2.71
15	5.27	5.27	5.12	5.10
17	6.50	6.50	6.50	6.52
19	7.35	7.35	7.33	7.34
21	2.65	2.65	2.65	2.67
22	0.94	0.94	0.91	0.93
23	1.21	1.21	1.20	1.21
24	1.06	1.06	1.06	1.06
25	0.90	0.90	0.89	0.88
26	-	-	-	-
27	2.10	2.10	2.12	2.14

## References and comments:

1. G. Höfle et al. Angew. Chem. Int. Ed. Engl. 1996, 35, 1567-1569.
2. Spectrum taken Nov. 15, 1996 from a sample of epothilone C isolated Oct./Nov. 1996.
3. Opposition Statement, p. 54-55.

## Conclusion:

Chemical shifts for C1-C8 and C15-C27 in EP-B-1186606 are identical with epothilone A, those for C9-C14 are identical with epothilone C.

Tab. 4  $^1\text{H}$ -NMR chemical shifts of epothilone D in  $\text{DMSO-D}_6$ 

H-Atoms	Epo B <sup>1</sup>	EP-B-1186606 18.11.96	GBF <sup>3</sup> 15.11.96	Kosan <sup>5</sup>
2a	2.35	2.35	2.32	2.34
2b	2.38	2.38	2.37	2.34
3	4.10	4.10	4.15	4.14
3OH	5.08	5.08	5.10	-
6	3.11	3.11	3.09	3.09
7	3.48	3.48	3.49	3.48
7OH	4.46	4.46	3.18	-
8	1.29	1.29	1.34	1.33
9a	1.14	1.14	1.15	1.15
9b	1.38	1.38	1.35	1.35
10a	1.14	1.14	1.02	1.02
10b	1.43	1.35	1.65	1.65
11a	1.31	1.75	1.76	1.75
11b	1.61	2.10	2.30	2.29
12	-	-	-	-
13	2.84	5.08	5.10	5.14
14a	2.05	2.30	2.12	2.12
14b	1.84	2.65	2.66	2.66
15	5.29	5.29	5.10	5.09
17	6.51	6.51	6.48	6.48
19	7.35	7.35	7.33	7.33
21	2.65	2.65	2.65	2.65
22	0.90	0.90	0.90	0.90
23	1.19	1.19	1.18	1.18
24	1.07	1.07	1.08	1.08
25	0.91	0.91	0.91	0.91
26	1.19	1.63	1.64	1.64
27	2.11	2.11	2.11	2.11

## References and comments:

1. G. Höfle et al. Angew. Chem. Int. Ed. Engl. 1996, 35, 1567-1569.
2. Spectrum taken Nov. 15, 1969 from a sample of epothilone D isolated Oct./Nov. 1996.
3. Recent sample of epothilone D.
4. Opposition Statement, p. 58-59.

## Conclusion:

Chemical shifts for C1-C8 and C15-C27 in EP-B-1186606 are identical with epothilone A, those for C9-C14 are identical with epothilone C.

Tab. 5  $^{13}\text{C}$ -NMR chemical shifts of epothilone isomers (Epos) with molecular mass  $m/z = 477$  in  $\text{CDCl}_3$

Nr.	Epo C Hardt <sup>1</sup>	trans-Epo C Schinzer <sup>2</sup>	trans-Epo C Danishefsky <sup>3</sup>	Epo D <sub>1</sub> Hardt <sup>1</sup>	Epo D <sub>2</sub> Hardt <sup>1</sup>	Max. delta > 2.0
1	220.6	219.9	219.9	217.0	216.8	3.6
2	170.4	170.5	170.5	169.7	170.4	
3	165.0	164.9	165.0	165.0	165.9	
4	152.1	152.1	152.0	152.2	152.3	
5	138.7	137.1	137.1	138.5	139.8	
6	133.5	134.3	134.4	137.7	137.5	4.2
7	125.0	125.7	125.7	120.7	120.5	4.8
8	119.5	119.8	119.8	121.1	119.2	
9	115.8	116.0	116.0	116.3	116.3	
10	78.5	76.6	77.6	78.8	80.8	2.3
11	74.2	75.8	75.8	77.2	74.3	3.0
12	72.4	72.4	72.4	67.7	69.7	4.7
13	53.4	52.5	52.5	52.5	48.6	4.8
14	41.8	43.6	43.6	46.5	48.4	6.6
15	39.3	38.9	38.8	30.6	39.9	
16	38.6	37.7	37.8	37.6	36.6	
17	31.8	36.2	36.2	32.3	32.7	4.4
18	31.5	32.4	32.5	31.8	32.2	
19	27.6	30.5	30.6	29.5	30.9	3.3
20	27.5	27.2	27.3	25.5	26.0	
21	22.7	21.0	21.0	22.1	23.6	
22	19.1	20.7	20.7	19.2	19.2	
23	18.7	19.1	19.0	16.6	17.1	2.1
24	15.9	16.4	16.4	15.5	15.4	
25	15.5	15.7	15.7	14.5	12.7	2.8
26	13.5	14.8	14.8	9.7	12.4	3.8

#### Refences and comments:

5. I. H. Hardt et al., *J. Nat. Prod.* **2001**, *64*, 847- 856.
6. D. Schinzer et al., *Chem. Eur. J.* **1999**, *5*, 2483- 2491.
3. PCT/US97/22381; D. Meng et al., *J. Am Chem. Soc.* **1997**, *119*, 10073- 10092.

#### Conclusion:

Chemical shifts for individual carbon atoms vary by 2.1 up to 6.6 ppm.

Tab. 6  $^{13}\text{C}$ -NMR chemical shifts of epothilone C in  $\text{CDCl}_3$ 

Nr.	Danishefsky Original <sup>1,2</sup>	Danishefsky Corrected <sup>3</sup>	Nicolaou <sup>4</sup>	Nicolaou <sup>5</sup>	Schinzer <sup>6</sup>	Hardt <sup>7</sup>	Max. delta
1	226.5	220.4	220.6	220.2	220.5	220.6	0.4
2	176.5	170.4	170.4	170.6	170.3	170.4	0.3
3	171.1	165.0	165.0	165.4	165.0	165.0	0.4
4	158.2	152.1	151.9	153.8	152.0	152.1	1.8 <sup>9</sup>
5	144.7	138.6	138.7	139.2	138.6	138.7	0.6
6	139.6	133.5	133.4	134.1	133.4	133.5	0.7
7	131.1	125.0	125.0	126.1	125.0	125.0	1.1
8	125.7	119.6	119.4	120.4	119.5	119.5	1.0
9	122.0	115.9	115.8	116.9	115.8	115.8	0.8
10	84.6	78.5	78.4	79.2	78.4	78.5	0.8
11	80.2	74.1	74.1	74.9	74.1	74.2	0.8
12	78.6	72.5	72.3	73.2	72.4	72.4	0.9
13	59.4	53.3	53.3	54.2	53.3	53.4	0.9
14	47.9	41.8	41.7	42.5	41.8	41.8	0.8
15	45.4	39.3	39.2	40.3	39.2	39.3	1.1
16	44.6	38.5	38.5	39.5	38.5	38.6	1.0
17	38.5	32.4	32.4	32.9	32.5	31.8	1.1
18	37.9	31.8	31.7	32.6	31.7	31.5	1.1
19	33.7	27.6	27.6	28.6	27.6	27.6	1.0
20	33.6	27.5	27.4	28.4	27.5	27.5	1.0
21	28.7	22.6	22.7	23.3	22.7	22.7	0.7
22	25.1	19.0	19.0	19.3	19.0	19.1	0.3
23	25.0	18.9	18.6	19.1	18.7	18.7	0.5
24	21.9	15.8	15.9	16.4	15.8	15.9	0.6
25	21.7	15.6	15.5	16.3	15.5	15.5	0.8
26	19.6	13.5	13.5	14.4	13.5	13.5	0.9

## References and comments:

1. PCT/US97/22381
2. D. Meng et al., J. Am Chem. Soc. 1997, 119, 10073- 10092.
3. Offset of 6.1 ppm.
4. K. C. Nicolaou et al., J. Amer. Chem. Soc. 119, 7960, 1997.
5. K. C. Nicolaou et al., J. Amer. Chem. Soc. 119, 7974, 1997.
6. D. Schinzer et al., Chem. Eur. J. 1999, 5, 2483- 2491.
7. I. H. Hardt et al., J. Nat. Prod. 2001, 64, 847- 856.

## Conclusion:

Chemical shifts for individual carbon atoms vary by 0.3 up to 1.1 ppm.

Tab. 7  $^1\text{H}$ - and  $^{13}\text{C}$ -NMR chemical shifts of epothilones C and D in  $\text{DMSO-D}_6$

Epothilone C					Epothilone D				
H - Atoms	GBF <sup>1</sup> 1. 9. 05	Kosan <sup>2</sup>	GBF <sup>1</sup> 1. 9. 05	Kosan <sup>2</sup>	C - Atoms	GBF <sup>1</sup> 1. 9. 05	Kosan <sup>2</sup>	GBF <sup>1</sup> 1. 9. 05	Kosan <sup>2</sup>
2a	2.34	2.35	2.34	2.34	1	170.08	170.6	170.08	170.3
2b	2.41	2.43	2.37	2.34	2	38.77	38.8	38.88	39.1
3	4.11	4.14	4.14	4.14	3	70.85	70.8	70.81	71.0
3OH	5.10	-	5.08	-	4	53.18	53.2	53.27	53.5
6	3.08	3.10	3.09	3.09	5	217.50	217.5	217.41	217.6
7	3.49	3.51	3.48	3.48	6	44.28	44.3	44.46	44.7
7OH	-	-	4.41	-	7	75.13	75.1	75.45	75.6
8	1.36	1.38	1.34	1.33	8	36.53	36.6	36.29	36.5
9a	1.03	1.05	1.15	1.15	9	27.57	27.6	29.91	30.1
9b	1.55	1.56	1.35	1.35	10	29.98	30.0	25.91	26.1
10a	1.16	1.19	1.01	1.02	11	27.57	27.6	31.57	31.8
10b	1.36	1.37	1.66	1.65	12	133.08	133.1	138.38	138.6
11a	1.89	1.90	1.76	1.75	13	124.54	124.5	120.28	120.5
11b	2.21	2.22	2.30	2.29	14	31.05	31.1	31.86	32.1
12	5.47	5.48	-	-	15	78.44	78.4	78.94	79.1
13	5.38	5.40	5.15	5.14	16	137.35	137.3	137.57	137.8
14a	2.15	2.14	2.12	2.12	17	118.72	118.7	118.78	119.0
14b	2.69	2.71	2.66	2.66	18	152.25	152.3	152.24	152.5
15	5.13	5.10	5.09	5.09	19	117.48	117.5	117.46	117.7
17	6.50	6.52	6.48	6.48	20	164.20	164.2	164.19	164.4
19	7.33	7.34	7.34	7.33	21	18.86	18.9	18.85	19.1
21	2.65	2.67	2.66	2.65	22	20.22	20.2	19.93	20.1
22	0.91	0.93	0.90	0.90	23	22.49	22.5	22.54	22.7
23	1.19	1.21	1.18	1.18	24	16.01	16.0	16.24	16.4
24	1.06	1.06	1.08	1.08	25	17.37	17.5	17.71	17.9
25	0.89	0.88	0.91	0.91	26	-	-	22.95	23.2
26	-	-	1.64	1.64	27	14.65	14.7	14.52	14.7
27	2.12	2.14	2.11	2.11					

References and comments:

- 1 New isolates from So ce90wild DSM 6773 (produced August 1 – 31, 2005).
- 2 Opposition Statement, p. 54-55.

**Conclusion:** All signals for GBF and Kosan samples are identical within the experimental error. The maximal shift differences of 0.52 and 0.22 ppm are observed for C-1.